

N O T I C E

THIS DOCUMENT HAS BEEN REPRODUCED FROM
MICROFICHE. ALTHOUGH IT IS RECOGNIZED THAT
CERTAIN PORTIONS ARE ILLEGIBLE, IT IS BEING RELEASED
IN THE INTEREST OF MAKING AVAILABLE AS MUCH
INFORMATION AS POSSIBLE

CR-160609

SPT



(NASA-CR-160609) CARDIOPULMONARY DATA
ACQUISITION SYSTEM. VERSION 2.0, VOLUME 2:
DETAILED SOFTWARE/HARDWARE DOCUMENTATION
Final Report (Technology, Inc., Houston,
Tex.) 195 p dC A09/MF A01

880-33084

Uncias
33201

CSSL 05B G3/52



TECHNOLOGY INCORPORATED

LIFE SCIENCES DIVISION

FINAL REPORT

ON THE

CARDIOPULMONARY DATA ACQUISITION SYSTEM



VERSION 2.0

VOLUME 2

DETAILED SOFTWARE/HARDWARE DOCUMENTATION

Prepared for the NASA Johnson Space Center
Cardiovascular Research Laboratory

MSC-18783



February 18, 1980

P. O. BOX 58827

HOUSTON, TEXAS 77058

FINAL REPORT
ON THE
CARDIOPULMONARY DATA ACQUISITION SYSTEM
VERSION 2.0
VOLUME 2
DETAILED SOFTWARE/HARDWARE DOCUMENTATION

Prepared for the NASA Johnson Space Center
Cardiovascular Research Laboratory

February 18, 1980

Contract NAS7-14880
Project 0160-20


TECHNOLOGY INCORPORATED
LIFE SCIENCES DIVISION
16821 Buccaneer, Suite 206
Houston, Texas 77058

Table of Contents

<u>Volume I (User's Guide)</u>	<u>Page</u>
Title Page -----	i
Approval Sheet -----	ii
Abstract -----	iii
CDAS Development -----	iv
Table of Contents -----	v
Overview of the Cardiopulmonary Data Acquisition System -----	1
Purpose and Philosophy -----	2
Hardware Configuration -----	3
Software Description -----	3
System Operation -----	5
Alternate Operation -----	9
Future Expansion -----	10
User's Manual -----	13
General Information -----	14
PIP (Peripheral Interchange Program) -----	19
Initializing a Floppy Disk -----	19
Getting a Floppy Directory (Contents) Listing -----	20
Deleting Files -----	22
Renaming a File -----	22
Copying Data Files (Making Disk Backups) -----	23
Terminating PIP -----	26
PRETEST (Pretest Subject & Test Information) -----	27
STRESS (Stress Test Data Collection Program) -----	29
REPORT (Test Summary Report & Data Editing) -----	39
PLOT (Test Data Plotting Program) -----	43
EDICON (Edit/Change Calibrations Reference Constants) ---	45
Appendix A: Trouble Shooting Chart -----	47
Appendix B: Sample Hard Copy Output -----	53
Appendix C: Equations Used in the Software -----	65
Appendix D: A/D Sampling Specifications & Channel Assignments ----	69
Appendix E: Notes on Operation of CDAS from Tape Playback -----	70
<u>Volume II (Detailed Software/Hardware Documentation)</u>	
Appendix F: Wiring and Timing Diagrams -----	1
General -----	2
LSI-11 Computer Control Panel&Interface Cables -----	9
Appendix G: Flowcharts of the Software -----	
PRETEST -----	
STRESS -----	
EDIT -----	
PLOT -----	
Appendix H: Program Listings -----	
PRETEST -----	
STRESS -----	
REPORT -----	
PLOT -----	
EDICON -----	
Appendix I: Format of the Floppy Disk File -----	

APPENDIX F
WIRING AND TIMING DIAGRAMS

CDAS - CARDIOPULMONARY LABORATORY
SIMPLIFIED HARDWARE CONFIGURATION



**JSC CARDIOPULMONARY LABOPATORY
STRESS TEST TAPE LOG FORM**

Analog Tape # _____ Subject Name: _____ Subject # _____
 Date: _____ DOY: _____ SS# _____ AGE: _____ HT.(cm) _____ WT.(kg) _____

CALIBRATION DATA

Analog Tape Channel	Data	Calib. Value (low/high)
3	VCG-X	} 1 mV
5	VCG-Y	
7	VCG-Z	
2	HEART RATE	60/180 bpm
4	{ ELEVATION	8/32%
	{ WORK LOAD	60/240 watts
6	SPIRO. VOLUME	0/8 liters
9	{ TREAD. SPEED	2/8 mph
	{ BIKE SPEED	20/80 rpm
11	SYSTOLIC BP	50/200mmHg
13	DIASTOLIC BP	50/200mmHg
10	O ₂	0/ _____ %
14	N ₂	0/ _____ %
12	CO ₂	0/ _____ %
1	TAPE REF. (Direct)	
8	IRIG-B TIME (Direct)	

CALIBRATION START TIMES:

ECG: _____
 Low: _____
 High: _____
 Ambient Air: _____

AMBIENT CONDITIONS:

O₂: _____ %
 N₂: _____ %
 CO₂: _____ %
 Temperature: _____ ° C
 Bar. Pres.: _____ mmHg

TEST TYPE

(Circle One) Bicycle Ergometer Treadmill (Modified Balke) Treadmill (Modified Bruce)

Begin Rest: _____ (10 min) Begin Exercise: _____

Begin Recovery: _____ (8 min) Begin Recovery: _____

Comments/Anomalous Conditions: _____

CDAS INPUT/OUTPUT CONNECTIONS

VCG, BPMS, Front Panel, and Analog Interface

BPMS SBP analog output	—————→	to analog interface SBP input
BPMS DBP analog output	—————→	to analog interface DBP input
BPMS Korothoff Sounds output	—————→	to analog interface K-Sounds input
BPMS Pressure analog output	—————→	(not used)
<hr/>		
Spirometer controller output	—————→	to analog interface spir. vol. input
<hr/>		
Treadmill speed & elevation outputs	} —————→	to front panel
Bicycle ergometer speed & WL outputs		treadmill/bicycle switch
<hr/>		
Treadmill/Bicycle switch outputs	—————→	to analog interface speed & work load/ elevation inputs
<hr/>		
LSI-11 D/A outputs (0 or 1)	—————→	to ergometer work load control input
<hr/>		
VCG system X lead output	} —————→	to front panel VCG switch
VCG Y lead output		
VCG Z lead output		
<hr/>		
Front panel VCG trigger select switch output	—————→	to ECG trigger inputs of BPMS and cardiotachometer
<hr/>		
Cardiotachometer output	—————→	to analog interface HR input and to front panel DPM switch
<hr/>		
All analog interface outputs	} —————→	to DPM/Scope Monitor switches
Front panel BNC DPM/Scope Inputs		
<hr/>		
DPM/Scope Monitor switches outputs	—————→	to front panel DPM select switch and oscilloscope vertical inputs
Front panel DPM switch outputs	—————→	to DPM
<hr/>		
All analog interface outputs	} —————→	to patch panel (see Chart)
Mass spectrometer interface outputs		
<hr/>		
T.C. generator IRIG-B code out	—————→	to analog tape chan. 8 record input
Tape recorder chan. 8 reprod. output	—————→	to T.C. generator/reader code in

Abbreviations:

VCG	Vectorcardiograph
BPMS	Blood Pressure Measuring System (Sphygmometrograph)
WL	Work Load (of bicycle ergometer)
D/A	Digital to Analog output from LSI-11
T.C.	Time Code
DPM	Digital Panel Meter

CDAS Relay Patch Panel, Analog Tape Recorder, LSI-11 A/D Inputs


Signal Name	Signal Source	Patch Panel Relay Pole No.	Tape Rec. Playback Ch. (for Relay Input #2)	Relay Output to Buffer Amp. Input #	Buffer Amp. Output to A/D Input #	Strip Chart Chan. No. (for Tape Output)
Spir. Vol.	Spiro. Contr.	0	6	10	0	5
O ₂	Mass	1	10	1	1	
N ₂	Spectrometer	2	14	2	2	
CO ₂	Interface	3	12	3	3	
Ht. Rate	Analog Interface Outputs	4	2	4	4	
ML/Elev.		5	4	5	5	6
Speed		6	9	6	6	7
SBP		7	11	7	7	
DBP		8	13	8	8	
(Spare)		9				
(Spare)		10				
24v for RT/PB lights	Pwr. Supply	11				
Tape Spd.Ref.	Tape Rec.	<div style="display: flex; align-items: center; justify-content: center;"> } <div> not switched </div> </div>	1*			1
VCG-X	VCG		3*			2
VCG-Y	System		5*			3
VCG-Z			7*			8**
TRIG-B	Time Code Generator		8*			
Time Code						

NOTES:

- (1) Except for patch panel relay pole #11:
- All relay inputs are from outputs of analog and mass spectrometer interfaces
- All relay #2 inputs are from analog tape recorder playback (reproduce) outputs.

- (2) Patch panel relay pole #11 is used to control the Real Time-Playback lights, which indicate the status of the latching relay. The pole of the switch is connected to the 24v. supply, the #1 contact goes to the Real Time light, and the #2 contact goes to the Playback light.
- (3) All cable shields (signal returns) are switched along with the cable center conductors, and all isolated from each other at the patch panel.
- (4)* These signals are not switched by the patch panel relay.
- (5)** Slow code, not IRIG-B, is sent to strip chart recorder channel 8.

COAX LATCHING RELAY PATCH PANEL



LSI-11 COMPUTER CONTROL PANEL AND INTERFACE CABLES

THIS PAGE PLANK NOT FILMED

CDAS MICROCOMPUTER HARDWARE REQUIREMENTS

<u>Supplier</u>	<u>Model No.</u>	<u>Description</u>
MDB	MLSI-KD11-HA	LSI-11 CPU
MDB	MLSI-KEV11	Extended Arithmetic Option
MDB	MLSI-MSV11-DD	32K x 16 bit RAM, on-board refresh
MDB	MLSI-DRV11C	Parallel Interface
MDB	MLSI-SMU	System Monitoring Unit/Line Time Clock/ Front Panel Controls/Bus Terminator
MDB	MLSI-BPA84	Backplane - Card Guide Assembly
MDB	MLSI-BA11-000	Enclosure
MDB	MDB-250-T-5/12	Power Supply
MDB	MLSI-DLV11	Asynchronous Interface
MDB	MLSI-KW11P	Programmable Real-Time Clock
ADAC	1030-16PD-A- 3PGA-2-A-P	Analog Data Acquisition and Control System
Data Systems	DSD 210-L11-2	Dual Flexible Disk System
DEC	LA36	Decwriter II Hard Copy Terminal

CVDAS Input/Output Signals - Parallel Interface

Input Register - Address 167774

Input Register Lines	Function
IN 00	CALIB. Switch
IN 01	REST Sw.
IN 02	EXERCISE Sw.
IN 03	RECOVERY Sw.
IN 04	PAUSE Sw.
IN 05	END TEST Sw.
IN 06	FVC Sw.
IN 07	(Spare Switch)
IN 08	} not used
IN 09	
IN 10	
IN 11	
IN 12	
IN 13	
IN 14	
IN 15	

Output Register - Address 167772

Output Register Lines	Function
OUT 00	CALIB. Light
OUT 01	REST Light
OUT 02	EXERCISE Light
OUT 03	RECOVERY Light
OUT 04	PAUSE Light
OUT 05	END TEST Light
OUT 06	FVC Light
OUT 07	(Spare Light)
OUT 08*	TEST TYPE 0
OUT 09*	TEST TYPE 1
OUT 10	Decr. TM Elev.
OUT 11	Incr. TM Elev.
OUT 12	Incr. TM Speed
OUT 13	Decr. TM Speed
OUT 14	not used
OUT 15	not used

*See below

Control/Status Register - Address 167770

Control Lines (used as outputs)

CSR 1	CSR 0	Operating Mode
0	0	Normal (Operate)
0	1	Low Calib.
1	0	High Calib.
1	1	Ambient Air/Standby

Test Type Control Lines

OUT 09	OUT 08	Selected Test Type
0	0	Treadmill
0	1	Bicycle Ergometer
1	0	LBNP
1	1	Other

LSI-11 SERIAL INTERFACE CABLE

Computer 40-Pin Board Connector		Agile A1 Terminal 25-Pin EIA Connector	
Signal	Pin No.	Signal	Pin No.
EIA Received Data In	8	Protective Ground	1
Transmitted Data	6	Transmitted Data	2
Signal Ground	1,2	Received Data	3
Data Set Ready	22	Signal Ground	7
Data Terminal Ready	26	Data Terminal Ready	20
		Carrier Detect	8

LSI-II ANALOG/DIGITAL CONVERTER CABLE

Signal Name	A/D Board Connector C1 Pin No.	Analog Cable 50-Pin Connector J1 Pin No.	Coax Cable BNC Connector Designation
A/D Ch. 0 in	31	16	0
A/D Ch. 1 in	29	15	1
A/D Ch. 2 in	27	14	2
A/D Ch. 3 in	25	13	3
A/D Ch. 4 in	23	12	4
A/D Ch. 5 in	21	11	5
A/D Ch. 6 in	19	10	6
A/D Ch. 7 in	17	9	7
A/D Ch. 8 in	32	33	8
A/D Ch. 9 in	30	32	9
A/D Ch. 10 in	28	31	10
A/D Ch. 11 in	26	30	11
A/D Ch. 12 in	24	29	12
A/D Ch. 13 in	22	28	13
A/D Ch. 14 in	20	27	14
A/D Ch. 15 in	18	26	15
D/A Ch. 1 out	37	46	D1
D/A Ch. 1 ret.	38	45	(D1 Shield)
D/A Ch. 2 out	39	44	D2
D/A Ch. 2 ret.	40	43	(D2 Shield)
Amp. Lo in (Source Ret.)	33	50	(A/D Input Shields)
Power Return (Chassis Grd.)	34	49	N.C.

Related Equipment

Seperate reports are available on the following items, which are used with the LSI-11 CVDAS system for tests in the JSC Cardiovascular Laboratory:

Microprocessor - Based Cardiotachometer

Analog Interface and Control System for the
JSC Cardiovascular Laboratory

CIA	CIS	CI6
		•+5 •DOWN UP• •FAST SLOW• •IFAST ISLOW• •IUP IDOWN•

CI1	CI2	CI3
•SW •SW •SW •SW •SW •SW	•SW •SW •SW •SW •SW •SW	•PVC •END •PA •REL •PA •PAC •CAL

GROUP B


A	B	C	D	E	F	G	H	I	J
Z4 7414	Z5 7486	Z6 7400	Z7 7438						
555 26									

GROUP A


A	B	C	D	E	F	G	H	I	J
Z1 7404	Z2 7404	Z3 7404							

SOLDOUT FRAME


HOLDOUT FRAME
2

		TECHNOLOGY INCORPORATED LIFE SCIENCES DIVISION HOLDOUT FRAME	
DRAWN BY	DATE	TITLE	
M. H. UTLEY	10/16/79		
DESIGN ENG.		CVDAS - CONTROL PANEL CIRCUIT BOARD	
PROD. ENG.		DWG. NO. TH7942-1A01	
		CHECKED 10/16/79	

C1




C2



WOLDOUT FRAME

C3



C4







EXHIBIT PRATT

2

			<p>TECHNOLOGY INCORPORATED LIFE SCIENCES DIVISION HOUSTON, TEXAS 77058</p>		
<p>DRAWN BY M.H. UTLEY</p>		<p>DATE 10/25/71</p>		<p>TITLE CVDAS - CONTROL PANEL COMPONENT LAYOUT</p>	
<p>DESIGN ENG.</p>					
<p>PROD. ENG.</p>				<p>DWG. NO. TH7942-1A02</p>	
				<p>REV. 1</p>	



REVISION PLANT



FOLDOUT FRAME

2



TECHNOLOGY INCORPORATED
LIFE SCIENCES DIVISION
HOUSTON, TEXAS

DRAWN BY
M. H. UTLEY

DATE
10-27-71

TITLE
CYDAS CONTROL PANEL
PUSHBUTTON SWITCH
DEBOUNCE CIRCUIT







DESIGN ENG.

DESIGNED BY



TH7742-1E01

DATE

1 OF 1




FOLDOUT FRAME



BOLDOUT FRAME

2

NOTES: (1) COMMON SIDE OF ALL LIGHTS IS
CONNECTED TO +24V SUPPLY OUTPUT.

			TECHNOLOGY INCORPORATED LIFE SCIENCES DIVISION HOUSTON, TEXAS 77058	
DRAWN BY M.H. UTLEY		DATE 10-25-79	TITLE CVDAS- CONTROL PANEL PUSHBUTTON LIGHT DRIVERS	
DESIGN ENG				
PROJ. ENG			DWG. NO. TH7942-1E02	
			SHEET 1 OF 1	

LSI-
PARALL
INTER
CONNE

CIRCUIT
BOARD
CONNECTORS

CHASSIS
CONNECTOR

INTERCONNECTING
RIBBON
CABLE

	C12	P12	J1		
CAL SWITCH OUT	6	6	1	1	TT
REST SWITCH OUT	5	5	2	2	LL
EX SWITCH OUT	4	4	18	18	HE
REC. SWITCH OUT	3	3	8	8	EE
PAUSE SWITCH OUT	2	2	3	3	KK
END SWITCH OUT	1	1	4	4	HH
ENC SWITCH OUT	9	9	6	6	EE
SPARE SWITCH OUT	10	10	7	7	CC
			9	9	Z
			10	10	V
			11	11	W
			12	12	V
			13	13	U
			15	15	P
			16	16	N
			17	17	M
GROUND 0			32	32	S

ROLLOUT FRAME
1

LSI-II
PARALLEL
INTERFACE
CONNECTOR

↓ J2

TING

TT IN 00

LL IN 01

HE IN 02

BB IN 03

KK IN 04

HH IN 05

EE IN 06

CC IN 07

Z IN 08

Y IN 09

W IN 10

V IN 11

U IN 12

P IN 13

N IN 14

M IN 15

S GROUND (COMMON, 31)

FOLDOUT FRAME
2



TECHNOLOGY INCORPORATED
LIFE SCIENCES DIVISION
HOUSTON, TEXAS 77058

DRAWN BY
M.H. UTLEY

DATE
10-26-79

TITLE
CVDAS CONTROL PANEL
I/O INTERCONNECT
PARALLEL INTERFACE INPUTS

DESIGN ENG

PROD ENG

DWG. NO.
TH 7942-1E03

SHEET
1 OF 1


LSI-11
PARALLEL
INTERFACE
CONNECTOR


CONNECTING
CABLE


	J1
C	OUT 00
K	OUT 01
RR, NN	OUT 02
U	OUT 03
L	OUT 04
N	OUT 05
R	OUT 06
T	OUT 07
Z	OUT 10
AA	OUT 11
BB	OUT 12
EE	OUT 13
HH	OUT 14
JJ	OUT 15
W	OUT 08
X	OUT 09
K	CSR 0 (ON CONV. J2)
DD	CSR 1
S	GROUND


FOLDOUT FRAME

2

			<p>TECHNOLOGY INCORPORATED LIFE SCIENCES DIVISION HOUSTON, TEXAS 77058</p>		
<p>DRAWN BY M.H. UTLEY</p>		<p>DATE 10/25/76</p>		<p>TITLE CVDAS- CONTROL PANEL I/O INTERCONNECT PARALLEL INTERFACE OUTPUTS</p>	
<p>DESIGN ENG.</p>		<p>PROJ. ENG.</p>		<p>DWG. NO. TH7942-1E04</p>	
<p></p>		<p></p>		<p>SHEET 1 OF 1</p>	




21



FOLDOUT FRAME

ADRIAL
ROLLER
SECTOR TREADMILL
CONTROLLER




5 ORG TO SPEED INTERFACE



2 YEL TO SPEED INTERFACE

6 GRN TO ELEV. INTERFACE

3 BLU TO ELEV. INTERFACE

NOT TO BE DRAWN
2

 TECHNOLOGY INCORPORATED LIFE SCIENCES DIVISION HOUSTON, TEXAS 77058		
DRAWN BY M.H.UTLEY	DATE 10/26/79	TITLE CYDAS-CONTROL PANEL TO TREADMILL INTERCONNECT
DESIGN ENG.		
PROJ. ENG.		DWG NO TH7942-1E06
		SHEET 1 OF 1





FOLDOUT FRAME

A YEL
 B ORG.
 C RED
 D BRW.
 E BLU
 F GRN
 G GRAY
 H VIO
 J BLK
 K WHT

FASTER ORG - PIN 5
 SLOWER YEL - PIN 2
 UP GRN - PIN 6
 DOWN BLU - PIN 3
 +5 VOLT VIO - PIN 1
 SPARE BRW - PIN 4


POLYESTER FILM
 2

		TECHNOLOGY INCORPORATED LIFE SCIENCES DIVISION HOUSTON, TEXAS 77058		
		DRAWN BY M. H. UTLEY	DATE 10-24-79	TITLE CVDAS- LAYOUT-TREADMILL CONTROLLER MOD. BOARD
DESIGN ENG.			PROJ. NO TH7942-1A03	
			SHEET 1 OF 1	



INCR


DECR



WELDOUT FRAME

2

NOTE 3(1) +5VDC IS CONNECTED THRU
MANUAL/AUTO SWITCH.



TECHNOLOGY INCORPORATED
LIFE SCIENCES DIVISION
HOUSTON, TEXAS 77058

DRAWN BY
M.H. UTLEY

DATE
10-24-79


TITLE
CVDAS-TREADMILL
CONTROLLER
MODIFICATION

DESIGN ENG.


PROJ. ENG.

DWG. NO.
TH7942-1E07

SHEET
1 OF 1




FOLDOUT FRAME
1



FOLDOUT

2

			TECHNOLOGY INCORPORATED LIFE SCIENCES DIVISION HOUSTON, TEXAS 77058	
DRAWN BY M.H. UTLEY		DATE 9-30-79	TITLE CVDAS- CONTROL PANEL POWER SUPPLY	
DESIGN ENG.				
PROJ. ENG.			DWG. NO. TH7942-1E08	
			SHEET 1 OF 1	

APPENDIX G
FLOWCHARTS OF THE SOFTWARE

Hierarchy Diagram of the CDAS PROGRAMS






Figure 6

Program: PRETEST


PRETEST asks the user the pretest questions about the test subject and stores them in the data file it creates.






Flow Chart of STRESS

STRESS is the main program that tests the CDAS front pannel switches; only those lite are functional.




SUBROUTINE DATA


Subroutine data performs all the disk I/O for CDAS one argument is passed, either OUT or IN




SUBROUTINE CAL
 Calibration is accomplished three phases:
 Phase 1: low calibration
 Phase 2: high calibration
 Phase 3: ambient calibration



SUBROUTINE CAL (Continued)




SUBROUTINE CAL CONTINUED



SUBROUTINE INITI


This subroutine is used by CAL, PFT and EX subroutines to reestablish the value in the ADCOM1 common area used by ADISR.





Flow Chart for SUBROUTINE PFT


SUBROUTINE EX

The exercise subroutine collects the subjects data during the REST, EXERCISE and RECOVERY phases for both treadmill or bicycle test.











**Routine PAUSE
for SUBROUTINE EX**


This routine is
called when the pause
button is pushed.






Subroutine WRKADJ


The work adjustment subroutine regulates the workload on the ergometer to meet the target heart rates.



Subroutine SNATCH


This routine prevents the A/D interrupt routine from updating the heart rate and workload data to maintain accuracy while it is being moved to the WRKADJ routine






SUBROUTINE TRSBUF

This subroutine disables the clock interrupts by altering the program status word then makes a copy of ADCOM1 and GAS30. ADCOM1 and GAS30 are zeroed and certain values in ADCOM1 are reset in their original value. The clock interrupts are again enabled.




SUBROUTINE BREATH

Subroutine BREATH is used to determine the oxygen consumption and carbon dioxide production in each breath exhaled by the subject and keep a 30 sec. sum.




SUBROUTINE RPT30

This routine prints the 30 second reports for the collected data and places the 1 minute reports in the patients data file.




SUBROUTINE REPORT


Prints the correct headers for the 30 sec reports depending if test is a treadmill or bicycle.





Routine ADISR


This program is the interrupt service routine which is triggered by the programmable clock 100 times a second to collect the analog signals from the lab equipment.
















SUBROUTINE LIGHT

This routine is used by subroutine CAL and EX to turn on the appropriate panel lights to cue the operator to the functions allowed at various times during a test.




PROGRAM REPORT

This program reproduces n number of copies of the one minute summary report; and also, allows the user to edit the stored data



EDIT SUBROUTINE

This subroutine allows the user to edit the data stored in the subjects data file stored on disk



APPENDIX H
PROGRAM LISTINGS

```

C *****
C *****
C *
C *
C *          TECHNOLOGY INCORPORATED
C *          LIFE SCIENCE DIVISION
C *
C *****
C *
C * PROGRAM NAME:..... PRETES
C * AUTHOR:..... OGDOTHY STEPHENS
C * DATE:..... 2/NOV/78
C * MODIFIED:..... ABHIJIT GADGIL (DEC/1979)
C *
C *-----*
C *
C * COMPUTER SYSTEM:..... DEC PDP 1103 "LSI-11"
C * OPERATING SYSTEM:..... RT-11 VOIC-03A
C *
C *-----*
C *
C * COMPILING SEQUENCE:
C *
C *      .R FORTRA <CR>
C *      *PRETES=PRETES/L <CR>
C *      **<CR>
C *
C *-----*
C *
C * RUN MODUAL LINKING SEQUENCE:
C *
C *      .R LINK <CR>
C *      *PRETES=PRETES,DX0:SYSLIB/F <CR>
C *      **<CR>
C *
C *-----*
C *
C * CALLING SEQUENCE:
C *
C *      .R PRETEST <CR>
C *
C *-----*
C *
C * PURPOSE:
C *
C *      ALLOWS THE USER TO CREATE A PATIENT DATA FILE ON DX1:
C *      AND ENTER IN THE PERTINENT SUBJECT INFORMATION.
C *
C *****
C *****
C CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C

```

```
C
C      PROGRAM: PRETEST
C
C      PURPOSE: QUESTIONS NECESSARY FOR TREADMILL
C
C      AND BICYCLE SUBJECT DATA
C
C
C
C
```

```
C*****
```

```
0001      REAL CRT(10)
0002      LOGICAL*1 YES,JAUM
0003      REAL*4 IMONTH,MONTH(12)
0004      DATA MONTH/'JAN','FEB','MAR','APR','MAY','JUN','JUL',
0005      &'AUG','SEP','OCT','NOV','DEC'/
0005      COMMON/X/MON(24)
0006      INTEGER NFILE(3),FILE(4),ICRT(1536)
0007      EQUIVALENCE (ICRT(50),CRT(1))
0008      EQUIVALENCE (ICRT(54),TMP),(ICRT(57),PRS)
0009      DATA FILE/3RDX1,3RXXX,3RXXX,3RDAT/
0010      DATA YES/1HY/
C      ZERO OUT THE DATA BLOCKS (6 OF THEM)
0011      DO 800 I=1,1536
0012      800      ICRT(I)=0
C * * * * * SET UP QUEUE ELEMENTS FOR DISK I/O
0013      CALL LOCK
0014      IF(IGSET(3).NE. 0) STOP 'QUEUE ELEMENT FAILURE'
0016      IGOOF=0
0017      TYPE 400
0018      400      FORMAT(' PATIENT FILENAME :',S)
0019      ACCEPT 410,NFILE
0020      410      FORMAT(3A2)
0021      CALL IRAD50(6,NFILE,FILE(2))
C * * * * * OPENING DISK FILE
0022      ICHAN=IGETC()
0023      IF(ICCHAN.LT.0) STOP 'NO AVAILABLE CHANNEL'
0025      IF(IENTER(ICCHAN,FILE,6).LT. 0) STOP 'FILE ALLOCATION FAILURE'
0027      TYPE 199
0028      199      FORMAT('          PRETEST QUESTIONS ')
0029      IF(IGOOF.GT.0) GO TO 121
0031      521      TYPE 21
0032      21      FORMAT(' 1. SUBJECT NAME: ',S)
0033      ACCEPT 22, (ICRT(I),I=26,39)
0034      22      FORMAT(14A2)
0035      IF(IGOOF.GT. 0) GO TO 121
0037      522      TYPE 1
0038      1      FORMAT(' 2. SUBJECT SEX(M=MALE,F=FEMALE): ',S)
0039      ACCEPT 2, ICRT(2)
0040      2      FORMAT(A2)
0041      IF(IGOOF.GT.0) GO TO 121
0043      523      TYPE 23
0044      23      FORMAT(' 3. SOCIAL SEC. NO.(NNN NN NNNN) ')
0045      ACCEPT 3, (ICRT(I),I=3,5)
0046      3      FORMAT(I3,1X,I2,1X,I4)
```

```

0047      IF(IG00F.GT.0) GO TO 121
0049  524      TYPE 33
0050  33      FORMAT(' 4. DATE OF BIRTH (DD-MMM-YY): ')
0051      ACCEPT 4, IDAY, IMONTH, IYR
0052  4      FORMAT( I2, 1X, A3, 1X, I2)
0053      DO 100 I=1, 12
0054      IF(IMONTH.NE.MONTH(I)) GO TO 100
0056      ICRT(22)=I
0057      ICRT(23)=IDAY
0058      ICRT(24)=IYR
0059      GO TO 300
0060  100     CONTINUE
0061      TYPE 200
0062  200     FORMAT(' ERROR IN DATE OF BIRTH. PLEASE RE-ENTER')
0063      GO TO 524
0064  300     CONTINUE
0065      IF(IG00F.GT.0) GO TO 555
0067  525     TYPE 65
0068  65      FORMAT(' 5. HEIGHT (CM): ', $)
0069      ACCEPT 5, CRT(1)
0070  5      FORMAT(F6.0)
0071      IF(IG00F.GT.0) GO TO 121
0073  526     TYPE 266
0074  266     FORMAT(' 6. WEIGHT (KG) : ', $)
0075      ACCEPT 6, CRT(2)
0076  6      FORMAT(F6.1)
0077      IF(IG00F.GT.0) GO TO 121
0079  527     TYPE 67
0080  67      FORMAT(' 7. TEST DATE (DD-MMM-YY): ')
0081      ACCEPT 7, IDAY, IMONTH, IYR
0082  7      FORMAT(I2, 1X, A3, 1X, I2)
0083      DO 500 I=1, 12
0084      IF(IMONTH.NE.MONTH(I)) GO TO 500
0086      ICRT(40)=I
0087      ICRT(41)=IDAY
0088      ICRT(42)=IYR
0089      GO TO 700
0090  500     CONTINUE
0091      TYPE 600
0092  600     FORMAT(' ERROR IN TEST DATE. PLEASE RE-ENTER')
0093      GO TO 527
0094  700     CONTINUE
0095      GO TO 555
0096  31      TYPE 68
0097  68      FORMAT(' 8. UNIQUE NO.(AAAAA): ')
0098      ACCEPT 8, ICRT(7), ICRT(8)
0099  8      FORMAT(A2, I4)
0100      IF(IG00F.GT.0) GO TO 121
0102  528     TYPE 69
0103  69      FORMAT(' 9. RETEST NO.: ', $)
0104      ACCEPT 9, ICRT(9)
0105  9      FORMAT(I2)
0106      IF(IG00F.GT.0) GO TO 121
0108  529     TYPE 70

```

```
0109 70      FORMAT(' 10. CPID NO.(NAN): ')
0110      ACCEPT 10, ICRT(10)
0111 10      FORMAT(I4)
0112      IF(IGOOF.GT.0) GO TO 121
0114 530      TYPE 71
0115 71      FORMAT(' 11. TEST MODE:(1=LIVE,2=TAPE) ',S)
0116      ACCEPT 11, ICRT(11)
0117 11      FORMAT(I2)
0118      IF(IGOOF.GT.0) GO TO 121
0120 531      TYPE 72
0121 72      FORMAT(' 12. TYPE OF TEST : (2=TREADMILL,1=ERGOMETER) ',S)
0122      ACCEPT 12, ICRT(12)
0123 12      FORMAT(I2)
0124      IF(IGOOF.GT.0) GO TO 121
0126 532      TYPE 73
0127 73      FORMAT(' 13. AMBIENT TEMP.(C) : ',S)
0128      ACCEPT 13, TMP
0129 13      FORMAT(F5.0)
0130      IF(IGOOF.GT.0) GO TO 121
0132 533      TYPE 74
0133 74      FORMAT(' 14. AMBIENT PRESSURE(MMHG): ',S)
0134      ACCEPT 14, PRS
0135 14      FORMAT(F6.0)
0136      IF(IGOOF.GT.0) GO TO 121
0138      IF(ICRT(12) .EQ. 2) GO TO 121
0140 534      TYPE 75
0141 75      FORMAT(' 15. TARGET HR @ 4 MIN.(BPM): ',S)
0142      ACCEPT 276, ICRT(15)
0143 276      FORMAT(I4)
0144      IF(IGOOF.GT.0) GO TO 121
0146 535      TYPE 277
0147 277      FORMAT(' 16. TARGET HR @ 8 MIN.(BPM): ',S)
0148      ACCEPT 278, ICRT(16)
0149 278      FORMAT(I4)
0150      IF(IGOOF.GT.0) GO TO 121
0152 536      TYPE 279
0153 279      FORMAT(' 17. TARGET HR @ 12 MIN.(BPM): ',S)
0154      ACCEPT 280, ICRT(17)
0155 280      FORMAT(I4)
0156      IF(IGOOF.GT.0) GO TO 121
0158 537      TYPE 281
0159 281      FORMAT(' 18. TARGET HR @ 16 MIN.(BPM): ',S)
0160      ACCEPT 282, ICRT(18)
0161 282      FORMAT(I4)
0162      IF(IGOOF.GT.0) GO TO 121
0164 538      TYPE 191
0165 191      FORMAT(' 19. MAX TEST ? (1=YES,2=NO): ',S)
0166      ACCEPT 19, ICRT(19)
0167 19      FORMAT(I2)
0168      IF(ICRT(19) .EQ. 2) GO TO 121
0170      ICRT(20)=220-ICRT(25)
0171      TYPE 20, ICRT(20)
0172 20      FORMAT(' 20. TARGET HR @ MAX: ',I4)
0173 50      FORMAT(I2)
```

```
0174 121      TYPE 122
0175 122      FORMAT(' ALL INFORMATION CORRECT ? (Y=YES) ')
0176      ACCEPT 222,JNUM
0177 222      FORMAT(A1)
0178      IF(JNUM.EQ. YES) GO TO 120
0180      IGOOF=IGOOF+1
0181      TYPE 322
0182 322      FORMAT(' TYPE QUESTION NUMBER OF INCORRECT ENTRY ')
0183      ACCEPT 323,INUM
0184 323      FORMAT(I2)
0185      IF(INUM.LE. 0) GO TO 338
0187      IF(INUM.GT.20) GO TO 338
0189      GO TO 490
0190 338      TYPE 339
0191 339      FORMAT(' * * * * * TRY AGAIN, WRONG ENTRY ! ! ! ! ! ')
0192      GO TO 121
0193 490      GO TO (521,522,523,524,525,526,527,31,528,529,530
*,531,532,533,534,535,536,537,538),INUM

C
C      CALCULATE AGE FROM DOB AND TEST DATE
C
C
C
C      MONTH CHECK * * * * *
C
C
0194 555      ICRT(25)=ICRT(42)-ICRT(24)
0195      IF (ICRT(22)-ICRT(40)) 666,664,663
0196 663      ICRT(25)=ICRT(25)-1
0197      GO TO 666
0198 664      IF (ICRT(23).GT.ICRT(41)) GO TO 663
0200 666      TYPE 667, ICRT(25)
0201 667      FORMAT(' AGE: ',I2)
0202      IF(IGOOF.GT.0) GO TO 121
0204      GO TO 31

C * * * * * WRITING TO DISK FILE
0205 120      IF(IWRITW(1536,ICRT,0,ICHAN).LT.0) STOP 'DISK WRITE ERROR'
C * * * * * CLOSING DISK FILE
0207      CALL UNLOCK
0208      CALL CLOSEC(ICHAN)
0209      STOP
0210      END
```


FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
CRT	000150	REAL*4 ARRAY (10)
MONTH	006006	REAL*4 ARRAY (12)
NFILE	006066	INTEGER*2 ARRAY (3)
FILE	006074	INTEGER*2 ARRAY (4)
ICRT	000006	INTEGER*2 ARRAY (1536)
YES	006104	LOGICAL*1 VARIABLE
JNUM	010202	LOGICAL*1 VARIABLE
IMONTH	010204	REAL*4 VARIABLE
TMP	000160	REAL*4 VARIABLE
PRS	000166	REAL*4 VARIABLE
I	010210	INTEGER*2 VARIABLE
LOCK	000000	INTEGER*2 PROCEDURE
IQSET	000000	INTEGER*2 PROCEDURE
IGOOF	010212	INTEGER*2 VARIABLE
IRAD50	000000	INTEGER*2 PROCEDURE
ICHAN	010214	INTEGER*2 VARIABLE
IGETC	000000	INTEGER*2 PROCEDURE
IENTER	000000	INTEGER*2 PROCEDURE
IDAY	010216	INTEGER*2 VARIABLE
IYR	010220	INTEGER*2 VARIABLE
INUM	010222	INTEGER*2 VARIABLE
IWRITW	000000	INTEGER*2 PROCEDURE
UNLOCK	000000	REAL*4 PROCEDURE
CLOSEC	000000	REAL*4 PROCEDURE

COMMON BLOCK /X/ LENGTH 000060

MON 000000 INTEGER*2 ARRAY (24)

TECHNOLOGY INCORPORATED
LIFE SCIENCE DIVISION

PROGRAM NAME:..... STRESS
AUTHOR:..... ROY A. REED
..... CHUCK MANN
..... WILLIAM G. CROSIER
DATE:..... 2/NOV/78

COMPUTER SYSTEM:..... DEC PDP 1103 "LSI-11"
OPERATING SYSTEM:..... RT-11 VOIC-03A

COMPILING SEQUENCE:

.R FORTRA <CR>
*STRESS=STRESS,INITI,DATA/U <CR>
*CAL=CAL <CR>
*PFT=PFT <CR>
*EX=EX,RPTHOR,WRKADJ,BREATH,HPT30,TIME <CR>
*LIGHT=LIGHT <CR>
*<CR>
.R MACRO <CR>
*ADISR=ADISR <CR>
*TRSEUF=TRSEUF <CR>
*SNATCH=SNATCH<CR>
* <CR>

RUN MODUAL LINKING SEQUENCE:

.R LINK <CR>
*STRESS=STRESS,CAL,PFT,EX,CX0:SYSLIB/F/C <CR>
*LIGHT,ADISR,TRSEUF,SNATCH <CR>
*<CR>

CALLING SEQUENCE:

.R STRESS <CR>

PURPOSE:

A MICROCOMPUTER - BASED DATA ACQUISITION AND CONTROL
SYSTEM FOR CARDIOPULMONARY STRESS TESTING.

C PROGRAM: PROCESS, VERSION 2.0 29 AUG. 1979
C FUNCTION: RUN CARDIOPULMONARY STRESS TEST
C

0001 INTEGER MAJREV,MINREV
0002 INTEGER DROUTB,DRINB,TESTYP
0003 INTEGER*4 IN,OUT
0004 COMMON /PDATA/ICATA(1536)
0005 DATA IN/'IN '//,OUT/'OUT '//
0006 COMMON/X/MO(24)
0007 DATA MO/'JA','N ','FE','B ','MA','R ','AP','R ',
@'MA','Y ','JU','N ','JU','L ','AU','G ','SE','P ',
@'OC','T ','NO','V ','DE','C '

0008 EQUIVALENCE (ICATA(13),MAJREV),(ICATA(14),MINREV)

C-- SETS THE MAJOR AND MINOR REVISION NUMBERS
C FOR THIS VERSION OF CVDAS

0009 MAJREV=2
0010 MINREV=0

C
C--DRV11 PARALLEL INTERFACE ADDRESSES:

0011 DROUTB= "167772 ! " OUTPUT BUFFER (KATHY'S DRAI,DRRI)
0012 DRINB= "167774 ! " INPUT BUFFER (KATHY'S DRAO,DRBO)

C THE LOW ORDER TWO BITS OF DRCSR WILL BE DINGLED
C BY PROCES FOR CALIBRATION
C

C
C--LOCAL PARAMETER INITIALIZATION:

C SET ALL PUSHBUTTON BIT PATTERNS
0013 ICALPB="1
0014 IFVCPB="100
0015IRSTPB="2
0016 IENDT="40

C
C--FOLLOWING SUBROUTINE DATA INPUTS SUBJECT DATA AND CERTAIN OTHER INFO.

0017 CALL DATA(IN)
C SWITCH ANALOG INTERFACE TO APPROPRIATE TEST TYPE
0018 TESTYP=1
0019 IF(ICATA(12).EQ.2)TESTYP=0
0021 IF(ICATA(12).EQ.3)TESTYP=2
0023 TESTYP=256*TESTYP
0024 CALL IPOKE(DROUTB,TESTYP)
0025 IF(ICATA(43).GT.0)GO TO 1000
C TURN ON CAL LIGHT
0027 CALL LIGHT("41)
0028 200 IRES=IPEEK(DRINB)
0029 IF(IENDT.EQ.IRES)GO TO 2000
0031 IF(ICALPB .NE. IRES) GO TO 200
0033 300 CALL CAL
D TYPE 1
D1 FORMAT(' AFTER CAL')
0034 IDATA(43)=IDATA(43)+1
0035 CALL DATA(OUT)
0036 1000 CONTINUE

```
D      TYPE 2
D2     FORMAT(' AFTER 1000')
C
C--TURN ON APPROPRIATE LIGHTS.
C
0037   CALL LIGHT("143)
0038   IRES=0
0039   IRES=IPEEK(DRINE)
0040   IF(IRES.EQ.IENDT)GO TO 2000
C
C--CHECK WHICH ONE IF SO AND DO WHAT IS REQUESTED:
C
0042   IF(ICALPB .EQ.IRES) GO TO 300   !REPEAT CALIBRATION
0044   IF(IFVCPB .EQ.IRES) CALL PFT     ! FORCED VITAL CAPACITY TEST
0046   IF(IRSTPB .EQ.IRES)CALL EX      !EXERCISE TEST
D      TYPE 3
D3     FORMAT(' AFTER EX')
0048   IF(IRSTPB .NE. IRES) GO TO 1000
0050 2000 CONTINUE
D      TYPE 4
D4     FORMAT(' AFTER 2000')
0051   CALL IPOKE(DROUTB,0)
0052   CALL DATA(OUT)
0053   STOP
0054   END
```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
------	--------	------------

DROUTB	000024	INTEGER*2 VARIABLE
DRINB	000026	INTEGER*2 VARIABLE
TESTYP	000030	INTEGER*2 VARIABLE
IN	000006	INTEGER*4 VARIABLE
OUT	000012	INTEGER*4 VARIABLE
ICALPB	000032	INTEGER*2 VARIABLE
IFVCPB	000034	INTEGER*2 VARIABLE
IRSTPB	000036	INTEGER*2 VARIABLE
IENDT	000040	INTEGER*2 VARIABLE
DATA	000000	REAL*4 PROCEDURE
IPOKE	000000	INTEGER*2 PROCEDURE
LIGHT	000000	INTEGER*2 PROCEDURE
IRES	000042	INTEGER*2 VARIABLE
IPEEK	000000	INTEGER*2 PROCEDURE
CAL	000000	REAL*4 PROCEDURE
PFT	000000	REAL*4 PROCEDURE
EX	000000	REAL*4 PROCEDURE

COMMON BLOCK /PDATA/ LENGTH 006000

DATA	000000	INTEGER*2 ARRAY (1536)
MAJREV	000030	INTEGER*2 VARIABLE
MINREV	000032	INTEGER*2 VARIABLE

COMMON BLOCK /X/ LENGTH 000060

MD	000000	INTEGER*2 ARRAY (24)
----	--------	----------------------

C-- SUBROUTINE INITI (VERSION 2.0), USED TO SET THE TIMERS, COUNTERS,
C-- & FLAGS USED IN THE ADISR INTERRUPT ROUTINE.
C--

```
0001 SUBROUTINE INITI
0002 INTEGER BPCNT,BPCNT1,BPTIME
0003 INTEGER CWLCNT,AHRCNT,AHR,CWL,CSPD,CSBP,CDBP
0004 INTEGER AHR1,CWL1,CSPD1,CSBP1,CDBP1
0005 INTEGER COUNT,CNT10,CNT30,SECS,SECS30,GASDSP
0006 INTEGER FLRCNT,FLRFLG,EXHCNT,EXHALE
0007 COMMON /ADCOM1/IDUM1(34)
0008 COMMON/TIMCNT/COUNT,CNT10,CNT30,FLRCNT,FLRFLG,EXHCNT
0009 EQUIVALENCE (IDUM1(1),GASCSP),(IDUM1(10),AHRCNT)
0010 EQUIVALENCE (IDUM1(2),EXHALE)
0011 EQUIVALENCE (IDUM1(14),AHR),(IDUM1(15),AHR1)
0012 EQUIVALENCE (IDUM1(16),CWL),(IDUM1(17),CWL1)
0013 EQUIVALENCE (IDUM1(18),CWLCNT),(IDUM1(28),SECS)
0014 EQUIVALENCE (IDUM1(19),BPTIME),(IDUM1(22),BPCNT1)
0015 EQUIVALENCE (IDUM1(20),CSPD),(IDUM1(21),CSPD1)
0016 EQUIVALENCE (IDUM1(23),CSBP),(IDUM1(24),CSBP1)
0017 EQUIVALENCE (IDUM1(25),CDBP),(IDUM1(26),CDBP1)
0018 EQUIVALENCE (IDUM1(27),BPCNT)
0019 EQUIVALENCE (IDUM1(29),SECS30),(IDUM1(32),AHL)
0020 EQUIVALENCE (IDUM1(33),AHL1),(IDUM1(34),AWLCNT)
0021 EXHALE=0
0022 BPTIME=550
0023 BPCNT1=0
0024 BPCNT=10
0025 AHR=0
0026 AHR1=0
0027 CWL=0
0028 CWL1=0
0029 CSPD=0
0030 CSPD1=0
0031 CSBP=0
0032 CSBP1=0
0033 CDBP=0
0034 CDBP1=0
0035 COUNT=10
0036 CNT10=10
0037 CNT30=30
0038 EXHCNT=0
0039 CWLCNT=0
0040 AHRCNT=0
0041 SECS=0
0042 SECS30=0
0043 GASDSP=0
0044 AHL=0
0045 AHL1=0
0046 AWLCNT=0
0047 RETURN
0048 END
```

FORTAN IV STORAGE MAP

NAME OFFSET ATTRIBUTES

COMMON BLOCK /ADCOM1/ LENGTH 000106

IDUM1	000000	INTEGER*2	ARRAY (34)
GASDSP	000000	INTEGER*2	VARIABLE
AHRCNT	000022	INTEGER*2	VARIABLE
EXHALE	000002	INTEGER*2	VARIABLE
AHR	000032	INTEGER*2	VARIABLE
AHR1	000034	INTEGER*2	VARIABLE
CWL	000036	INTEGER*2	VARIABLE
CWL1	000040	INTEGER*2	VARIABLE
CWLCNT	000042	INTEGER*2	VARIABLE
SECS	000066	INTEGER*2	VARIABLE
BPTIME	000044	INTEGER*2	VARIABLE
BPCNT1	000052	INTEGER*2	VARIABLE
CSPD	000046	INTEGER*2	VARIABLE
CSPD1	000050	INTEGER*2	VARIABLE
CSBP	000054	INTEGER*2	VARIABLE
CSBP1	000056	INTEGER*2	VARIABLE
CDBP	000060	INTEGER*2	VARIABLE
CDBP1	000062	INTEGER*2	VARIABLE
BPCNT	000064	INTEGER*2	VARIABLE
SECS30	000070	INTEGER*2	VARIABLE
AWL	000076	REAL*4	VARIABLE
AWL1	000100	REAL*4	VARIABLE
AWLCNT	000102	REAL*4	VARIABLE

COMMON BLOCK /TIMCNT/ LENGTH 000014

COUNT	000000	INTEGER*2	VARIABLE
CNT10	000002	INTEGER*2	VARIABLE
CNT30	000004	INTEGER*2	VARIABLE
FLRCNT	000006	INTEGER*2	VARIABLE
FLRFLG	000010	INTEGER*2	VARIABLE
EXHCNT	000012	INTEGER*2	VARIABLE

```
0001      SUBROUTINE DATA(ISTAT)
0002      LOGICAL*1 IFILE(12)
0003      INTEGER FILE(4)
      C
0004      INTEGER*4 IOSTAT,CUT
      C
      C
0005      COMMON /PCDATA/ IPATNT(1536) ! THIS IS PATIENT DATA.
      C
      C
0006      DATA OUT/'OUT '/
0007      DATA IFILE/'D','X','1',' ',' ',' ',' ',' ',' ',' ','D','A','T'/
      C
0008      IF(ISTAT(2).NE.0) STOP 'QUEUE ALLOCATION FAILURE'
0010      IF(ISTAT.EQ.CUT)GO TO 10
0012      TYPE 1025
0013 1025  FORMAT('S  PLEASE TYPE PATIENT DATA FILE NAME (UP TO 6 SYMBOLS):'
0014          ACCEPT 1030,(IFILE(I),I=4,9)
0015 1030  FORMAT(6A1)
      C
      C
      C
0016 10  CONTINUE
      C
      C--CONVERT ASCII FILE NAME TO RADIX 50
      C
0017      IF(ISTAT.NE.CUT) CALL IRAC50(12,IFILE,FILE)
      C
      C--OPEN FILE:
      C
0019      ICHAN =IGETC()
0020      IF (ICHAN .LT. 0) STOP 'NO AVAILABLE CHANNEL'
0022      IF (LOOKUP(ICHAN,FILE) .LT. 0) STOP 'FILE ALLOCATION FAILURE'
      C
      C--NOW READS OR WRITES DATA:
      C
0024      IF(ISTAT.EQ.CUT) GO TO 1050
0026      IF(IREADW(1536,IPATNT,0,ICHAN).LT.0) STOP 'DISK READ ERROR'
0028      GO TO 1060
0029 1050  CONTINUE
0030      IF(IWRITE(1536,IPATNT,0,ICHAN).LT.0) STOP 'DISK WRITE ERROR'
0032 1060  CONTINUE
0033      CALL CLOSEC(ICHAN)
0034      RETURN
0035      END
```


FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
IFILE	000016	LOGICAL*1 ARRAY (12)
FILE	000032	INTEGER*2 ARRAY (4)
IOSTAT	000014	INTEGER*4 PARAMETER VARIABLE
OUT	000042	INTEGER*4 VARIABLE
IQSET	000000	INTEGER*2 PROCEDURE
I	000330	INTEGER*2 VARIABLE
IRAD50	000000	INTEGER*2 PROCEDURE
ICHAN	000332	INTEGER*2 VARIABLE
IGETC	000000	INTEGER*2 PROCEDURE
LOOKUP	000000	INTEGER*2 PROCEDURE
IHEADW	000000	INTEGER*2 PROCEDURE
IWRITW	000000	INTEGER*2 PROCEDURE
CLOSEC	000000	REAL*4 PROCEDURE

COMMON BLOCK /PDATA/ LENGTH 006000

IPATNT 000000 INTEGER*2 ARRAY (1536)

```

C-- SUBROUTINE CALIBRATION, VERSION 2.0
C-- WRITTEN BY: ROY A. REED
C-- DATE OCT 13, 1978
C
C   REWRITTEN BY: WILLIAM G. CROSIER
C   DATE: AUG 24, 1979
C   MODIFIED BY: ABHIJIT GADGIL
C   DATE: DEC, 1979
C--
0001 SUBROUTINE CAL
C--
C--
C-- SET UP VARIABLES
C--
0002 LOGICAL*1 IANS, YES, NO
0003 INTEGER SEX, SSNO(3), UNIGNO(2), RETNO
0004 INTEGER CPID, DATE(3), NAME(14)
0005 INTEGER GASBUF, ADCTBF, IDATA
0006 INTEGER GASDSP, EXHALE, SVTHRS, FVC, PFT
0007 INTEGER HRTHRS, CWLCNT, AHRCNT
0008 INTEGER BPTIME, EPCNT, BPCNT1
0009 INTEGER LFLSHR, FILE(4)
0010 INTEGER PCSR, PCCB, PCS, PCC
0011 INTEGER DROUTB, CRINE, DRSH
0012 INTEGER*4 CSPD, AHR, CWL, CSBP, CDBP
0013 REAL SLPS(9), SLFST(2), SEP
0014 REAL AMBTMP, ANBFFS
C   REAL AO2, AN2, AC02
C   REAL HCSV, LCSV, FCCX, LCOX, HCN2, LCN2, HCCO2, LCCO2
C   REAL HCHR, LCHR, HCWL, LCWL, HCAS, LCAS, HCSBP, LCSBP, HCCBP, LCDBP
0015 REAL RDATA(24), CALBUF(22)
C   REAL LSVOL, LO2, LN2, LC02
C   REAL HSVOL, HO2, HN2, HC02
C   REAL LHR, LWL, LAS, LSBP, LDBP
C   REAL MHR, HWL, HAS, HSBP, HDBP
0016 REAL REFVAL(2,9), CALDAT(2,9), TEMP, CNT
0017 REAL*8 LHA, LOW, -IGH, AMB, CHAN(9)
C--
C-- SET UP COMMON AREAS
C--
0018 COMMON /ADCOM/JCHAN(9)
0019 COMMON /PDATA/ICATA(1536)
0020 COMMON /GASCOM/GASBUF(2000)
0021 COMMON /ADCOM1/ADCTBF(34)
0022 COMMON /X/MD(24)
C--
C-- SET UP EQUIVALENCE STATEMENTS
C--
0023 EQUIVALENCE (ADCTBF(1), GASDSP), (ADCTBF(2), EXHALE)
%, (ADCTBF(3), SVTHRS), (ADCTBF(4), FVC)
%, (ADCTBF(8), PFT), (ADCTBF(10), AHRCNT), (ADCTBF(11), HRTHRS)
%, (ADCTBF(14), AHR)
%, (ADCTBF(16), CWL)
%, (ADCTBF(18), CWLCNT), (ADCTBF(19), BPTIME)

```

```

%,(ADCTBF(20),CSPD)
%,(ADCTBF(22),BPCNT1),(ADCTBF(23),CSBP)
%,(ADCTBF(25),CDSP)
%,(ADCTBF(27),BPCNT)
%,(ADCTBF(30),LFLSHR)
0024      EQUIVALENCE (IDATA(2),SEX),(IDATA(3),SSNO(1))
%,(IDATA(7),UNIGNO(1)),(IDATA(9),RETNQ)
%,(IDATA(10),CPID)
%,(IDATA(26),NAME(1))
%,(IDATA(40),DATE(1))
0025      EQUIVALENCE (IDATA(73),RCATA(1))
0026      EQUIVALENCE (IDATA(54),AMBTMP),(IDATA(57),AMBPRS)
C          EQUIVALENCE(CALBUF(1),HCSV),(CALBUF(2),LCSV)
C          %,(CALBUF(3),HCCX),(CALBUF(4),LCOX)
C          %,(CALBUF(5),HCN2),(CALBUF(6),LCN2)
C          %,(CALBUF(7),HCCO2),(CALBUF(8),LCCO2)
C          %,(CALBUF(9),HCHR),(CALBUF(10),LCHR)
C          %,(CALBUF(19),HCSBP),(CALBUF(20),LCSBP)
C          %,(CALBUF(21),HCSBP),(CALBUF(22),LCSBP)
0027      EQUIVALENCE (CALBUF,REFVAL)
C--
C--      SET UP DATA VALUES
C--
0028      DATA SLPS/.0025,.0300,.0245,.0025,.0720,.0757,.0249,.0623,.0624/
0029      DATA SLPST/0.01,0.0025/
0030      DATA SEP/.02/
0031      DATA YES/1HY/,NO/1HN/
0032      DATA DRQUT8/"167772/,DRINE/"167774/,DRSR/"167770/
0033      DATA PCSR/"172540/,PCCB/"172542/
0034      DATA PCS/"113/,PCC/"144/
C          DATA MD/"JA",N,"FE",B,"MA",R,"AP",R,
C          @ "MA",Y,"JU",N,"JU",L,"AU",G,"SE",P,
C          @ "OC",T,"NO",V,"DE",C /
0035      DATA LOW,HIGH,AMB/" LOW ",HIGH " ,AMBIENT "/
0036      DATA CHAN/"SP. VOL.",OXYGEN " ,NITROGEN",CO2 " ,
@ "HT.RATE",WL/ELEV.",SPEED " ,SYS. BP",DIAS. BP"/
0037      IF(IDATA(12).EQ.2) SLPS(6)=SLPST(1)
0039      IF(IDATA(12).EQ.2) SLPS(7)=SLPST(2)
C--
C--      READ IN THE CALIBRATION FILE, WHICH CONTAINS THE STANDARD DATA
C--      USED IN CALIBRATION OF SLOPES AND INTERCEPTS.
C--
0041      CALL IRADSO(12,'DXOCALCONCAT',FILE)
0042      ICHAN=IGETC()                                !GET A I/O CHANNEL
0043      IF(ICCHAN.LT.0)STOP 'CHANNEL ALLOCATION FAILURE IN CAL'
0045      IF(LOOKUP(ICCHAN,FILE).LT.0)STOP 'FILE ALLOCATION FAILURE IN CAL'
0047      IF(IREADW(44,CALBUF,1,ICCHAN).LT.0) STOP 'DISK READ ERROR IN CAL'
0049      IF(IDATA(12).EQ.2)GO TO 5
0051      REFVAL(1,7)=CALBUF(15)                        !BIKE SPEED          HIGH REF. VALI
0052      REFVAL(2,7)=CALBUF(16)                        !                      LOW REF. VALUE
0053      GO TO 6
0054      5      CONTINUE
0055      REFVAL(1,6)=CALBUF(13)                        !TREADMILL ELEVATION    HIGH
0056      REFVAL(2,6)=CALBUF(14)                        !                      LOW

```

```

0057      REFVAL(1,7)=CALPUF(17)      !TREADMILL SPEED      HIGH
0058      REFVAL(2,7)=CALBUF(18)      !                      LOW
0059      6      CONTINUE
0060      REFVAL(1,8)=CALBUF(19)      !SBP      HIGH
0061      REFVAL(2,8)=CALBUF(20)      !                      LOW
0062      REFVAL(1,9)=CALBUF(21)      !DBP      HIGH
0063      REFVAL(2,9)=CALBUF(22)      !                      LOW
0064      7      CONTINUE
0065      IF(IDATA(11).EQ.1) GO TO 10
      C-      INITIALIZE PROGRAMMABLE GAIN CODE BITS IN A/D COMMAND
      C-      WORD ARRAY (LATER USED BY INTERRUPT SERVICE ROUTINE)
0067      IAD="10
0068      DO 9 ICH=1,9
0069      JCHAN(ICH)=IAD
0070      9      IAD=IAD+"400
0071      10     CONTINUE
      C--
      C      INITIALIZE PARAMETERS
      C--
0072      PFT=0
0073      FVC=0
0074      HRTHRS=-2048
0075      SVTHRS=-2048
0076      LFLSHR="1      !SET CAL LIGHT TO FLASH
      C
      C      COLLECT LOW CAL DATA
      C
0077      MODE=1      !      SET LOW CAL MODE
0078      ISUB=2      !      USE SECOND ROW OF CALDAT ARR
0079      DELAY=5.0      !      SET DELAY AT 5 SECONDS
0080      LHA=LOW      !      SET ALPHA IDENTIFIER FOR CAL
0081      ASSIGN 162 TO IRETRN      !      ASSIGN RETURN LABEL
0082      GO TO 5000      !      GO TO CAL DATA COLLECTION/AVERAGING
      C
      C      COLLECT HIGH CAL DATA
      C
0083      162     MODE=2      !      SET HIGH CAL MODE
0084      ISUB=1      !      USE FIRST ROW OF CALDAT ARRAY
0085      DELAY=5.0      !      SET DELAY AT 5 SECONDS
0086      LHA=HIGH      !      SET ALPHA IDENTIFIER FOR CAL MODE
0087      ASSIGN 262 TO IRETRN      !      ASSIGN RETURN LABEL
0088      GO TO 5000      !      GO TO DATA COLLECTION/AVERAGING ROUTINE
      C--
      C-- CALCULATE THE SLOPES AND Y INTERCEPTS FOR ALL 9 CHANNELS
0089      262     K=7
0090      DO 265 ICHAN=1,9
0091      RDATA(K)=0.0
0092      DENOM=CALDAT(1,ICHAN)-CALDAT(2,ICHAN)
0093      IF(DENOM.NE.0.0)RDATA(K)=(REFVAL(1,ICHAN)-REFVAL(2,ICHAN))/DENOM
0095      K=K+1
0096      RDATA(K)=REFVAL(2,ICHAN)-RDATA(K-1)*CALDAT(2,ICHAN) !Y-INT.
0097      265     K=K+1
      C
      C      COLLECT AMBIENT AIR DATA
      C

```

```

0098      MODE=3              !      SET AMBIENT AIR MODE
0099      ISUB=1              !      USE FIRST ROW OF CALDAT ARRAY
0100      DELAY=20.0          !      SET DELAY 20 SECONDS
0101      LHA=AMB             !      SET ALPHA IDENTIFIER FOR CAL MODE
0102      ASSIGN 370 TO IRETRN !      ASSIGN RETURN LABEL
0103      GO TO 5000          !      GO TO DATA COLLECTION/AVERAGING ROUTINE

C--
C-- PH20=VAPOR PRESSURE OF WATER AT AMBIENT
C--      PRESSURE AND TEMPERATURE.
C--
0104 370 PH20=.0369*AMBTMP*AMBTMP-.4012*AMBTMP+10.76
C--
C-- STPD=FACTOR TO ADJUST VOLUME AT AMBIENT CONDITIONS
C--      TO VOLUME AT STANDARD TEMPERATURE AND PRESSURE (DRY).
C--      NORMAL BODY TEMPERATURE IS 310 DEGREES KELVIN
C--
0105      BTPS=310.0/(AMBTMP+273.0)
0106      STPD=(273.16/(AMBTMP+273.16))*((AMBPRS-PH20)/760.)

C--
C-- CALCULATE THE PERCENTS
C--
0107      RDATA(3)=CALDAT(1,2)*RDATA(9)+RDATA(10)      !AMBIENT 02
0108      RDATA(4)=CALDAT(1,3)*RDATA(11)+RDATA(12)     !AMBIENT N2
0109      RDATA(5)=CALDAT(1,4)*RDATA(13)+RDATA(14)     !AMBIENT CO2
0110      RDATA(1)=STPD                                  !STPD FACTOR
0111      RDATA(2)=BTPS                                  !BTPS FACTOR

C--
C--
C-- SWITCH ANALOG INTERFACE TO NORMAL MODE
C--
0112      CALL IPOKE(DRSR,0)

C--
C-- WRITE REPORT
C--
0113      IMM=2*DATE(1)
0114      IMM1=IMM-1
0115 999 TYPE 1000
0116 1000 FORMAT(/,32X,'CALIBRATION REPORT',/36X,'CDAS - V2.0')
0117      TYPE 1010,CPID,(UNIGNO(I),I=1,2),RETNO,DATE(2),MD(IMM1),
      @MD(IMM),DATE(3)
0118      TYPE 1020,(NAME(I),I=1,14),SEX,(SSNO(I),I=1,3)
0119 1010 FORMAT(/,5X,'CPID NO.: ',15,28X,'UNIQUE NO.: 'A2,I4,/,
      X5X,'RETEST NO.: ',12,28X,'DATE: ',12,1X,2A2,I2,/)
0120 1020 FORMAT(5X,'SUBJECT NAME: '14A2,1X,'SEX: ',A2,/,
      X5X,'SOC. SEC. NO.: ',13,'-',12,'-',14)
0121      TYPE 1030
0122 1030 FORMAT(/,5X,'AMBIENT CONDITIONS:',23X,'CALIBRATION FACTORS:'
      X,/,40X,'CHANNEL SLCPE Y-INT.')
0123      TYPE 1040,AMBTMP,RDATA(7),RDATA(8)
      X,AMBPRS,RDATA(9),RDATA(10)
      X,(RDATA((I-9)/2),RDATA(I),RDATA(I+1),I=11,19,2)
      X,(RDATA(I),I=21,24)
0124 1040 FORMAT(/,5X,'TEMP.: ',F4.1,' C',16X,'SP.VOL',5X,F9.5,3X,F8.3,
      X/,5X,'PRES.: ',F5.1,' MM HG',12X,'02 ',F9.5,3X,

```

```

XF0.3,/,5X,'STPD FACTOR: ',F5.3,16X,'N2',F9.5,3X,F8.3
X/,5X,'BTP9 FACTOR: ',F5.3,16X,'CO2',F9.5,3X,F8.3
X/,5X,'O2: ',F4.1,'%',17X,'HR',F9.5,3X,F8.3
X/,5X,'N2: ',F4.1,'%',17X,'ELEV/WL',F9.5,3X,F8.3
X/,5X,'CO2: ',F4.2,'%',17X,'SPEED',F9.5,3X,F8.3
X/,40X,'SBP',F9.5,3X,F8.3
X/,40X,'DBP',F9.5,3X,F8.3

0125 IRA=1
0126 IF(IDATA(11).EQ.2)IRA=1.414
0128 LINE=29
0129 DO 3000 IC=1,9
0130 IS=(IC-1)*2+7
0131 SMAX=SLPS(IC)*IRA+(SLPS(IC)*IRA*SEP)
0132 SMIN=SLPS(IC)*IRA-(SLPS(IC)*IRA*SEP)
0133 IF(RDATA(IS).LT.SMAX.AND.RDATA(IS).GT.SMIN) GO TO 3000
0135 PERER=((RDATA(IS)-SLPS(IC))/SLPS(IC))*100.0
0136 LINE=LINE+2
0137 TYPE2010, CHAN(IC),PERER
0138 2010 FORMAT(/1X,A8,' CALIB. OUT OF RANGE BY ',F6.1,'%')
0139 3000 CONTINUE
0140 LINE=66-LINE
0141 DO 3100 IL=1,LINE
0142 3100 TYPE 3150
0143 3150 FORMAT(1X)
C
C SET TIME OF DAY TIMER IF SYSTEM CLOCK SHOWS LESS THAN
C 5 HOURS (18000 SECS) SINCE MIDNIGHT OR SYSTEM BOOT, OR
C IF WE ARE RUNNING PLAYBACK MODE (SINCE SYSTEM CLOCK TIME
C WILL NOT BE TEST TIME).
C
0144 IF (SECNDS(0.).LT.18000..OR.IDATA(11).EQ.2) CALL TIMSET
C--
C-- SET LIGHTS ON FRONT PANNEL
C--
0146 CALL LIGHT("1")
C--
C-- ALL DONE RETURN TO PROCES
C--
0147 RETURN
C
C ROUTINE TO COLLECT & AVERAGE DATA FROM A/D CONVERTER
C FOR CALIBRATION
C
C INITIALIZE COUNTERS, TIMERS, & FLAGS
0148 5000 CALL INITI
0149 BPTIME=0
0150 BPCNT=40
0151 BPCNT1=40
C IS TEST BEING RUN IN REAL TIME(LIVE) ?
0152 IF (IDATA(11).EQ.1) GO TO 5010
0154 5006 TYPE 5007, LHA
0155 5007 FORMAT (/ ' IS TAPE READY FOR ',A7,' CALIB. (Y OR N) ?',$(
0156 ACCEPT 5008,IANS
0157 5008 FORMAT(A1)

```

```

0158      IF (IANS,NE,YES) GO TO 5006
0160 5010  CONTINUE
0161      CALL LIGHT ("1)                ! TURN ON CAL LIGHT
0162      CALL IPOKE(DRSR,MODE)          ! SEND CAL SIGNAL TO LAB EQUIP.
0163      TI=SECNDS(0.0)
0164 5015  DELTA=SECNDS(TI)
0165      IF (DELTA,LT,DELAY) GO TO 5015 ! WAIT FOR SIGNALS TO STABILIZE
0166      CALL IPOKE (PCCB,PCC)           ! SET PROGRAMMABLE CLOCK FOR 100
0167      CALL IPOKE (PCSR,PCS)           ! START CLOCK FOR DATA ACQUISITI
0168 5020  IF (BPCNT1,GT,0) GO TO 5020 ! WAIT UNTIL 4 SEC. OF DATA IS COLLECTED
0169      CALL IPOKE (PCSR,0)             ! STOP CLOCK
0171      CALL LIGHT ("1)                ! TURN ON CAL LIGHT
0172      C DO 5051 I=1, GASDSP,8
      C TYPE 5050, I, (GASBUF(J),J=I,I+7) ! TYPE GAS BUFFER CONTENTS
0173 5050  FORMAT (/1X,I4,' '),8(1X,I6))
0174 5051  CONTINUE
      C
      C CALCULATE AVG. VALUE FOR EACH OF THE 9 CHANNELS
      C
0175      DO 5080 K=1,9
0176 5080  CALDAT(ISUB,K)=0.0              ! INITIALIZE VALUES FOR SUMMING
0177      SAMP=FLOAT(GASDSP)/4.0
0178      DO 6010 K=1,4
0179      DO 6000 IAV=1, GASDSP,4          ! CHAN. 1-4 ARE GAS CHANNELS
0180 6000  CALDAT(ISUB,K)=CALDAT(ISUB,K) + FLOAT(GASBUF(IAV+K-1))
0181 6010  CALDAT(ISUB,K)=CALDAT(ISUB,K)/SAMP
      D TYPE 6050,LHA,GASDSP,SAMP,(CALDAT(ISUB,K),K=1,4)
0182 6050  FORMAT(/1X,A8,' CAL VALUES'/' GASDSP      SAMP      SP VOL',6X,'02',
      @ 8X,'N2',7X,'CC2'/I10,5F10.2)
      C
      C CALCULATE AVERAGES FOR CHAN. 5-9
      C
0183      IF(IAJFLT(AHR,TEMP).EQ.-2) TYPE 9999
0185      IF(AHRCNT,NE.0) CALDAT(ISUB,5)=TEMP/FLOAT(AHRCNT)
0187      IF(IAJFLT(CWL,TEMP).EQ.-2) TYPE 9999
0189      IF(CWLCNT,NE.0) CALDAT(ISUB,6)=TEMP/FLOAT(CWLCNT)
0191      IF(IAJFLT(CSPD,TEMP).EQ.-2) TYPE 9999
0193      IF(CWLCNT,NE.0) CALDAT (ISUB,7)=TEMP/FLOAT(CWLCNT)
0195      CNT=FLOAT(BPCNT-BPCNT1)
0196      IF(CNT,EQ.0.0) GO TO 6055
0198      IF(IAJFLT(CSBP,TEMP).EQ.-2) TYPE 9999
0200      CALDAT(ISUB,8)=TEMP/CNT
0201      IF(IAJFLT(CDBP,TEMP).EQ.-2) TYPE 9999
0203      CALDAT(ISUB,9)=TEMP/CNT
0204 6055  CONTINUE
      D TYPE 6060,AHRCNT,CWLCNT,BPCNT,BPCNT1,CNT,(CALDAT(ISUB,K),K=5,9)
0205 6060  FORMAT (/ ' AHRCNT CWLCNT BPCNT BPCNT1      CNT'/4I7,F7.2/
      @ 5X,'HR',8X,'WL',6X,'SPEED',6X,'SBP',7X,'DBP'/5F10.2)
      C
0206      GO TO IRETRN                  ! RETURN
      C
0207 9999  FORMAT(/,1X,'OVERFLOW IN INTEGER*4 TO REAL CONVERSION')
0208      END

```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
FILE	000014	INTEGER*2 ARRAY (4)
SLPS	000024	REAL*4 ARRAY (9)
SLPST	000070	REAL*4 ARRAY (2)
CALBUF	000100	REAL*4 ARRAY (22)
REFVAL	000100	REAL*4 ARRAY (2,9) VECTORED
CALDAT	000230	REAL*4 ARRAY (2,9) VECTORED
CHAN	000340	REAL*8 ARRAY (9)
IAN5	002676	LOGICAL*1 VARIABLE
YES	000454	LOGICAL*1 VARIABLE
NO	000455	LOGICAL*1 VARIABLE
PCSR	000464	INTEGER*2 VARIABLE
PCCB	000466	INTEGER*2 VARIABLE
PCS	000470	INTEGER*2 VARIABLE
PCC	000472	INTEGER*2 VARIABLE
DROUTB	000456	INTEGER*2 VARIABLE
DRINB	000460	INTEGER*2 VARIABLE
DRSR	000462	INTEGER*2 VARIABLE
SEP	000450	REAL*4 VARIABLE
TEMP	002700	REAL*4 VARIABLE
CNT	002704	REAL*4 VARIABLE
LHA	002710	REAL*8 VARIABLE
LOW	000474	REAL*8 VARIABLE
HIGH	000504	REAL*8 VARIABLE
AMB	000514	REAL*8 VARIABLE
IRAD50	000000	INTEGER*2 PROCEDURE
ICHAN	002720	INTEGER*2 VARIABLE
IGETC	000000	INTEGER*2 PROCEDURE
LOOKUP	000000	INTEGER*2 PROCEDURE
IREADW	000000	INTEGER*2 PROCEDURE
IAD	002722	INTEGER*2 VARIABLE
ICH	002724	INTEGER*2 VARIABLE
MODE	002726	INTEGER*2 VARIABLE
ISUB	002730	INTEGER*2 VARIABLE
DELAY	002732	REAL*4 VARIABLE
IRETRN	002736	INTEGER*2 VARIABLE
K	002740	INTEGER*2 VARIABLE
DENOM	002742	REAL*4 VARIABLE
PH20	002746	REAL*4 VARIABLE
BTPS	002752	REAL*4 VARIABLE
STPD	002756	REAL*4 VARIABLE
IPOKE	000000	INTEGER*2 PROCEDURE
IMM	002762	INTEGER*2 VARIABLE
IMM1	002764	INTEGER*2 VARIABLE
I	002766	INTEGER*2 VARIABLE
IRA	002770	INTEGER*2 VARIABLE
LINE	002772	INTEGER*2 VARIABLE
IC	002774	INTEGER*2 VARIABLE
IS	002776	INTEGER*2 VARIABLE
SMAx	003000	REAL*4 VARIABLE
SMin	003004	REAL*4 VARIABLE
PERER	003010	REAL*4 VARIABLE
IL	003014	INTEGER*2 VARIABLE

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
------	--------	------------

SECNDS	000000	REAL*4 PROCEDURE
TIMSET	000000	REAL*4 PROCEDURE
LIGHT	000000	INTEGER*2 PROCEDURE
INITI	000000	INTEGER*2 PROCEDURE
TI	003016	REAL*4 VARIABLE
DELTA	003022	REAL*4 VARIABLE
SAMP	003026	REAL*4 VARIABLE
FLOAT	000000	REAL*4 PROCEDURE
IAV	003032	INTEGER*2 VARIABLE
IAJFLT	000000	INTEGER*2 PROCEDURE

COMMON BLOCK /ADCOM/ LENGTH 000022

JCHAN 000000 INTEGER*2 ARRAY (9)

COMMON BLOCK /PDATA/ LENGTH 006000

IDATA	000000	INTEGER*2 ARRAY (1536)
SEX	000002	INTEGER*2 VARIABLE
SSNO	000004	INTEGER*2 ARRAY (3)
UNIGNO	000014	INTEGER*2 ARRAY (2)
RETNO	000020	INTEGER*2 VARIABLE
CPID	000022	INTEGER*2 VARIABLE
NAME	000062	INTEGER*2 ARRAY (14)
DATE	000116	INTEGER*2 ARRAY (3)
RDATA	000220	REAL*4 ARRAY (24)
AMBTMP	000152	REAL*4 VARIABLE
AMBPRS	000160	REAL*4 VARIABLE

COMMON BLOCK /GASCOM/ LENGTH 007640

GASBUF 000000 INTEGER*2 ARRAY (2000)

COMMON BLOCK /ADCOM1/ LENGTH 000104

ADCTBF	000000	INTEGER*2 ARRAY (34)
GASDSP	000000	INTEGER*2 VARIABLE
EXHALE	000002	INTEGER*2 VARIABLE
SVTHRS	000004	INTEGER*2 VARIABLE
FVC	000006	INTEGER*2 VARIABLE
PFT	000016	INTEGER*2 VARIABLE
AHRCNT	000022	INTEGER*2 VARIABLE
HRTHRS	000024	INTEGER*2 VARIABLE
AHR	000032	INTEGER*4 VARIABLE
CWL	000036	INTEGER*4 VARIABLE
CWLCNT	000042	INTEGER*2 VARIABLE
BPTIME	000044	INTEGER*2 VARIABLE
CSPD	000046	INTEGER*4 VARIABLE
BPCNT1	000052	INTEGER*2 VARIABLE
CSBP	000054	INTEGER*4 VARIABLE
CDBP	000060	INTEGER*4 VARIABLE

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
------	--------	------------

BPCNT	000064	INTEGER*2 VARIABLE
-------	--------	--------------------

LFLSHR	000072	INTEGER*2 VARIABLE
--------	--------	--------------------

COMMON BLOCK /X/	LENGTH	000060
------------------	--------	--------

MD	000000	INTEGER*2 ARRAY (24)
----	--------	----------------------

```

C-- FORCED VITAL CAPACITY SUBROUTINE, VERSION 2.0
C--
0001 SUBROUTINE PFT
C--
C-- SET UP VARIABLES
C--
0002 INTEGER AGE
C-- !PATIENTS AGE AND HEIGHT, TO BE USED IN
C-- !THE CALCULATIONS OF THE PREDICTED VALUES.
0003 INTEGER CPID
C-- !PATIENTS CARDIOPULMONARY IDENTITY NUMBER.
0004 INTEGER NAME(14)
C-- !PATIENT'S NAME.
0005 INTEGER SSNO(3)
C-- !PATIENT'S SOCIAL SECURITY NUMBER.
0006 INTEGER DATE(3)
C-- !DATE THE TEST WAS CONDUCTED.
0007 INTEGER RETNO
C-- !PATIENT'S RETEST NUMBER.
0008 INTEGER UNIGNO(2)
C-- !PATIENT'S UNIQUE NUMBER.
0009 INTEGER M,F,SEX
C-- !SUBJECTS SEX TO BE USED IN CALCULATIONS.
0010 INTEGER ADCTBF(34)
C-- !COMMON AREA FOR A/D COMMAND VARIABLES
0011 INTEGER GASBUF(2000)
C-- !A/D COMMON DATA BUFFER.
0012 INTEGER GASDSP
C-- !GAS BUFFER DISPLACEMENT
0013 INTEGER SVTHRS
C-- !SPIROMETER VOLUME THRESHOLD.
0014 INTEGER FVC
C-- !FLAG (=1 FOR FORCED VITAL CAPACITY
C !EQUAL TO 0 FOR ALL OTHERS).
0015 INTEGER PFT
C-- !FLAG =1
0016 INTEGER PUSHES
C-- !CONTAINS 1'S FOR BUTTONS THAT ARE
C !PUSHED.
0017 INTEGER PCSR
C-- !ADDRESS OF THE PROGRAMABLE CLOCK
C !STATUS REGISTER.
0018 INTEGER PCCB
C-- !ADDRESS OF THE PROGRAMABLE CLOCK
C !COUNTER BUFFER.
0019 INTEGER EXHALE
C-- !DONE FLAG SET BY THE A/D INTERRUPT
C !ROUTINE.
0020 INTEGER PCS
C-- !HAS THE VALUE TO BE PLACED INTO THE
C !CLOCK STATUS REGISTER.
0021 INTEGER PCC
C-- !HAS THE VALUE TO BE PLACED INTO THE
C !CLOCK COUNTER BUFFER.

```

```
0022      INTEGER OLD
C--      !VARIABLE TO HOLD THE OLD VALUE OF EXHAST.
0023      INTEGER DROUTB,CRINB
C--      !VARIABLES THAT CONTAIN THE ADDRESS OF THE
C--      !DR11 PARRELL INTERFACE CARD INPUT AND OUTPUT BUFFERS.
0024      REAL HIGHT,BPTS
C--      !PATIENT'S HEIGHT IN CM.
0025      REAL SLOPE,YNTRCP
C--      !VARIABLES THAT CONTAIN INFORMATION PASSED
C--      !THROUGH GENCOM COMMON BLOCK TO CALCULATE THE
C--      !SPIROMETER THRESHOLD.
C--
C--      SET UP THE COMMON AREAS.
C--
0026      COMMON/PDATA/IDATA(1536)
C--      !SET UP THE COMMON BUFFER FOR PATIENT DATA
C--      !ALSO, CONTAINES THE PATIENTS INFORMATION.
0027      COMMON /GASCON/GASELF
C--      !DATA BUFFER COMMON AREA.
0028      COMMON /ADCOM1/ADCTBF
C--      !COMMON AREA TO PASS COMMAND VARIABLES AND
C--      !CONSTANTS.
0029      COMMON/X/MD(24)
0030      EQUIVALENCE (ICATA(85),SLOPE)
C--      !CONTAINES THE SLOPE FROM CALIBRATION ROUTINE.
0031      EQUIVALENCE (ICATA(87),YNTRCP)
C--      !CONTAINES THE Y-INTERCEPT FROM THE CALIBRATION ROUTINE.
0032      EQUIVALENCE(ICATA(3),SSNC(1))
C--      !THESE THREE DATA WORDS CONTAIN
C--      !THE PATIENT'S SOCIAL SECURITY NUMBER.
0033      EQUIVALENCE (IDATA(7),UNIGNC(1))
C--      !THIS DATA WORD CONTAINES THE PATIENT'S
C--      !UNIQUE NUMBER.
0034      EQUIVALENCE (ICATA(9),RETNOC)
C--      !THIS DATA WORD CONTAINES THE SUBJECT'S
C--      !RETEST NUMBER.
0035      EQUIVALENCE (ICATA(10),CPID)
C--      !THIS DATA WORD CONTAINES THE PATIENT'S
C--      !CARDIOPULMUNARY NUMBER.
0036      EQUIVALENCE(IDATA(40),DATE(1))
C--      !THESE THREE DATA WORDS CCNTAIN THE
C--      !PATIENT'S TEST DATE.
0037      EQUIVALENCE (IDATA(26),NAME(1))
C--      !THE 14 DATA WORDS CONTAIN THE
C--      !PATIENT'S NAME.
0038      EQUIVALENCE (ICATA(2),SEX)
CC      !THIS DATA WORD CONTAINES THE PATIENTS SEX.
0039      EQUIVALENCE (ICATA(25),AGE)
C--      !THIS DATA WORD CONTAINES THE PATIENTS AGE.
0040      EQUIVALENCE (ICATA(50),HIGHT)
C--      !THIS DATA WORD CONTAINES THE PATIENTS HEIGHT.
0041      EQUIVALENCE (ICATA(75),BPTS)
0042      EQUIVALENCE (ADCTBF(1),GASCSP)
C--      GAS BUFFER DISPLACEMENT
```

```

0043      EQUIVALENCE (ADCTBF(2),EXFALE)
C--      !DONE FLAG SET BY THE A/D INTERRUPT ROUTINE.
0044      EQUIVALENCE (ADCTBF(3),SVTHRS)
C--      !SPIROMETER VOLUME THRESHOLD, MUST BE CALCULATED.
0045      EQUIVALENCE (ADCTBF(4),FVC)
C--      !FLAG EQUAL TO -1 FOR FORCED VITAL CAPACITY.
0046      EQUIVALENCE (ADCTBF(7),PUSHES)
C--      !CONTAINS 1'S FOR BUTTONS THAT ARE
C        !PUSHED BY PROCESS.
0047      EQUIVALENCE (ADCTBF(8),PFT)
C--      !FLAG FOR THE PFT TEST, MUST BE SET TO -1.
0048      EQUIVALENCE (ADCTBF(30),LFLSHR)
C--      !CONTAINS LIGHT TO FLASH
0049      EQUIVALENCE (IDATA(129),RFVC)
C--      !CONTAINS THE FORCED VITAL CAPACITY IN THE
C--      !PATIENTS DATA BUFFER, WHICH WILL BE STORED
C--      !ON THE FLOPPY DISK.
0050      EQUIVALENCE (IDATA(131),FEV1)
C--      !CONTAINS THE FORCED EXPIRATORY VOLUME IN 1 SEC.
0051      EQUIVALENCE (IDATA(133),PEFR)
C--      !CONTAINS THE PEAK EXPIRATORY FLOW RATE.
0052      EQUIVALENCE (IDATA(135),FEF)
C--      !CONTAINS THE FORCED EXPIRATORY FLOW BETWEEN 25% AND 75% MAX FVC.
0053      EQUIVALENCE (IDATA(137),FEF212)
C--      !CONTAINS THE MEANS FORCED EXPIRATORY FLOW BETWEEN 200ML. AND 1200
0054      EQUIVALENCE (IDATA(141),F1FP)
C--      !CONTAINS THE % FEV1/FVC.
0055      EQUIVALENCE (IDATA(145),PFVC)
C--      !CONTAINS THE PREDICTED FORCED VITAL CAPACITY.
0056      EQUIVALENCE (IDATA(147),OPFVCP)
C--      !CONTAINS THE % OBSERVED FVC/ PREDICTED FVC.
0057      EQUIVALENCE (IDATA(149),PFEV1)
C--      !CONTAINS THE PREDICTED FORCED EXPIRATORY VOLUME IN 1 SEC.
0058      EQUIVALENCE (IDATA(151),CPFEV1)
C--      !CONTAINS THE % OBSERVED FEV1 / PREDICTED FEV1.
C--
C--      INITIALIZE PARAMETERS
C--
0059      DATA M/'M'/,F/'F'/
C--      !SET VARIABLES M AND F EQUAL TO THEIR ACSII VALUES.
0060      DATA DROUTB/'167772'/
C--      !DATA OUT BUFFER FOR THE PARRELL INTERFACE CARD.
0061      DATA DRINB/'167770'/
C--      !DATA IN BUFFER FOR THE PARRELL INTERFACE CARD.
0062      DATA PCSR/'172540'/
C--      !OCTAL ADDRESS OF THE PROGRAMABLE CLOCK(KW11-P)
C        !STATUS REGISTER.
0063      DATA PCCB/'172542'/
C--      !OCTAL ADDRESS OF THE PROGRAMABLE CLOCK
C        !COUNTER BUFFER.
0064      DATA PCC/'144'/
C--      !SET CLOCK TO GET 100 INTERRUPTS/SEC.
0065      DATA PCS/'113'/
C        DATA MD/'JA', 'N ', 'FE', 'B ', 'MA', 'R ', 'AP', 'W ',

```

```

C      @'MA','Y ','JU','N ','JU','L ','AU','G ','SE','P ',
C      @  'OC','T ','NO','V ','CE','C '
C--    !OCTAL VALUE TO BE LOADED INTO THE CLOCK STATUS
C      !REGISTER. BITS 0,1,3, AND 6 ARE TURNED ON.
C      !    BIT 6 - ALLOWS DCNE TO CAUSE AN INTERRUPT.
C      !    BIT 3 - SELECTS REPEATED INTERRUPT MODE.
C      !    BIT 1 - SET CLOCK RATE TO 10KHZ.
C      !    BIT 0 - STARTS THE CLOCK.
0066    FVC=-1
C--    !SET A/D FOR FORCED VITAL CAPACITY DATA ACQUISITION.
0067    PFT=-1
C--    !INITILIZE THE PFT FLAG.
0068    SVTHRS=IFIX((0.1-YNTRCP)/SLOPE)
C--    !CALCULATE THE SPIROMETER THRESHOLD .
C--
C--    SET UP CLOCK STATUS AND COUNTER BUFFER.
C--
0069    CALL IPOKE(DRCUTE,0)
C--    !INITILIZE THE DRV11 OUT BLFFER TO 0.
0070    10    CONTINUE
0071    CALL INITI
C--    !SET TIMERS TO INITIAL VALLES.
0072    EXHALE=0
0073    OLD=0
C--    !SET FLAGS TO 0
0074    CALL LIGHT ("100)
0075    LFLSHR="100
C--    !SET THE BIT TO CAUSE THE FVC LIGHT TO FLASH.
0076    CALL IPOKE(PCCB,PCC)
C--    !SET CLOCK COUNTER FOR 100 INTERRUPTS PER. SECOND.
0077    CALL IPOKE(PCSR,PCS)
C--    !TURN THE CLOCK ON. IPOKE PUTS THE VALUES INTO THE
C      !SPECIFIED ADDRESS.
C
C--    CHECK FOR THE DONE FLAG OR USER TERMINATION
C--
0078    TI=SECNDS(0.0)
C--    !SET START TIME TO 0 FOR ELAPSED TIME CHECK.
0079    100    DELTA=SECNDS(TI)
C--    !GET THE ELAPSED TIME.
0080    IF(EXHALE.EQ.0.AND.OLD.EQ.-1) GO TO 1000
C--    !TEST FOR DONE FLAG FROM THE A/D.
C--    !BRANCH TO RETURN ON SUCESSFUL COMPLETION.
0082    OLD=EXHALE
C--    !SET OLD EQUAL TO THE CURRENT VALUE OF EXHALE.
0083    IF(DELTA.LT.10.0) GO TO 100
C--    !CHECK FOR 10 SECONDS OF ELAPSED TIME. IF NOT GREATER THAN
C--    !10 BRANCH TO 10 TO RECHECK.
0085    TYPE 200
0086    200    FORMAT(/,1X,'FORCED VITAL CAPACITY TEST TIMED OUT',
X/,1X,'PLEASE RETEST SUBJECT',64(/))
0087    CALL IPOKE(PCSR,0)
C--    !TURN THE CLOCK OFF.
0088    GO TO 2200

```

```
C-- !BRANCH TO 10, FVC TEST TIMED OUT, TRY IT AGAIN.
0089 1000 CONTINUE
0090      CALL IPOKE(PCSR,0)
C-- !TURN THE CLOCK OFF.
D      TYPE 999,(GASBUF(1),I=1,GASDSP)
0091 999  FORMAT(/,1X,8(2X,I6))
0092      IF(GASDSP.GT.100)GO TO 1050
0094 1010 TYPE 1001
0095 1001  FORMAT(/,1X,'EXPIRATION BREATH TOO SHORT;',
X/,1X,'PLEASE RETEST SUBJECT',64(/))
D      TYPE 1002
0096 1002  FORMAT(/,'S TYPE IN THE S. VOL. THRESHOLD ')
D      ACCEPT 1003,SVTFRS
0097 1003  FORMAT(I6)
0098      GO TO 2200
0099 1050 CONTINUE
0100      CALL LIGHT("100")
C--
C-- GET THE FVC.
C--
0101      GASAVG=0
C-- INITILIZE THE GAS AVERGER TO 0.
0102      IPKVL=-2047
C-- !SET PEKK VALUE TO 0
D      DO 1100 IPK=1,GASDSP
0103      IF(GASBUF(IPK).LT.IPKVL)GO TO 1100
0104      IPKVL=GASBUF(IPK)
0106      IPEAK=IPK
0107      IPEAK=IPK
0108 1100 CONTINUE
C-- !THIS DO LOOP GET MAX PEEK VALUE
0109      IF(IPEAK.LT.100) GO TO 1010
0111      DO 1150 IAVG=IPEAK-3,IPEAK
C-- SET UP TO GET 4 VALUES
0112 1150 GASAVG=GASAVG+FLOAT(GASBUF(IAVG))
D      TYPE 1151,IPEAK,GASAVG
0113 1151  FORMAT(/,1X,' IPEAK  GASAVG ',I5,5X,F10.2)
0114      RFVC=((GASAVG/4.)*SLOPE+YNTRCP)*8PTS
C-- !GET THE AVERAGE OF THE FVC
C-- !AND NORMALIZE THE FVC.
C--
C-- GET THE FORCED EXPIRATORY VOLUMN IN 1 SECOND.
C--
0115      FEV1=(FLOAT(GASBUF(100))*SLOPE+YNTRCP)*8PTS
C-- !THE FEV1 IS EQUAL TO THE 100 TM.
C-- !DATA WORD. NORMALEZED.
C--
C-- GET THE PEAK EXPIRATORY FLCW RATE (PEFR)
C--
0116      PEFR=-2047
C-- !SET THE INITIAL VALUE TO ZERO.
0117      DO 1200 IP=1,IPEAK
C-- !ANALYZE 1 SEC. CF DATA.
0118      TAN=FLOAT(GASBUF(IP+4)-GASBUF(IP))/0.04
C-- !CALCULATE THE TANGANT FOR EACH TIME INTERVAL.
```

c-2

```
C-- 11 DATA POINT REPRESENTS 10 MSEC. REAL TIME.
0119 IF(TAN.GT.PEFR)PEFR=TAN
C-- 1FIND THE MAXIMUM TANGANT VALUE,
C-- 1THIS IS THE PEFR.
0121 1200 CONTINUE
0122 PEFR=(PEFR*SLOPE+YNTRCP)*8PTS
C-- 1STORE THE PEFR VALUE IN THE CORRECT BUFFER LOCATION.
C--
C-- GET THE FORCED EXPIRATORY FLOW BETWEEN 25% AND 75% MAX FVC.
C-- 1(FORMERLY MMEF)
C--
C-- FEF 25% - 75%
C--
0123 ITI75=1
0124 ITI25=1
C-- 1SET INITIAL TIMES FOR THE 25% AND 75%
C-- 1VALUES TO 0
0125 FVC75=(RFVC*0.75)
C-- 1CALCULATE 75% OF FVC.
0126 FVC25=(RFVC*0.25)
C-- 1CALCULATE 25% OF FVC.
0127 DO 1300 IF=1,IPEAK
C-- 1SET UP TO ANALYZE 1 SEC. CF TIME.
0128 RNUM=FLOAT(GASBLF(IF))*SLOPE+YNTRCP
0129 IF(RNUM.LE.FVC75)ITI75=IF
C-- 1GET THE TIME WHEN THE FVC IS 75% MAXIMUM.
0131 IF(RNUM.LE.FVC25)ITI25=IF
C-- 1GET THE TIME WHEN THE FVC IS 25% MAXIMUM.
0133 1300 CONTINUE
D TYPE 1301,FVC75,FVC25,ITI75,ITI25
0134 1301 FORMAT(/,1X,' FVC75 FVC25 ',2(2X,F10.2),
X/,1X,' ITI75 ITI25 ',2(7X,15))
0135 FEF=(FLOAT(GASBLF(ITI75)-GASBLF(ITI25))*SLOPE)/
X(FLOAT(ITI75-ITI25)*0.01)
0136 FEF=FEF*8PTS
C-- 1CALCULATE THE FEF 25% - 75%
C-- 1BY CALCULATING THE SLOPE.
C-- 1STORE THE EFE VALUE IN THE CORRECT BUFFER LOCATION.
C--
C-- GET THE MEANS FORCED EXPIRATORY FLOW BETWEEN 200 ML. AND 12000 ML.
C-- 1OF THE FORCED VITAL CAPACITY, (FEF 200 - 1200)
C-- 1(FORMERLY MEFR).
0137 ITI12=1
0138 ITI2=1
C-- 1SET THE INITIAL TIMES-FOR THE FVC'S AT 200
C-- 1AND 1200 MLS. TO 0.
0139 DO 1400 ITI=1,IPEAK
C-- 1SET UP TO ANALYZE 1. SEC CF REAL TIME.
0140 RNUM=FLOAT(GASBLF(ITI))*SLOPE+YNTRCP
0141 IF(RNUM.LE.0.2)ITI2=ITI
C-- 1GET THE TIME WHEN THE FVC IS EQUAL TO 200 ML.
0143 IF(RNUM.LE.1.2)ITI12=ITI
C-- 1GET THE TIME WHEN THE FVC IS EQUAL TO 1200 ML.
0145 1400 CONTINUE
```



```

D      TYPE 1401,ITI12,ITI2
0146  1401  FORMAT(/,1X,' ITI12  ITI2 ',2(5X,I5))
0147      FEF212=(FLOAT(GASBUF(ITI12)-GASBUF(ITI2))*SLOPE)/
      X(FLOAT(ITI12-ITI2)*0.01)
0148      FEF212=FEF212*8PTS
C--    !CALCULATE THE FEF 200 = 1200 MLS.
C--    !BY CALCULATING THE SLOPE BETWEEN THEM.
C--    !STORE THR FEF212 VALUE IN THE CORRECT BUFFER LOCATION.
C--
C--    PREDICTED FORCED VITAL CAPACITY AND FORCED VITAL
C--    CAPACITY AT 1 MIN. ARE CALCULATED USING THE FORMULAS REFERENCED
C--    BY:
C--    REUBEN M. CHERNIACK, MC.
C--    IN:
C--    PULMONARY FUNCTION TESTING
C--    W.B. SAUNDERS, PP. 243, 1977.
C--
C--    GET THE % FEV1/FVC
C--
0149      F1FP=FEV1/RFVC*100.0
C--    !CALCULATE THE % FEV1/FVC
C--    !STORE THE F1FP IN THE CORRECT BUFFER LOCATION
C--
C--    GET THE PREDICTED FVC
C--
0150      IF(SEX.EQ.M)PFVC=0.06584*HIGHT-0.02954*AGE-5.12451
C--    !FORMULA FOR CALCULATING THE
C--    !PREDICTED FVC FOR MALES.
0152      IF(SEX.EQ.F)PFVC=0.04071*HIGHT-0.02147*AGE-2.56958
C--    !FORMULA FOR CALCULATING THE
C--    !PREDICTED FVC FOR FEMALES.
C--    !STORE THE PFVC VALUE IN THE CORRECT BUFFER LOCATION.
C--
C--    GET THE PERCENT CF THE OBSERVED FVC TO THE PREDICTED FVC
C--
0154      OPFVCP=RFVC/PFVC*100.0
C--    !CALCULATE THE RATIO % OF THE
C--    !OBSERVED TO PREDICTED FVC.
C--    !STORE THE OPFVCP VALUE IN THE CORRECT BUFFER LOCATION.
C--
C--    GET THE PREDICTED FEV1.
C--
0155      IF(SEX.EQ.M)PFEV1=0.0425*HIGHT-0.03509*AGE-2.59946
C--    !FORMULA FOR CACLLATING THE
C--    !PREDICTED FEV1 FOR MALES.
0157      IF(SEX.EQ.F)PFEV1=0.04071*HIGHT-0.02147*AGE-2.56958
C--    !FORMULA FOR CALCULATING THE
C--    !PREDICTED FEV1 FOR FEMALES.
C--    !STORE THE PFEV1 VALUE IN THE CORRECT BUFFER LOCATION.
C--
C--    GET THE PERCENT OBSERVED FEV1 TO THE PREDICTED FEV1.
C--
0159      OPFEV1=FEV1/PFEV1*100.0
C--    !CALCULATE THE % OF THE OBSERVED

```

```

C-- !FEV1 DIVIDED BY THE PREDICTED FEV1.
C-- !STORE THE OPFEV1 VALUE IN THE CORRECT BUFFER LOCATION.
C--
C-- PRINT OUT REPORT
C--
0160      IMM=2*(DATE(1))
0161      IMM1=IMM-1
0162      TYPE 2000,CPID,(UNIGNO(I),I=1,2),RETNO,DATE(2),
      @MD(IMM1),MD(IMM),DATE(3)
0163      TYPE 2050,(NAME(I),I=1,14),SEX,(SSNO(I),I=1,3)
0164      2000  FORMAT(/,32X,'FVC REPORT',/36X,'CCAS - V2.0',//
      $SX,'CPID NO.: ',I5,28X,'UNIQUE NO.: ',A2,I4,/,
      $SX,'RETEST NO.: ',I2,28X,'DATE: ',I2,I1X,2A2,I2,/)
0165      2050  FORMAT($X,'SUBJECT NAME: ',14A2,I1X,'SEX: ',A2,/,
      $SX,'SOC. SEC. NO.: ',I3,'-',I2,'-',I4)
0166      TYPE 2100,RFVC,FEV1,PEFR,FEF,FEF212,F1FP,PFVC,OPFVCP,
      %PFEV1,OPFEV1
0167      2100  FORMAT(////,$X,'FVC (L) ',F8.2,/,
      $SX,'FEV1 (L) ',F8.2,/,
      $SX,'PEFR (L/SEC) ',F8.2,/,
      $SX,'MMFR (L/SEC) ',F8.2,/,
      $SX,'MEFR (L/SEC) ',F8.2,/,
      $SX,'% FEV1/FVC ',F8.2,/,
      $SX,'% PRED FVC (L) ',F8.2,/,
      $SX,'% OBS/PRED FVC ',F8.2,/,
      $SX,'% PRED FEV1 (L) ',F8.2,/,
      $SX,'% OBS/PRED FEV1 ',F8.2,/,
      $SX,'ALL PFT DATA BTPS',////////////////////////////////)
0168      2200  CALL LIGHT('103)      !SET LIGHTS TO IDLE STATUS
0169      RETURN
0170      END

```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
M	000014	INTEGER*2 VARIABLE
F	000016	INTEGER*2 VARIABLE
PCSR	000024	INTEGER*2 VARIABLE
PCCB	000026	INTEGER*2 VARIABLE
PCS	000032	INTEGER*2 VARIABLE
PCC	000030	INTEGER*2 VARIABLE
OLD	001560	INTEGER*2 VARIABLE
DROUTB	000020	INTEGER*2 VARIABLE
DRINB	000022	INTEGER*2 VARIABLE
IFIX	000000	INTEGER*2 PROCEDURE
IPOKE	000000	INTEGER*2 PROCEDURE
INITI	000000	INTEGER*2 PROCEDURE
LIGHT	000000	INTEGER*2 PROCEDURE
TI	001562	REAL*4 VARIABLE
SECNDS	000000	REAL*4 PROCEDURE
DELTA	001566	REAL*4 VARIABLE
GASAVG	001572	REAL*4 VARIABLE
IPKVL	001576	INTEGER*2 VARIABLE
IPK	001600	INTEGER*2 VARIABLE
IPEAK	001602	INTEGER*2 VARIABLE
IAVG	001604	INTEGER*2 VARIABLE
FLOAT	000000	REAL*4 PROCEDURE
IP	001606	INTEGER*2 VARIABLE
TAN	001610	REAL*4 VARIABLE
ITI75	001614	INTEGER*2 VARIABLE
ITI25	001616	INTEGER*2 VARIABLE
FVC75	001620	REAL*4 VARIABLE
FVC25	001624	REAL*4 VARIABLE
IF	001630	INTEGER*2 VARIABLE
RNUM	001632	REAL*4 VARIABLE
ITI12	001636	INTEGER*2 VARIABLE
ITI2	001640	INTEGER*2 VARIABLE
ITI	001642	INTEGER*2 VARIABLE
IMM	001644	INTEGER*2 VARIABLE
IMM1	001646	INTEGER*2 VARIABLE
I	001650	INTEGER*2 VARIABLE

COMMON BLOCK /PDATA/ LENGTH 006000

IDATA	000000	INTEGER*2 ARRAY (1536)
SLOPE	000250	REAL*4 VARIABLE
YNTRCP	000254	REAL*4 VARIABLE
SSNO	000004	INTEGER*2 ARRAY (3)
UNIGNO	000014	INTEGER*2 ARRAY (2)
RETNO	000020	INTEGER*2 VARIABLE
CPID	000022	INTEGER*2 VARIABLE
DATE	000116	INTEGER*2 ARRAY (3)
NAME	000062	INTEGER*2 ARRAY (14)
SEX	000002	INTEGER*2 VARIABLE
AGE	000060	INTEGER*2 VARIABLE
HIGHT	000142	REAL*4 VARIABLE

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES	
BPTS	000224	REAL*4	VARIABLE
RFVC	000400	REAL*4	VARIABLE
FEV1	000404	REAL*4	VARIABLE
PEFR	000410	REAL*4	VARIABLE
FEF	000414	REAL*4	VARIABLE
FEF212	000420	REAL*4	VARIABLE
F1FP	000430	REAL*4	VARIABLE
PFVC	000440	REAL*4	VARIABLE
OPFVCP	000444	REAL*4	VARIABLE
PFEV1	000450	REAL*4	VARIABLE
OPFEV1	000454	REAL*4	VARIABLE

COMMON BLOCK /GASCOM/ LENGTH 007640

GASBUF 000000 INTEGER*2 ARRAY (2000)

COMMON BLOCK /ADCOM1/ LENGTH 000104

ADCTBF	000000	INTEGER*2	ARRAY (34)
GASDSP	000000	INTEGER*2	VARIABLE
EXHALE	000002	INTEGER*2	VARIABLE
SVTHRS	000004	INTEGER*2	VARIABLE
FVC	000006	INTEGER*2	VARIABLE
PUSHES	000014	INTEGER*2	VARIABLE
PFT	000016	INTEGER*2	VARIABLE
LFLSHR	000072	INTEGER*2	VARIABLE

COMMON BLOCK /X/ LENGTH 000060

MD 000000 INTEGER*2 ARRAY (24)

SUBROUTINE ADISR

VERSION: 2

DATE: 08/24/79

AUTHOR: BILL CHOSIER

PURPOSE: SAMPLE SELECTED A/D CHANNELS, CHECK FRONT SWITCH PANEL (THRU PARALLEL INTERFACE), FLICKER OR FLASH SELECTED FRONT PANEL LIGHTS, AND MAINTAIN ONE-SECOND AND 30-SECOND COUNTERS FOR TIMING APPROPRIATE TEST ACTIVITIES.

SUMMARY OF ROUTINE:

ADISR IS CALLED EACH TIME THE KKV11-A PROGRAMMABLE CLOCK REQUESTS AN INTERRUPT. THE INTERRUPT RATE SHOULD BE SET BY THE MAIN PROGRAM AT THE START OF EACH TEST AT 100 HZ.

EACH TIME ADISR IS CALLED, A SCAN

COUNTER (COUNT) IS DECREMENTED AND ANALOG CHANNEL 0 (SPIROMETER VOLUME) IS SAMPLED AND COMPARED TO A THRESHOLD VOLUME (SVTHRS) WHICH IS AN INTEGER VALUE PREVIOUSLY CALCULATED BY THE MAIN PROGRAM WHICH CORRESPONDS TO A PARTICULAR AMOUNT OF SPIROMETER DISPLACEMENT SUCH AS 0.2 LITERS. IF THE SAMPLED VOLUME IS LESS THAN THE THRESHOLD, NOTHING IS DONE WITH THE SAMPLE, AND THE ROUTINE STORES 0 IN THE EXHALATION FLAG (EXHALE) AND THEN GOES TO READ THE FRONT PANEL SWITCHES.

HOWEVER, IF THE SAMPLED VOLUME IS GREATER THAN OR EQUAL TO THE THRESHOLD, "EXHALE" IS SET TO -1, THE SAMPLE IS STORED IN THE GAS BUFFER (GASBUF) AND ITS POINTER IS INCREMENTED. (THE POINTER IS GENERATED BY ADDING THE STARTING ADDRESS OF GASBUF TO A RELATIVE DISPLACEMENT, GASDSP, ACCESSIBLE TO THE MAIN PROGRAM. AT THE END OF THE GAS SAMPLING THE UPDATED VALUE OF GASDSP IS STORED AGAIN. THIS ELIMINATES THE NEED TO ADD AN INDEX VALUE TO THE BUFFER STARTING ADDRESS FOR EACH GAS VALUE.)

IF AN FVC IS NOT IN PROGRESS, THE PROGRAM THEN SAMPLES ANALOG CHANNELS 1 (O2), 2 (N2), AND 3 (CO2), AND STORES EACH SAMPLE IN SUCCESSIVE LOCATIONS IN GASBUF. AFTERWARDS, THE UPDATED GAS BUFFER DISPLACEMENT POINTER IS STORED IN GASDSP. IT THEN POINTS TO THE NEXT AVAILABLE LOCATION IN GASBUF.

NEXT, THE SCAN COUNTER IS CHECKED TO SEE IF IT EQUALS 0. IF NOT, THE PROGRAM

RESTORES THE REGISTERS AND RETURNS. IF COUNT IS = 0, THEN A NUMBER OF ACTIVITIES TAKE PLACE AS FOLLOWS:

1. THE A/D IS TRIGGERED TO SAMPLE CH. 4 (HR).
2. COUNT IS RE-INITIALIZED TO 10. THIS DETERMINES HOW FREQUENTLY (RELATIVE TO THE INTERRUPT RATE) ALL OF THESE OTHER ACTIVITIES TAKE PLACE.
3. THE FRONT PANEL SWITCHES ARE READ, AND BITS FOR SWITCHES THAT WERE JUST PUSHED, WHICH WERE NOT PREVIOUSLY RECOGNIZED BY THE MAIN PROGRAM, ARE SET IN "PUSHES". THE MAIN PROGRAM SHOULD ACKNOWLEDGE RECEIPT OF THIS INFORMATION BY SAVING THE CONTENTS OF PUSHES AND THEN IMMEDIATELY CLEARING PUSHES BY STORING 0 IN IT.
4. THE PFT FLAG IS CHECKED AND IF NOT=0, THE PROGRAM BRANCHES TO END.
5. HEART RATE (HR) IS SAMPLED AND IF IT IS GREATER THAN HRTHRS, AN INTEGER VALUE CORRESPONDING TO A LOWER HEART RATE THRESHOLD AND PREVIOUSLY DETERMINED BY THE MAIN PROGRAM, THEN THE SAMPLED VALUE IS ADDED TO THE CURRENT HR (CHR) AND AVERAGE HR (AHR) ACCUMULATORS, AND THEIR CORRESPONDING COUNTERS (CHRCNT & AHRCNT) ARE INCREMENTED.
6. WORK LOAD (WL) OR ELEVATION (ELEV) AND SPEED ARE NEXT

SAMPLED. THE SAMPLED DATA IS ADDED TO THE APPROPRIATE ACCUMULATORS (CWL & CSPD) AND THE WL-ELEV-SPEED COUNTER (CWL CNT) IS INCREMENTED.

7. BPTIME, THE BP TIMER, IS DECREMENTED, AND AT THE START OF THE 55-TH SECOND OF EACH MINUTE AN INTERNAL BP COUNTER IS SET TO COLLECT SBP AND DBP SAMPLES DURING THE NEXT 10 PASSES THRU THIS PORTION OF ADISR (EVERY 0.1 SEC). IN ORDER TO WORK PROPERLY, BPTIME MUST BE INITIALIZED BY THE MAIN PROGRAM AT THE BEGINNING OF THE EXERCISE TEST, JUST BEFORE THE PROGRAMMABLE CLOCK IS STARTED, TO $10 \times (55 - \text{SECS})$, WHERE SECS IS THE NUMBER OF SECONDS PAST THE BEGINNING OF THE CURRENT MINUTE. OTHERWISE THE BLOOD PRESSURE CHANNELS WILL NOT BE SAMPLED AT THE CORRECT TIME.

THE SBP AND DBP VALUES ARE ADDED TO THEIR RESPECTIVE ACCUMULATORS, CSBP AND CDBP. THE MAIN PROGRAM SHOULD DIVIDE THE SUMS BY 10 (RATHER THAN BY A VARIABLE COUNTER AS FOR HR, WL-ELEV, AND SPEED) IN ORDER TO GET THE AVERAGE BLOOD PRESSURES FOR THE MINUTE.

8. AFTER SPEED IS SAMPLED (AND CONDITIONALLY, SBP & DBP), A 0.1 SECOND COUNTER IS DECREMENTED, AND EVERY 1 SECOND SECS IS INCREMENTED BY ONE. ANOTHER 1-SECOND COUNTER IS ALSO DECREMENTED, AND EVERY 30 SECS. SEC30 IS INCREMENTED BY 1. SECS & SEC30 MAY BE USED BY OTHER ROUTINES TO TIME OTHER ACTIVITIES AND MAY EITHER BE LEFT TO RUN CONTINUOUSLY, OR THEY MAY BE RESET TO ZERO PERIODICALLY BY OTHER ROUTINES. THESE TIMERS WILL NOT WORK PROPERLY UNLESS THE KMW11-A IS SET TO PROVIDE INTERRUPTS EVERY 10 MSEC.

NOTES: 1. THE ACCUMULATORS FOR CURRENT HR (CHR), AVERAGE HR (AHR), WL-ELEV (CWL), SPEED (CSPD), SBP (CSBP), AND DBP (CDBP) ARE DOUBLE-WORD INTEGERS, ALTHOUGH THEIR COUNTERS ARE SINGLE-WORD INTEGERS. EACH ACCUMULATOR AND ITS RESPECTIVE COUNTER (EXCEPT FOR SBP & DBP) MUST BE RESET TO 0 BY THE MAIN PROGRAM EACH TIME ONE DESIRES TO START A NEW AVERAGING INTERVAL. THIS IS NORMALLY DONE EACH TIME SEC30 IS INCREMENTED, EXCEPT FOR CURRENT HR, WHICH MAY BE RESET EACH TIME SECS IS INCREMENTED.

2. GENERAL REGISTERS USED:

- R0 -- POINTS TO A/D COMMAND & STATUS REGISTER
- R1 -- POINTS TO ONE OF SAMPO THRU SAMPB (A/D TRIGGER COMMAND WORDS)
- R2 -- POINTER FOR SAMPLED A/D DATA AND OTHER CALCULATED PARAMETERS
- R3 -- USED FOR SCAN COUNTER (RESET EACH TIME FRONT PANEL SWITCHES AND SLOW A/D SAMPLES ARE READ), ACKNOWLEDGED FRONT PANEL SWITCHES, AND SIGN EXTENSIONS OF A/D SAMPLES FOR HR, WL-ELEV, SPEED, SBP, AND DBP.
- R4 -- CONTAINS A/D SAMPLED DATA TEMPORARILY (EXCEPT O2, N2, AND CO2), AND FRONT PANEL SWITCH CONDITIONS
- R5 -- CONTAINS MASK FOR CHECKING A/D DONE FLAG.

```

1      ; TITLE ADISR
2      ; PROGRAM: ADISR.MAC
3      ; VERSION: 2.0
4      ; ANALOG TO DIGITAL CONVERTER INTERRUPT SERVICE ROUTINE FOR CDAS
5      ; AUTHOR: WILLIAM G. CROSBIE
6      ; DATE: 24 AUG 79
7
8
9
10     ;
11     ;
12     ;
13     ;
14     ;
15     ;
16     ;
17     ;
18     ;
19     ;
20     ;
21     ;
22     ;
23     ;
24     ;
25     ;
26     ;
27     ;
28     ;
29     ;
30     ;
31     ;
32     ;
33     ;
34     ;
35     ;
36     ;
37     ;
38     ;
39     ;
40     ;
41     ;
42     ;
43     ;
44     ;
45     ;

```

.TITLE ADISR
 PROGRAM: ADISR.MAC
 VERSION: 2.0
 ANALOG TO DIGITAL CONVERTER INTERRUPT SERVICE ROUTINE FOR CDAS
 AUTHOR: WILLIAM G. CROSBIE
 DATE: 24 AUG 79

.MCALL .REGCEP,..V2..
 ..V2..
 .REGCEP

ACSR=176770 ; A/D CONTROL & STATUS REGISTER
 ACIA=176772 ; A/D INPUT DATA BUFFER
 DAC1=176760 ; D/A #1 OUTPUT DATA BUFFER
 DAC2=176762 ; D/A #2 OUTPUT DATA BUFFER
 CRCLTR=167772 ; FRONT PANEL OUTPUT BUFFER (LIGHTS)
 CHIAN=167774 ; FRONT PANEL INPUT BUFFER (SWITCHES)
 CLKCSR=172540 ; PROGRAMMABLE CLOCK CONTROL & STATUS REG.

ADISR: MCV R0,-(SP) ; SAVE REGISTERS ON STACK
 MCV R1,-(SP) ;
 MCV @ACSR,R0 ; PUT ADDRESS OF A/D CSR IN R0

MCV @SAMPC,R1 ; PUT A/D TRIGGER WORD POINTER IN R1

MCV (R1)+,(R0) ; TRIGGER A/D FOR CHAN. 0 (SPIR. VOL.)
 MCV R2,-(SP) ; SAVE OTHER REGISTERS
 MCV R3,-(SP)
 MCV R4,-(SP)
 MCV R5,-(SP)
 MCV @200,R5 ; STORE A/D DONE MASK IN R5

MCV @GASPLF,R2 ; PUT GAS DATA BUFFER STARTING ADR. IN R2

MCV GASDSP,R3 ; PUT GAS DISPLACEMENT POINTER (WORD COUNT) IN R3

CMP R3,#2000. ; SEE IF GAS DATA BUFFER IS FULL

BMI GETDSP ; IF NOT, GO CALCULATE DISPLACEMENT POINTER.
 CLR R3 ; RESET R3 TO 0.
 CLR GASDSP ; RESET GAS DISPLACEMENT TO 0.

IAC GASCVF ; INCREMENT GAS DISPLACEMENT OVERFLOW FLAG.

GETDSP: ASL R3 ; MULT. R3 BY 2 (CONVERT WORD TO BYTE COUNT)
 ACC R3,R2 ; & ADD REL. DISPLACEMENT TO GENERATE POINTER
 MCV COUNT,R3 ; LOAD SCAN COUNTER

DEC R3 ; DECREMENT COUNTER
 MCV R3,CCLAT ; & STORE IT (LEAVE IT IN R3 FOR LATER USE)

LOOP0: BIT (R0),R5 ; SEE IF A/D CONVERSION IS DONE
 BEC LOOPC ; IF NOT, WAIT
 MCV @ACIA,R4 ; PUT A/D SAMPLE IN R4

TST EXHALE ; EXHALATION IN PROGRESS DURING LAST SAMPLE?

```

06 00114 001431      PEG      NOEXH      ; IF NOT, GO TO NOEXH
07 00116 020467      CMP      R4,SVTHRS   ; IS SAMPLE .GT. THRESHOLD?
      000004
08 00122 100013      SPL      HIGH1      ; IF SO, GO RESET LOCNT & CHECK FVC FLAG
09 00124 005367      DEC      LOCAT      ; IF SAMPLE .LE. THRESHOLD, DECREMENT LOCNT
      000574
10 00130 002013      BGE      FVCTST      ; IF LOCNT IS STILL .GE. 0, GO TO FVCTST
11 00132 005067      CLM      EXHALE      ; OTHERWISE, SET EXHALE FLAG = 0 (FALSE)
      000002
12 00134 005267      INC      EXHCA1      ; INCREMENT END OF BREATH EXHALATION COUNTER
      000012
13 00142 012767      MCV      #5,LOCAT    ; RESET LOCNT TO 5
      000005
      000550
14 00150 000456      BR       READSH      ; GO READ FRONT PANEL SWITCHES
15 00152 012767 HIGH1: MCV      #5,LOCNT    ; RESET LOCNT TO 5
      000005
      000540
16 00160 005767 FVCTST: TST      FVC      ; IS THIS AN FVC TEST?
      000006
17 00164 001430      BGE      EXH        ; IF NOT, GO TO EXH
18 00166 010422 STORSH: MCV      R4,(R2)+   ; IF FVC, STORE SP. VOL. SAMPLE IN GAS BUFFER
19 00170 012767      MCV      #=-1,EXHALE  ; SET EXHALE = -1 (TRUE)
      177777
      000002
20 00176 000443      BR       READSH      ; GO READ FRONT PANEL SWITCHES
21 00200 020467 NOEXH: CMP      R4,SVTHRS   ; IS SAMPLE .GT. THRESHOLD?
      000004
22 00204 100004      SPL      HIGH      ; IF SO, GO TO HIGH
23 00206 012767      MCV      #5,HICNT    ; OTHERWISE, RESET HICNT TO 5
      000005
      000506
24 00210 000434      BR       READSH      ; GO READ FRONT PANEL SWITCHES
25 00216 005767 HIGH:  TST      FVC        ; IS THIS AN FVC TEST?
      000006
26 00222 001361      BGE      STORSH      ; IF SO, GO STORE SPIR. VOL. SAMPLE
27 00224 005367      DEC      HICNT      ; DECREMENT HICNT
      000472
28 00230 002026      BGE      READSH      ; IF HICNT IS STILL .GE. 0, GO TO READSH
29 00232 012767      MCV      #=-1,EXHALE  ; OTHERWISE, SET EXHALE = -1 (TRUE)
      177777
      000002
30 00240 012767      MOV      #5,HICNT    ; RESET HICNT TO 5
      000005
      000454
31 00246 010422 EXH:   MCV      R4,(R2)+   ; STORE SPIR. VOLUME SAMPLE IN GAS BUFFER
32 00250 012110      MCV      (R1)+,(R0)   ; TRIGGER A/D FOR CH. 1 (O2)
33 00252 031005 LOOP1: BIT      (R0),R5    ; SEE IF CONVERSION DONE
34 00254 001776      BEC      LOOP1      ; IF NOT WAIT
35 00256 013722      MCV      #=ACIA,(R2)+ ; STORE A/D DATA (O2)
      176772
36 00262 012110      MCV      (R1)+,(R0)   ; TRIGGER A/D FOR CH. 2 (N2)
37 00264 031005 LOOP2: BIT      (R0),R5    ; SEE IF CONVERSION DONE
38 00266 001776      BEC      LOOP2      ; IF NOT WAIT
39 00270 013722      MCV      #=ACIA,(R2)+ ; STORE A/D DATA (N2)
      176772
40 00274 012110      MOV      (R1)+,(R0)   ; TRIGGER A/D FOR CH. 3 (CO2)

```



```

01 00274 031005 LOOP3: BIT      (R0),R5      ; SEE IF CONVERSION DONE
02 00300 001776 SEC      LOOP3          ; IF NOT WAIT
03 00302 013722 MCV      @RACIA,(R2)+      ; STORE A/D DATA (C02)
                                176772
04 00304 012701 READSH: MCV      @SAMP4,R1    ; PUT CH. 4 A/D TRIGGER WORD POINTER IN R1
                                000010
05 00312 162702 SBR      @GASBUF,R2        ; SUBTRACT GAS BUF. STARTING ADDR. FROM
                                000000
                                ; GAS POINTER
06
07 00314 006202 ASR      R2                ;
08 00320 010267 MCV      R2,GASBUF          ; STORE NEW GAS DATA BUFFER DISPLACEMENT
                                000000
09 00324 005703 TST      R3                ; SCAN COUNTER = 0 ?
10 00326 003165 SET      FINISH          ; IF NOT, GO TO END
11 00330 012767 MCV      @R10,COUNT        ; RE-INITIALIZE SCAN COUNTER
                                000012
                                000000
12
13 00334 012702 MCV      @R0R18,R4          ; READ FRONT PANEL SWITCHES
14 00336 000010 MCV      @PFT,R2                ; PUT DATA POINTER IN R2
                                000010
15
16 00340 000000 MCV      @MASK,R2          ; PUT POINTER TO MASK IN R2
17 00342 000000 BIC      (R2)+,R4          ; MASK OUT UNUSED BITS
18 00344 000000 CCM      R4                ; COMPLEMENT OF CURRENT SWITCH CONDITIONS
19 00346 000000 MCV      (R2),R3          ; PUT PREVIOUS SWITCH CONDITIONS IN R3
20 00348 000000 BIC      R4,R3            ; CLEAR BITS IN R3 FOR SWITCHES THAT
21 00350 000000 ; ARE OFF (NOT PUSHED).
22 00352 000000 CCM      R4                ; RESTORE CURRENT SWITCH CONDITIONS IN R4
23 00354 000000 BIC      R3,R4            ; CLEAR BITS IN R4 FOR PREVIOUSLY RECOGNIZED
24 00356 000000 ; SWITCH PUSHERS. NOW R4 HAS 1'S ONLY FOR NEWLY
25 00358 000000 ; PUSHED SWITCHES, R0:
26 00360 000000 ; REPRESENT SWITCHES AND (.NOT.PREVS)
27 00362 000000 BIC      R4,R3            ; SET BITS IN R3 FOR NEW SWITCH PUSHERS, R0
28 00364 000000 ; R3 NOW HAS PRESENT SWITCH CONDITIONS
29 00366 000000 ; (SAME AS @R0R18)
30 00368 000000 MCV      R3,(R2)+          ; STORE PRESENT SWITCH CONDITIONS IN PREVS
31 00370 000000 BIS      R4,(R2)+          ; STORE NEW SWITCH PUSHER INFO IN PUSHERS
32 00372 000000 TST      (R2)+            ; IS PFT IN PROGRESS ?
33 00374 000000 BAE      FLICKR          ; IF SO, GO TO END
34 00376 000000 MCV      (R1)+,(R0)        ; TRIGGER A/D FOR CH. 4 (MR)
35 00378 000000 IAC      (R2)+            ; INCREMENT CURRENT MR POINTER
36 00380 000000 IAC      (R2)+            ; INCREMENT 30- SEC AVG. MR COUNTER
37 00382 031005 LOOP4: BIT      (R0),R5      ; SEE IF A/D CONVERSION DONE
38 00384 001776 SEC      LOOP4          ; IF NOT WAIT
39 00386 013704 MCV      @RACIA,R4          ; PUT MR A/D SAMPLE IN R4
                                176772
40 00388 006703 SXT      R3                ; EXTEND SIGN INTO R3
41 00390 012110 MCV      (R1)+,(R0)        ; TRIGGER A/D FOR CH. 5 (WL-ELEV)
42 00392 022204 CMP      (R2)+,R4          ; MR THRESHOLD .GT. SAMPLE ?
43 00394 100407 BHI      ADDMR          ; SKIP IF NOT
44 00396 005367 DEC      CMOCAT          ; DECREMENT MR COUNTERS SINCE MR
                                000020
45 00398 005367 DEC      AMRCAT          ; SAMPLE WAS INVALID (.LT. THRESHOLD).
                                000022
46 00400 012702 MCV      @CWL,R2          ; PUT POINTER TO CWL IN R2
                                000036
47 00402 000406 BR      LCOPI          ; GO READ WL-ELEV SAMPLED VALUE
48 00404 000422 ADDMR: ACC      R4,(R2)+    ; ACCUMULATE CURRENT MR SUM

```

127	0410	005512	ACC	(R2)	: ADD CARRY
128	0416	000322	ACC	R3,(R2)+	: ADD HIGH-ORDER PART OF CWR
129	0420	000422	ACC	R4,(R2)+	: ACCUMULATE 30-SEC AVG. WR SUM
130	0422	005503	ACC	R3	: ADD CARRY
131	0424	000322	ACC	R3,(R2)+	: ADD HIGH-ORDER PART OF AMR
132	0426	031005	LOOP5: BIT	(R0),R5	: SEE IF A/D CONVERSION DONE
133	0430	001776	BEC	LOOP5	: IF NOT WAIT
134	0432	013704	MOV	SRAC1A,R4	: PUT WL-ELEV A/D SAMPLE IN R4
		176772			
135	0434	006703	EXT	R3	: EXTEND SIGN INTO R3
136	0440	012110	MOV	(R1)+,(R0)	: TRIGGER A/D FOR CH. 6 (SPEED)
137	0442	000447	ACC	R4,ANL	: ACCUMULATE CURRENT WL SUM
		000476			
138	0446	005567	ACC	ANL1	: ADD CARRY
		000100			
139	0452	000367	ACD	R3,ANL1	: ADD HIGH-ORDER PART OF ANL
		000100			
140	0456	005267	IAC	ANL,CAT	: INCREMENT CURRENT COUNTER
		000102			
141	0462	000422	ACC	R4,(R2)+	: ACCUMULATE 30 SEC AVG WL-ELEV SUM
142	0464	005503	ACC	R3	: ADD CARRY
143	0466	000322	ACC	R3,(R2)+	: ADD HIGH-ORDER PART OF CML
144	0470	005222	IAC	(R2)+	: INCREMENT 30 SEC AVG WL-ELEV COUNTER
145	0472	005322	DEC	(R2)+	: DECREMENT RP TIMER 8
146	0474	100006	BPL	LOOP6	: SKIP IF NOT TIME TO START BP READING
147	0476	012767	MOV	R400..BPTIME	: RESET BP TIMER TO 60 SEC = 10 SAMP/SEC
		001130			
		000044			
148	0500	016767	MOV	BPCAT,BPCNT1	: RESET BP COUNTER FOR BPCNT SAMPLES
		000064			
		000052			
149	0512	031005	LOOP6: BIT	(R0),R5	: SEE IF A/C CONVERSION DONE
150	0514	001776	BEC	LOOP6	: IF NOT WAIT
151	0516	013704	MOV	SRAC1A,R4	: PUT A/D SAMPLE (SPEED) IN R4
		176772			
152	0522	006703	EXT	R3	: EXTEND SIGN IN TO R3
153	0524	000422	ACD	R4,(R2)+	: ACCUMULATE CURRENT SPEED SUM
154	0526	005503	ACC	R3	: ADD CARRY
155	0530	000322	ACC	R3,(R2)+	: ADD HIGH-ORDER PARTS
156	0532	005322	DEC	(R2)+	: DECREMENT BP COUNTER
157	0534	100422	BMI	TIMER	: SKIP IF ALREADY HAVE BPCNT SAMPLES
158	0536	012110	MOV	(R1)+,(R0)	: TRIGGER A/D FOR CH. 7 (SRP)
159	0540	031005	LOOP7: BIT	(R0),R5	: SEE IF A/D CONVERSION DONE
160	0542	001776	BEC	LOOP7	: IF NOT WAIT
161	0544	013704	MOV	SRAC1A,R4	: PUT SRP A/D SAMPLE IN R4
		176772			
162	0550	006703	EXT	R3	: EXTEND SIGN INTO R3
163	0552	012110	MOV	(R1)+,(R0)	: TRIGGER A/D FOR CH. 8 (DBP)
164	0554	000422	ACC	R4,(R2)+	: ACCUMULATE CURRENT SRP SUM
165	0556	005503	ACC	R3	: ADD CARRY
166	0560	000322	ACC	R3,(R2)+	: ADD HIGH-ORDER PARTS
167	0562	031005	LOOP8: BIT	(R0),R5	: SEE IF A/D CONVERSION DONE
168	0564	001776	BEC	LOOP8	: IF NOT WAIT
169	0566	013704	MOV	SRAC1A,R4	: PUT SAMPLE (DBP) IN R4
		176772			
170	0572	006703	EXT	R3	: EXTEND SIGN INTO R3
171	0574	000422	ACC	R4,(R2)+	: ACCUMULATE CURRENT DBP SUM

```

172 0576 005503      ACC      R3      ; ADD CARRY
173 0600 000322      ACC      R3,(R2)+ ; ADD HIGH-ORDER PARTS
174 0602 005367 TIMER: DEC      CNT10  ; DECREMENT 0.1 SECOND COUNTER
      000002*
175 0604 003015      BGT      FLICKR   ; SKIP IF 1 SEC. HAS NOT ELAPSED
176 0610 012767      MCV      @10..CNT10 ; RESET CNT10 FOR 10 COUNTS
      000012*
      000002*
177 0616 005267      IAC      SECS     ; INCREMENT SECONDS COUNTER/FLAG
      000006*
178 0622 005367      DEC      CNT30    ; DECREMENT 30-SEC COUNTER
      000004*
179 0624 001005      PAB      FLICKR   ; SKIP IF 30 SECS HAVE NOT ELAPSED
180 0630 012767      MCV      @30..CNT30 ; RESET CNT30 FOR 30 COUNTS
      000036*
      000004*
181 0636 005267      IAC      SECS30   ; INCREMENT 30-SEC. COUNTER/FLAG
      000070*
182
183 0642 005367 FLICKR: DEC      FLPCAT  ; DECREMENT FLICKER COUNTER
      000006*
184 0644 003015      BGT      FINISH   ; SKIP IF NOT TIME TO FLICKER LIGHTS
185 0650 012767      MCV      @2..FLPCAT ; RESET FLPCAT TO CHANGE LIGHT STATUS
      000002*
      000006*
186
187 0656 005167      CCN      FLRFLG   ; EVERY 0.2 SEC. (FLICKER RATE= 2.5HZ)
      000010* ; COMPLEMENT FLICKER FLAG
188 0662 100008      BCL      LITECA   ; EVERY OTHER TIME, GO TURN LIGHTS ON
189 0664 006737      BIC      LFLSR,200R0TB ; TURN OFF APPROPRIATE LIGHTS
      000072*
      167772
190 0672 000403      BR       FINISH   ; GO TO END
191 0674 006737 LITEON: BIS      LFLSR,200R0TB ; TURN ON APPROPRIATE LIGHTS
      000072*
      167772
192
193 0702 012605 FINISH: MCV      (SP)+,R5 ; RESTORE REGISTERS
194 0704 012604      MCV      (SP)+,R4
195 0706 012603      MCV      (SP)+,R3
196 0710 012602      MCV      (SP)+,R2
197 0712 012601      MCV      (SP)+,R1
198 0714 012600      MCV      (SP)+,R0
199 0716 000002      RTI
200
201 0720 000005 LOCAT: .ACRD      5
202 0722 000005 LOCAT: .ACRD      5
203
204 000000* .CSECT TIMCNT
205 0000 000012 COUNT: .ACRD     10. ; SCAN COUNTER
206 0002 000012 CNT10: .ACRD     10. ; 0.1 SEC. COUNTER
207 0004 000036 CNT30: .ACRD     30. ; 1 SEC. COUNTER
208 0006 000000 FLPCNT: .ACRD     0 ; FLICKER COUNTER (FOR LIGHTS)
209 0010 000000 FLRFLG: .ACRD     0 ; FLICKER FLAG
210 0012 000000 EXMCNT: .ACRD     0 ; EXHALATION COUNTER.
211
212 000000* .CSECT AOCCH

```

213	0000	000020	SAMP0:	.ACRD	20	:	CH. 0 (SP. VOL) A/D START COMMAND
214	0002	000420	SAMP1:	.ACRD	420	:	CH. 1 (O2)
215	0004	001020	SAMP2:	.ACRD	1020	:	CH. 2 (A2)
216	0006	001420	SAMP3:	.ACRD	1420	:	CH. 3 (CO2)
217	0010	002020	SAMP4:	.ACRD	2020	:	CH. 4 (HR)
218	0012	002420	SAMP5:	.ACRD	2420	:	CH. 5 (HL-ELEV)
219	0014	003020	SAMP6:	.ACRD	3020	:	CH. 6 (SPEED)
220	0016	003420	SAMP7:	.ACRD	3420	:	CH. 7 (SRP)
221	0020	004020	SAMP8:	.ACRD	4020	:	CH. 8 (DBP)
222						:	
223		000000		.CSECT	ADCCM1		
224	0000	000000	GASCP:	.ACRD	0	:	GAS BUFFER DISPLACEMENT (SUBSCRIPT)
225	0002	000000	EXHALE:	.ACRD	0	:	FLAG (=0 FOR INHALATION, =1 FOR EXHALATION)
226						:	
227	0004	000000	SVTHRS:	.ACRD	0	:	SPIROMETER VOL. THRESHOLD
228	0006	000000	FVC:	.ACRD	0	:	FLAG (=1 FOR FVC, 0 ALL OTHERS)
229	0010	177400	MASK:	.ACRD	177400	:	MASKS OUT UNUSED BITS ON PARALLEL INTERFACE INPUT LINES (UNUSED SWITCHES)
230						:	
231	0012	000000	PREVSN:	.ACRD	0	:	CONTAINS 1'S FOR SWITCHES THAT HAVE ALREADY BEEN RECOGNIZED AS BEING PUSHED
232						:	
233	0014	000000	PUSHES:	.ACRD	0	:	CONTAINS 1'S FOR BUTTONS THAT ARE PUSHED
234	0016	000000	PFT:	.ACRD	0	:	FLAG (=1 FOR ALL PFT'S, 0 OTHERWISE)
235	0020	000000	CHRCNT:	.ACRD	0	:	
236	0022	000000	ANRCNT:	.ACRD	0	:	
237	0024	000000	PRTHRS:	.ACRD	0	:	
238	0026	000000	CHR:	.ACRD	0	:	
239	0030	000000	CHR1:	.ACRD	0	:	
240	0032	000000	ADR:	.ACRD	0	:	
241	0034	000000	ADR1:	.ACRD	0	:	
242	0036	000000	CALL:	.ACRD	0	:	
243	0040	000000	CPL:	.ACRD	0	:	
244	0042	000000	CALCNT:	.ACRD	0	:	
245	0044	001046	BPTIME:	.ACRD	550.	:	
246	0046	000000	CSPC:	.ACRD	0	:	
247	0050	000000	CSPC1:	.ACRD	0	:	
248	0052	000000	BPCNT1:	.ACRD	0	:	
249	0054	000000	CSP:	.ACRD	0	:	
250	0056	000000	CSP1:	.ACRD	0	:	
251	0060	000000	CCSP:	.ACRD	0	:	
252	0062	000000	CCSP1:	.ACRD	0	:	
253	0064	000012	HCNT:	.ACRD	10.	:	
254	0066	000000	SECS:	.ACRD	0.	:	
255	0070	000000	SECS30:	.ACRD	0	:	
256	0072	000000	LFLSHR:	.ACRD	0	:	
257	0074	000000	GASCVF:	.ACRD	0	:	
258	0076	000000	ANL:	.ACRD	0	:	
259	0100	000000	ANL1:	.ACRD	0	:	
260	0102	000000	ANLCNT:	.ACRD	0	:	
261						:	
262		000000		.CSECT	GASCCM		
263	0000		GASBUF:	.BLKW	2000.	:	GAS DATA BUFFER
264						:	
265		000000		.ASECT			
266		000104		.#104			
267	0104	000000		.ACRD	ADISB		KN11-P INTERRUPTS VECTOR TO ADISB
268	0106	000200		.ACRD	200		DISABLE ALL OTHER INTERRUPTS.
269		000001		.ENC			

ADISR RT-11 MACRO VM02-12 29-FEB-80 14:22:30 PAGE 14
SYMBOL TABLE

ADDH	000412R		ACIA	= 176772		ADISR	000000R	
ADSR	= 176770		ADR	000032R	004	AMRCNT	000022R	004
AMP1	000034R	004	AML	000076R	004	AMLCNT	000102R	004
AML1	000100R	004	BPCAT	000064R	004	BPCAT1	000052R	004
BPTIME	000044R	004	CCBP	00006CR	004	CCBP1	000062R	004
CNR	000026R	004	CHRCAT	00002CR	004	CHRI	000030R	004
CLPCSR	= 172540		CAT10	000002R	002	CNT30	000004R	002
COLNT	000000R	002	CSEF	000054R	004	CSBP1	000056R	004
CSPC	000046R	004	CSFC1	000050R	004	CHL	000036R	004
CHLCNT	000042R	004	CHL1	000040R	004	DAC1	= 176760	
DAC2	= 176762		CRIB	= 167774		DGUTB	= 167772	
EXM	000246R		EXHALE	000002R	004	EXMCNT	000012R	002
FINISH	000702R		FLICKR	000442R		FLRCNT	000008R	002
FLRFLG	000010R	002	FVC	000006R	004	FVCTST	000160R	
GASEUF	000000R	005	GASCSP	00000CR	004	GASQVF	000074R	004
GETCSP	000062R		MICAT	000722R		HIGH	000216R	
HIGH1	000152R		MRTPR8	000024R	004	LFLSHR	000072R	004
LITEON	000674R		LCCNT	000720R		LOOP0	000100R	
LOOP1	000252R		LCCF2	000264R		LCCP3	000276R	
LOOP4	000354R		LCCF5	000426R		LOOP6	000512R	
LOCP7	000540R		LCCF8	000562R		MASK	000010R	004
NOEXM	000200R		PC	=20000C7		PFT	000016R	004
PNEVSA	000012R	004	PLSRES	000014R	004	READSH	000306R	
R0	=2000000		R1	=20000C1		R2	=2000002	
R3	=2000003		R4	=20000C4		R5	=2000005	
SAMP0	000000R	003	SAMP1	000002R	003	SAMP2	000004R	003
SAMP3	000006R	003	SAMP4	000010R	003	SAMP5	000012R	003
SAMP6	000014R	003	SAMP7	000016R	003	SAMP8	000020R	003
SECS	000066R	004	SECS30	000070R	004	SP	=2000006	
STORSV	000166R		SVTPRS	000004R	004	TIMER	000602R	
...V2	= 000001							
.ABS.	000110	000						
	000724	001						
TIMCNT	000014	002						
ADCCM	000022	003						
ADCOM1	000104	004						
GASCOM	007640	005						

ERRORS DETECTED: 0
FREE CCRE: 17567, NCRDS

.LP:=ADISR

```

C-- EXERCIZE SUBROUTINE, VERSION 2.0
C--
C   WRITTEN BY: RCY A. REED
C   DATE:      OCT. 1978
C
C   REWRITTEN BY: WILLIAM G. CROSIER
C   DATE:      AUG. 24, 1979
C   MODIFIED BY: ABHIJIT GADGIL
C   DATE:      DEC. 1979
C
0001   SUBROUTINE EX
C--
C-- SET UP VARIABLES
C--
0002   LOGICAL*1 RETI(8), RETIX(8)
C-- !CONTAINES THE REAL TIME
0003   INTEGER EXHCNT, TIMBLF(6)
0004   INTEGER SECS, BPCNT, BPCNT1, BPTIME, SECS30
C-- !30 SECOND TIMER, FROM ADISR
0005   INTEGER*4 OUT
0006   INTEGER AGE
C-- !PATIENTS AGE, TO BE USED IN
C-- !THE CALCULATIONS OF THE PRECICTD VALUES.
0007   INTEGER CPID
C-- !PATIENTS CARDICFULMUNARY IDENTITY NUMBER.
0008   INTEGER NAME(14)
C-- !PATIENT'S NAME.
0009   INTEGER SSNO(3)
C-- !PATIENT'S SOCIAL SECURITY NUMBER.
0010   INTEGER DATE(3)
C-- !DATE THE TEST WAS CONCLCTED.
0011   INTEGER RETNG
C-- !PATIENT'S RETEST NUMBER.
0012   INTEGER UNIGNC(2)
C-- !PATIENT'S UNICLE NUMBER.
0013   INTEGER M,F,SEX
C-- !SUBJECTS SEX TO BE USED IN CALCLLATIONS.
0014   INTEGER ADCTEF(34)
C-- !COMMON AREA FOR A/D COMMAND VARIABLES
0015   INTEGER GASRUF(2000)
C-- !A/D COMMON DATA BLFFER.
0016   INTEGER GASDSP
C-- !GAS RUFFER DISPLACEMENT
0017   INTEGER EXHALE
0018   INTEGER SVTHRS
C-- !SPIROMETER VOLUME THRESHCLD.
0019   INTEGER FVC
C-- !FLAG (=1 FOR FORCED VITAL CAPACITY
C   !EQUAL TO 0 FOR ALL OTHERS).
0020   INTEGER PFT
C-- !FLAG =-1 FOR PULMUNARY FLACTION TESTS
C   !      =-0 FOR ALL OTHERS
0021   INTEGER HRTHRS
C-- !CONTAINES THE HEART RATE THRESHCLD.

```

```

C      !PUSHED.
0022      INTEGER PCSR
C--      !ADDRESS OF THE PROGRAMABLE CLOCK
C      !STATUS REGISTER.
0023      INTEGER PCCB
C--      !ADDRESS OF THE PROGRAMABLE CLOCK
C      !COUNTER BUFFER.
0024      INTEGER PCS
C--      !HAS THE VALUE TO BE PLACED INTO THE
C      !CLOCK STATUS REGISTER.
0025      INTEGER PCC
C--      !HAS THE VALUE TO BE PLACED INTO THE
C      !CLOCK COUNTER BUFFER.
0026      INTEGER OLD
C--      !VARIABLE TO HOLD THE OLD VALUE OF EXFAST.
0027      INTEGER DROUTB,CRINE
C--      !VARIABLES THAT CONTAIN THE ADDRESS OF THE
C--      !DW11 PARALLEL INTERFACE CARD INPUT AND OUTPUT BUFFERS.
0028      REAL HIGHT
C--      !PATIENT'S HEIGHT IN CM.
0029      REAL WT
0030      REAL SLPYCP(18)
C--      !VARIABLES THAT CONTAIN INFORMATION PASSED
C--      !THROUGH ADCOM1 COMMON BLOCK TO CALCULATE THE
C--      !SPIROMETER THRESHOLD (SLOPE & Y-INTERCEPT
C--      !ARRAY -- CALIB. FACTORS FOR ALL 9 CHANNELS
0031      REAL*8 SECTN,RST,EX,REC      !ALPHA LABELS FOR SECTION OF TEST
C--      SET UP THE COMMON AREAS.
C--
0032      COMMON/GAS30/VC2,VCC2,VCL,IERTH
0033      COMMON/TIMCNT/TIMBUF
0034      COMMON/PDATA/IDATA(1536)
C--      !SET UP THE COMMON BUFFER FOR PATIENT DATA
C--      !ALSO, CONTAINS THE PATIENTS INFORMATION.
0035      COMMON /GASCCN/GASBUF
C--      !DATA BUFFER COMMON AREA.
0036      COMMON /ADCOM1/ADCTFF
C--      !COMMON AREA TO PASS COMMAND VARIABLES AND
C      !CONSTANTS TO 8FROM A/D INTERRUPT SERVICE ROUTINE
0037      EQUIVALENCE (TIMBUF(6),EXFCNT)
C--      !CONTAINS THE END OF BREATH FLAG.
0038      EQUIVALENCE (IDATA(85),SLPYCP(1))
C--      ! CONTAINS THE SLOPE AND Y INTERCEPTS.
0039      EQUIVALENCE(IDATA(3),SSNC(1))
C--      !THESE THREE DATA WORDS CONTAIN
C--      !THE PATIENT'S SOCIAL SECURITY NUMBER.
0040      EQUIVALENCE (IDATA(7),UNIGNC)
C--      !THIS DATA WORD CONTAINS THE PATIENT'S
C--      !UNIQUE NUMBER.
0041      EQUIVALENCE (IDATA(9),RETNC)
C--      !THIS DATA WORD CONTAINS THE SUBJECT'S
C--      !RETEST NUMBER.
0042      EQUIVALENCE (IDATA(10),CPIC)
C--      !THIS DATA WORD CONTAINS THE PATIENT'S

```

```

C-- !CARDIOPULMONARY NUMBER.
0043 EQUIVALENCE(ICATA(40),DATE(1))
C-- !THESE THREE DATA WORDS CONTAIN THE
C-- !PATIENT'S TEST DATE.
0044 EQUIVALENCE(ICATA(26),NAME(1))
C-- !THE 14 DATA WORDS CONTAIN THE
C-- !PATIENT'S NAME.
0045 EQUIVALENCE(ICATA(2),SEX)
CC !THIS DATA WORD CONTAINS THE PATIENTS SEX.
0046 EQUIVALENCE(ICATA(25),AGE)
C-- !THIS DATA WORD CONTAINS THE PATIENTS AGE.
0047 EQUIVALENCE(ICATA(50),HEIGHT)
C-- !THIS DATA WORD CONTAINS THE PATIENTS HEIGHT.
0048 EQUIVALENCE(ICATA(52),WT)
0049 EQUIVALENCE(ICATA(72),KCNT)
0050 EQUIVALENCE(ACCTBF(19),BPTIME),(ACCTBF(22),BPCNT1)
0051 EQUIVALENCE(ACCTBF(27),BPCNT)
0052 EQUIVALENCE(ACCTBF(1),GASCSF)
C-- GAS BUFFER DISPLACEMENT
0053 EQUIVALENCE(ACCTBF(2),EXHALE)
0054 EQUIVALENCE(ACCTBF(3),SVTHRS)
C-- !SPIROMETER VOLUME THRESHOLD, MUST BE CALCULATED.
0055 EQUIVALENCE(ACCTBF(4),FVC)
C-- !FLAG EQUAL TO -1 FOR FORCED VITAL CAPACITY.
0056 EQUIVALENCE(ACCTBF(7),FLASHES)
C-- !CONTAINS 1'S FOR BUTTONS THAT ARE
C !PUSHED BY OPERATOR.
0057 EQUIVALENCE(ACCTBF(8),PFT)
C-- !FLAG FOR THE PFT TEST, MUST BE SET TO -1 FOR PFT'S
0058 EQUIVALENCE(ACCTBF(11),PRTHRS)
0059 EQUIVALENCE(ACCTBF(16),WL),(ACCTBF(20),AS)
0060 EQUIVALENCE(ACCTBF(29),SECS30),(ACCTBF(28),SECS)
C-- !CONTAINS THE ACISR 30 SEC. TIMER.
0061 EQUIVALENCE(ACCTBF(30),LFLSHR)
C-- !CONTAINS THE BITS OF THE LIGHTS TO FLASH.
C--
C-- INITIALIZE PARAMETERS
C--
0062 DATA RST, EX, REC/'* REST *', 'EXERCISE', 'RECOVERY'/
0063 DATA M/'M'/,F/'F'/
C-- !SET VARIABLES M AND F EQUAL TO THEIR ASCII VALUES.
0064 DATA DRCUTB/"167772/
C-- !DATA OUT BUFFER FOR THE PARALLEL INTERFACE CARD.
0065 DATA DRINB/"167774/
C-- !DATA IN BUFFER FOR THE PARALLEL INTERFACE CARD.
0066 DATA PCSR/"172540/
C-- !OCTAL ADDRESS OF THE PROGRAMABLE CLOCK(KW11-P)
C !STATUS REGISTER.
0067 DATA PCCB/"172542/
C-- !OCTAL ADDRESS OF THE PROGRAMABLE CLOCK
C !COUNTER BUFFER.
0068 DATA PCC/"144/
C-- !SET CLOCK TO GET 100 INTERRUPTS/SEC.
0069 DATA PCS/"113/

```



```

C-- !OCTAL VALUE TO BE LOADED INTO THE CLOCK STATUS
C !REGISTER. BITS 0,1,3, AND 6 ARE TURNED ON.
C ! BIT 6 - ALLOWS DCNE TO CAUSE AN INTERRUPT.
C ! BIT 3 - SELECTS REPEATED INTERRUPT MODE.
C ! BIT 1 - SET CLOCK RATE TO 10KHZ.
C ! BIT 0 - STARTS THE CLOCK.
0070 DATA OUT/'OUT '/'
0071 SVTHRS=IFIX((0.2-SLPYCP(2))/SLPYCP(1))
0072 HRTHRS=IFIX((35.-SLPYCP(10))/SLPYCP(9))
0073 FVC=0
C-- !INITIALIZE FVC FLAG.
0074 PFT=0
C-- !INITILIZE THE PFT FLAG.
D TYPE 5
C5 FORMAT(' BEFORE GLELE')
0075 IF(IGSET(5).NE.0) STOP 'GLELE ALLLOCATION FAILURE'
C--
0077 100 CONTINUE
C--
C-- SET UP THE REST PART
C--
0078 SECTN=RST
0079 CALL LIGHT('2)
0080 LITES='26
0081 ICT=61
0082 LFLSHR='2 !SET REST LIGHT TO FLASH
0083 IDATA(64)=0
0084 IDATA(68)=0
0085 KCNT=0 !INITILIZE RECCRD COUNTER
0086 ASSIGN 130 TO IRETHN !ASSIGN RETURN LABEL FOR START
0087 GO TO 520 !GO TO INITIALIZE & START ROUTINE.
0088 130 CONTINUE
0089 IF(EXHCNT.EQ.0)GO TO 135
D TYPE 3,GASDSP,EXHCNT
C3 FORMAT(2X,'GASDSP=',I3,'EXHCNT=',I3)
0091 CALL BREATH
D TYPE 4,V02,VCC2,VCL,IERTH
C4 FORMAT(2X,'V02=',F6.1,'VCC2=',F6.1,'VOL=',F6.1,'IERTH=',I3)
0092 135 CONTINUE
0093 IF(SECS30.EQ.0)GO TO 145
0095 WL=0.
0096 AS=0.
0097 CALL TRSHUF
0098 KCNT=KCNT+1
0099 CALL RPT30
0100 145 CONTINUE
0101 IRES=IPEEK(DRINE) !CHECK CONTROL PANEL PUSHBUTTON
0102 IF(IRES.EQ.'2)GC TO 100 !FESTART REST
0104 IF(IRES.EQ.'20)GO TO 160 !PAUSE
0106 IF(IRES.EQ.'4)GC TO 200 !START EXERCISE
0108 IF(KCNT.GE.20)GC TO 200 !MAX RECORDS REACHED
0110 GO TO 130 !GET NEXT 30.SEC RECORD
0111 160 ASSIGN 145 TO IWAIT1 !ASSIGN RETURN LABEL FOR PAUSE
0112 GO TO 562 !GO TO PAUSE ROUTINE

```

```

0113 200 CONTINUE
0114 CALL IPCKE(PCSR,0) !STOP CLOCK.
0115 IDATA(64)=(KCNT+1)/2
0116 CALL DATA (OUT)
0117 IF(IRES.EQ.4) GC TO 220
C WAIT UNTIL EX BLTCK IS PLSFD
0119 205 IF(IPEEK(DRINE).NE."4") GC TO 205
0121 207 IF(IPEEK(DRINE).NE.0) GC TO 207 !WAIT TILL EX RELEASED
C--
C-- SET UP EXERCIZE PART
0123 220 CALL IPCKE(PCSR,0)
0124 SECTN=EX
0125 CALL LIGHT("4")
0126 LITES="34
0127 LFLSHR="4
0128 ICT=65
C--
0129 IDATA(68)=0
0130 KCNT=0 !RECORD COUNT
0131 NREC=2*(60-IDATA(64))
0132 ILAST=0
0133 ASSIGN 230 TO IRETRN !ASSIGN RETURN LABEL FOR START ROUTINE
0134 GO TO 520 !GO TO INITIALIZE & START ROUTINE
0135 230 CONTINUE
0136 IF(SECS.EQ.ILAST)GC TO 233
0138 CALL WRKADJ
0139 ILAST=SECS
0140 233 CONTINUE
0141 IF(EXMCNT.EQ.0)GC TO 235
0143 CALL BREATH
0144 235 CONTINUE
0145 IF(SECS30.EQ.0)GC TO 245
0147 CALL TRSBUP
0148 KCNT=KCNT+1
0149 CALL RPT30
0150 245 CONTINUE
0151 IRES=IPEEK(DRINE) !CHECK USER RESPONSE
0152 IF(IRES.EQ."4")GC TO 220 !RESTART EXERCISE
0154 IF(IRES.EQ."20")GO TO 260 !PAUSE
0156 IF(IRES.EQ."10")GO TO 300 !START RECOVERY
0158 IF(KCNT.GE.NREC)GC TO 300 !EXERCISE DONE, MAX # OF RECORDS
0160 GO TO 230 !GET NEXT RECORD
0161 260 ASSIGN 245 TO IWAIT1 !ASSIGN RETURN LABEL FOR PAUSE ROUTINE
0162 GO TO 562 !GO TO PAUSE ROUTINE
0163 300 CONTINUE
0164 CALL IPCKE(PCSR,0)
0165 IDATA(68)=(KCNT+1)/2
0166 CALL DATA(OUT)
0167 IF(KCNT.GE.NREC)GC TO 900 !MAX REACHED, GO TO END
C--
C-- SET UP RECOVERY PART
C--
0169 310 CONTINUE
0170 CALL IPCKE(PCSR,0)

```

```

0171      MREC=2*(60-IDATA(64)-IDATA(68))
0172      KCNT=0
0173      SECTN=REC
0174      CALL LIGHT("10)
0175      LITES="70
0176      LFLSHR="10
0177      ICT=69
0178      ASSIGN 330 TO IRETRN      !ASSIGN RETURN LABEL FOR START ROUTINE
0179      GO TO 520      !GO TO INITIALIZE & START ROUTINE.
0180 330    CONTINUE
0181      IF(EXHCNT.EQ.0)GO TO 335
0183      CALL BREATH
0184 335    CONTINUE
0185      IF(SECS30.EQ.0)GO TO 345
0187      WL=0.
0188      AS=0.
0189      CALL IRSBUF
0190      KCNT=KCNT+1
0191      CALL RPT30
0192 345    CONTINUE
0193      IRES=IFEEK(DRINE)      !GET USER RESPONSE
0194      IF(IRES.EQ."10)GO TO 310      !RESTART RECOVERY
0196      IF(IRES.EQ."20)GO TO 360      !PAUSE
0198      IF(IRES.EQ."40)GO TO 400      !RECOVERY DONE (END OF TEST)
0200      IF(KCNT.GE.MREC) GO TO 400      !MAX RECORDS
0202      GO TO 330      !GET NEXT RECCRD
0203 360    ASSIGN 345 TO IWAIT1      !ASSIGN RETURN LABEL FOR PAUSE
0204      GO TO 562      !GO TO PAUSE ROUTINE
0205 400    CONTINUE
0206      CALL IPCKE(PCSR,0)      !TURN OFF CLOCK
0207 900    IDATA(72)=(KCNT+1)/2
           C      TYPE 1
           C1     FORMAT(' BEFORE LIGHT')
0208      CALL LIGHT(0)      !TURN OFF ALL LIGHTS
0209      RETURN

           C
           C      ROUTINE TO INITIALIZE & START EACH SECTION OF EXERCISE TEST
0210 520    IPAUSE=-1      !INITIALIZE PAUSE FLAG
0211      CALL TIMRD(IMR,IMIN,ISEC)
0212      IDATA(ICT)=IMR
0213      IDATA(ICT+1)=IMIN
0214      IDATA(ICT+2)=ISEC
0215 522    CONTINUE
0216      ISEC=IDATA(63)
0217      IF(IPAUSE.EQ.-1 .AND. SECTN.EQ.REC) GO TO 524
0219 523    CALL INITI
           C      TYPE 9,EXHCNT
           C9     FORMAT(2X,'EXHCNT=',I3)
0220      BPTIME=550 - (10*ISEC)
0221      IF(BPTIME.LT.0) BPTIME=BPTIME + 600
0223      V02=0.0
0224      VC02=0.0
0225      VOL=0.0
0226      IBRTH=0

```

```

      D      TYPE 10, V02, VCC2, VCL, IBRT+
      D10     FORMAT(2X, 3F7.1, I3)
0227  524    CALL IPCKE(PCCR, PCC)          !SET PROGRAMMABLE CLOCK TO 100 HZ
0228        CALL IPCKE(PCSR, PCS)          !SET STATUS IN CLOCK (START)
0229        CALL LIGHT(LITES)             !TURN ON APPROPRIATE PANEL LIGHTS
0230        IF(IPAUSE.NE.-1) GO TO 526
0232        IF(SECTN.EQ.RST) CALL RPT+CR
0234        TYPE 525, SECTN, IHR, IMIN, ISEC
0235  525    FORMAT(/5X, ' * * * ', AB, ' * * * '          REAL TIME: ', I2, ':',
@I2, ':', I2/)
0236  526    IF(IPAUSE.EQ.0) TYPE 527, JHR, JMIN, JSEC
0238  527    FORMAT(/' CONTINUE AT ', I2, ':', I2, ':', I2/)
0239  528    IF(IPEEK(DRINH).NE.0) GO TO 528      !WAIT TILL SWITCHES ARE RELEAS
0241        GO TO IRETRN                          !RETURN

      C
      C      PAUSE ROUTINE
      C
0242  562    IF(IPAUSE.EQ.1) GO TO 570      !ENDING PAUSE PERIOD?
0244        CALL IPCKE(PCSR, 0)              ! (START OF PAUSE) STOP CLOCK
0245        IPAUSE=1
0246        TYPE 564
0247  564    FORMAT(/' * PAUSE *')
0248        CALL LIGHT (LITES)                !TURN ON APPROPRIATE LIGHTS
0249  565    IF(IPEEK(DRINH).NE.0) GO TO 565 !WAIT TILL SWITCHES ARE RELEASED
0251        GO TO IWAIT1                      !GO WAIT FOR ANOTHER BUTTON PUS
0252  570    IPAUSE=0                        !CAME HERE IF ENDING PAUSE
0253        CALL TIMRD(JHR, JMIN, JSEC)
0254        ISEC=JSEC
0255        GO TO 523                          !END CF PAUSE ROUTINE
0256        END

```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
RETI	000014	LOGICAL*1 ARRAY (8)
RETX	000024	LOGICAL*1 ARRAY (8)
OUT	000104	INTEGER*4 VARIABLE
M	000064	INTEGER*2 VARIABLE
F	000066	INTEGER*2 VARIABLE
PCSR	000074	INTEGER*2 VARIABLE
PCCB	000076	INTEGER*2 VARIABLE
PCS	000102	INTEGER*2 VARIABLE
PCC	000100	INTEGER*2 VARIABLE
OLD	000314	INTEGER*2 VARIABLE
DROUTB	000070	INTEGER*2 VARIABLE
OWINB	000072	INTEGER*2 VARIABLE
SECTN	000316	REAL*8 VARIABLE
RST	000034	REAL*8 VARIABLE
EX	000044	REAL*8 VARIABLE
REC	000054	REAL*8 VARIABLE
IFIX	000000	INTEGER*2 PROCEDURE
IQSET	000000	INTEGER*2 PROCEDURE
LIGHT	000000	INTEGER*2 PROCEDURE
LITES	000326	INTEGER*2 VARIABLE
ICT	000330	INTEGER*2 VARIABLE
IRETRN	000332	INTEGER*2 VARIABLE
BREATH	000000	REAL*4 PROCEDURE
TRSHUF	000000	REAL*4 PROCEDURE
RPT30	000000	REAL*4 PROCEDURE
IRES	000334	INTEGER*2 VARIABLE
IPEEK	000000	INTEGER*2 PROCEDURE
IWAIT1	000336	INTEGER*2 VARIABLE
IPOKE	000000	INTEGER*2 PROCEDURE
DATA	000000	REAL*4 PROCEDURE
MREC	000340	INTEGER*2 VARIABLE
ILAST	000342	INTEGER*2 VARIABLE
WRKADJ	000000	REAL*4 PROCEDURE
IPAUSE	000344	INTEGER*2 VARIABLE
TIMRO	000000	REAL*4 PROCEDURE
IHR	000346	INTEGER*2 VARIABLE
IMIN	000350	INTEGER*2 VARIABLE
ISEC	000352	INTEGER*2 VARIABLE
INITI	000000	INTEGER*2 PROCEDURE
RPTHDR	000000	REAL*4 PROCEDURE
JHF	000354	INTEGER*2 VARIABLE
JMIN	000356	INTEGER*2 VARIABLE
JSEC	000360	INTEGER*2 VARIABLE

COMMON BLOCK /GAS30/ LENGTH 000016

V02	000000	REAL*4	VARIABLE
V002	000004	REAL*4	VARIABLE
V0L	000010	REAL*4	VARIABLE
IHRTH	000014	INTEGER*2	VARIABLE

FORTRAN IV STORAGE MAP

NAME OFFSET ATTRIBUTES

COMMON BLOCK /TIMCNT/ LENGTH 000014

TIMBUF 000000 INTEGER*2 ARRAY (6)

EXHCNT 000012 INTEGER*2 VARIABLE

COMMON BLOCK /PDATA/ LENGTH 006000

IDATA 000000 INTEGER*2 ARRAY (1536)

SLPYCP 000250 REAL*4 ARRAY (18)

SSNO 000004 INTEGER*2 ARRAY (3)

UNIGNO 000014 INTEGER*2 ARRAY (2)

RETNO 000020 INTEGER*2 VARIABLE

CPJD 000022 INTEGER*2 VARIABLE

DATE 000116 INTEGER*2 ARRAY (3)

NAME 000062 INTEGER*2 ARRAY (14)

SEX 000002 INTEGER*2 VARIABLE

AGE 000060 INTEGER*2 VARIABLE

HIGHT 000142 REAL*4 VARIABLE

WT 000146 REAL*4 VARIABLE

KCNT 000216 INTEGER*2 VARIABLE

COMMON BLOCK /GASCON/ LENGTH 007640

GASBUF 000000 INTEGER*2 ARRAY (2000)

COMMON BLOCK /ADCOM1/ LENGTH 000104

ADCTBF 000000 INTEGER*2 ARRAY (34)

HPTIME 000044 INTEGER*2 VARIABLE

BPCNT1 000052 INTEGER*2 VARIABLE

BPCNT 000064 INTEGER*2 VARIABLE

GASDSP 000000 INTEGER*2 VARIABLE

EXHALE 000002 INTEGER*2 VARIABLE

SVTHRS 000004 INTEGER*2 VARIABLE

FVC 000006 INTEGER*2 VARIABLE

PUSHES 000014 REAL*4 VARIABLE

PFT 000016 INTEGER*2 VARIABLE

HRTHRS 000024 INTEGER*2 VARIABLE

WL 000036 REAL*4 VARIABLE

AS 000046 REAL*4 VARIABLE

SECS30 000070 INTEGER*2 VARIABLE

SECS 000066 INTEGER*2 VARIABLE

LFLSHR 000072 INTEGER*2 VARIABLE

```

C
C-----
0001  SUBROUTINE WAKACJ
C-----
C  AUTHOR: CHARLES MANN      CONSULTANTS: PAUL SCHACHTER
C  DATE:   JULY 26, 1978      STAN FINK
C  PURPOSE: TO ADJUST THE WORK LOAD ON THE BICYCLE SO THAT A
C  SUBJECT'S HEART RATE MATCHES A GIVEN TARGET VALUE
C  AT A GIVEN TARGET TIME.
C  SUBROUTINE: SNATCH (ASSEMBLY LANGUAGE)
C-----
0002  INTEGER DAC1,CAC2,CLTPUT,IACDAT(34),ICALIB(1536)
0003  INTEGER*4 IHRAD,IWLAD
0004  REAL HRARY(30),WLARY(30)
0005  INTEGER IHRD(4)
0006  COMMON /ADCOM1/ IACCAT
0007  COMMON /PDATA/ ICALIB
0008  EQUIVALENCE (IACDAT(28),ICCLNT)      !ONE SECOND TIMER
0009  EQUIVALENCE (ICALIB(101),HRSLOP)
0010  EQUIVALENCE (ICALIB(103),HPRINT)
0011  EQUIVALENCE (ICALIB(105),WLSLOP)
0012  EQUIVALENCE (ICALIB(107),WLINT)
0013  EQUIVALENCE (ICALIB(15),IHRD(1))    !TARGET HEART RATES
0014  DATA WL/0./
0015  DATA CONST1/.25/,CONST2/0./
0016  DATA CAC1/"176760"/,CAC2/"176762"/
0017  DATA ISTDND/0/,IHRDFT/1/,HRS/80./
0018  DATA HRARY/30*80./
C  CHECK FOR TREADMILL--THIS ROUTINE WORKS ONLY FOR THE BIKE
0019  IF (ICALIB(12).NE.1) RETURN
C  SHIFT HR INTO ARRAY
DO 10 I=1,29
0021  WLARY(I)=WLARY(I+1)
0022  HRARY(I)=HRARY(I+1)
0023  10  CONTINUE
C  GET THE HR & WL DATA STORED BY THE ACISR
0025  CALL SNATCH(IHRAD,IWLAD,IHRCNT,IWLCNT)
C  CHANGE A/D HEART RATE INTO BPM
0026  IF((IHRCNT.EQ.0).OR.(IWLCNT.EQ.0)) GO TO 13
0028  HRARY(30)=HRSLOP*(AJFLT(IHRAD)/FLCAT(IHRCNT))+HPRINT
0029  WLARY(30)=WLSLOP*(AJFLT(IWLAD)/FLCAT(IWLCNT))+WLINT
0030  ACTWL=WLARY(30)
0031  ACTHR=HRARY(30)
0032  13  IF (MOD(ICOUNT,5).NE.0) GO TO 999
0034  IF(ICCLNT.EQ.5) GO TO 110
0036  IF(ICOUNT.EQ.240) GO TO 120
0038  IF(ICCLNT.EQ.480) GO TO 120
0040  IF(ICOUNT.EQ.720) GO TO 120
0042  IF((240.LT.ICCLNT).AND.(ICCLNT.LT.300)) GO TO 110
0044  IF((480.LT.ICCLNT).AND.(ICCLNT.LT.540)) GO TO 110
0046  IF((720.LT.ICCLNT).AND.(ICCLNT.LT.780)) GO TO 110
0048  IF (ICCLNT.GT.960) GO TO 103
C  STANDARD FLOATING WL ROUTINE
0050  R0  HPSOLD=HRS

```

```

      C CALCULATE HRS(SMCCTH H R AVG.)
0051      DO 60 I=26,30
0052 60      HRS=HRS+HRARY(I)
0053      HRS=HRS/6.
      C CHECK TO SEE IF H R IS CLOSE ENOUGH
0054      IF (ABS(FLOAT(IHRC(IHRCPT))-HRS).LE.1.0) GO TO 999
      C DON'T LET H R CHANGE TOO QUICKLY
0056      IF (ABS(HRS-HRSCLD).GT.5.0) GO TO 999
      C CALCULATE NEW WL
0058      WLOLD=WL
0059      WLOLD=WL
      ISAVE CLD WL
      WL=WLOLD+CONST1*(FLCAT(IHRC(IHRCPT))-HRS)+CONST2*((HRS-HRSCLD)
      & /ABS(FLOAT(IHRC(IHRCPT))-HRS))
0060      GO TO 140
      C DETERMINE THE RESTING HEART RATE
0061 110      WLFREE=0.
      C CALCULATE THE AVERAGE RESTING H R FROM PRE-EXERCISE RECORDS
0062      IREST=ICALIB(64)
0063      IR64=MCD(IREST,2)
0064      IRECNO=(IREST/2)+IR64
0065      IF (IRECNO.GT.0) GO TO 111
0067      HRFREE=65.
0068      GO TO 113
0069 111      HRFREE=0.
0070      DO 112 I=1,IRECNO
0071 112      HRFREE=HRFREE+ICALIB(258+(I-1)*20)
0072      HRFREE=HRFREE/IRECNO
0073      IF (HRFREE.LT.45.) HRFREE=45.
0075      IF (HRFREE.GT.70.) HRFREE=70.
      C SET HRD
0077 113      IHRCPT=1
      C INITIALIZE WL
0078      WL=0.
      C GO TO STANDARD WL ROUTINE
0079      GO TO 80
      C END OF EXERCISE--DECREASE WL TO 10% OF END IN 90 SEC
0080 103      WL=WL*.85
0081      GO TO 140
0082 117      ISTDND=1
0083 118      IF (ISTDND.EG.1) GO TO 80
0085      GO TO 130
      C SET NEW DESIRED H R
0086 120      IHRCPT=IHRCPT+1
      C INITIALIZE FLAG TO HOLD A CONSTANT WL
0087      ISTDND=0
      C SAVE OLD WL & HR
0088      WLPAST=WLFREE
0089      HRPAST=HRFREE
      C CALCULATE THE CONSTANT WL VALUE
0090      WLFREE=0.
0091      HRFREE=0.
0092      DO 125 I=1,30
0093      WLFREE=WLFREE+WLARY(I)
0094 125      HRFREE=HRFREE+HRARY(I)
0095      HRFREE=HRFREE/30.

```



```

0096      WLFREE=WLFREE/30.
0097      SLWL=(WLFREE-WLFAST)/(HREFREE-HRPAST)
0098      WLCNST=WLFREE
0099      IF(SLWL.LE.0) GO TO 117
C   IF TARGET HR IS LESS THAN THE CURRENT AVE. HR OR CLOSE TO IT,
C   GO TO STANDARD CORRECTOR WL RATHER THAN A CONSTANT WL.
0101      IF ((IFIX(HREFREE)-IHRD(IHRCPT)).GT.-5) GO TO 117
0103      WLCNST=WLFREE+SLWL*(FLOAT(IHRD(IHRCPT))-HREFREE)
C   KEEP WORK LOAD CONSTANT
0104      130  WL=WLCNST
C   OUTPUT THE WL VALUE TO THE RICYCLE
C   CHECK FOR MAX WORKLOAD
0105      140  IF (WL.GT.300.) WL=300.
C   CHECK FOR MIN WORKLOAD
0107      IF (WL.LT.0) WL=0.
0109      OUTPUT=IFIX(6.144*WL)
0110      CALL IPCKE(DAC1,OUTPUT)
0111      CALL IPCKE(DAC2,OUTPUT)
0112      999  CONTINUE
CC      IF (MOD(ICOUNT,15).NE.0) RETURN
CC      TYPE 825,ICOUNT,ACTHR,HRS,IHRD(IHRCPT),WL,ACTWL,WLFREE,HREFREE
CC825  FORMAT(' TIME=',I4,3X,'IHR=',F4.0,3X,'HRS=',F4.0,3X,'HPRD=',I4,3X,
CC      *      'CAL WL=',F4.0,3X,'ACT WL=',F4.0,3X
CC      *      ',WLF=',F4.0,3X,'HREF=',F4.0)
0113      RETURN
0114      END

```

IBACKUP CHANNEL

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
HRARY	000014	REAL*4 ARRAY (30)
WLARY	000204	REAL*4 ARRAY (30)
DAC1	000410	INTEGER*2 VARIABLE
DAC2	000412	INTEGER*2 VARIABLE
OUTPUT	000440	INTEGER*2 VARIABLE
IHRAD	000442	INTEGER*4 VARIABLE
IWLAD	000446	INTEGER*4 VARIABLE
WL	000374	REAL*4 VARIABLE
CONST1	000400	REAL*4 VARIABLE
CONST2	000404	REAL*4 VARIABLE
ISTAND	000414	INTEGER*2 VARIABLE
IHRDPT	000416	INTEGER*2 VARIABLE
HRS	000420	REAL*4 VARIABLE
I	000452	INTEGER*2 VARIABLE
SNATCH	000000	REAL*4 PROCEDURE
IHRCNT	000454	INTEGER*2 VARIABLE
IWLCNT	000456	INTEGER*2 VARIABLE
AJFLT	000000	REAL*4 PROCEDURE
FLOAT	000000	REAL*4 PROCEDURE
ACTWL	000460	REAL*4 VARIABLE
ACTHR	000464	REAL*4 VARIABLE
MOD	000000	INTEGER*2 PROCEDURE
HRSOLD	000470	REAL*4 VARIABLE
ABS	000000	REAL*4 PROCEDURE
WLOLD	000474	REAL*4 VARIABLE
WLFREE	000500	REAL*4 VARIABLE
IREST	000504	INTEGER*2 VARIABLE
IR64	000506	INTEGER*2 VARIABLE
IRECNO	000510	INTEGER*2 VARIABLE
HRFREE	000512	REAL*4 VARIABLE
WLPAST	000516	REAL*4 VARIABLE
HRPAST	000522	REAL*4 VARIABLE
SLWL	000526	REAL*4 VARIABLE
WLCNST	000532	REAL*4 VARIABLE
IFIX	000000	INTEGER*2 PROCEDURE
IPOKE	000000	INTEGER*2 PROCEDURE

COMMON BLOCK /ADCON1/ LENGTH 000104

IADDAT	000000	INTEGER*2 ARRAY (34)
ICGUNT	000066	INTEGER*2 VARIABLE

COMMON BLOCK /PODATA/ LENGTH 006000

ICALIB	000000	INTEGER*2 ARRAY (1536)
HRSLOP	000310	REAL*4 VARIABLE
HRINT	000314	REAL*4 VARIABLE
WLSLOP	000320	REAL*4 VARIABLE
WLINT	000324	REAL*4 VARIABLE
IHRD	000034	INTEGER*2 ARRAY (4)

```

1 .....
2 .TITLE SNATCH
3 .....
4
5 .AUTHOR: CHARLES MARR
6 .DATE: SEPTEMBER 21, 1978
7 .PURPOSE: TO GET THE HR & WL DATA STORED IN COMMON BY THE INTERRUPT
8 .SERVICE ROUTINE AND TO INITIALIZE THESE LOCATIONS BACK TO
9 .ZERO.
10 .....
11
12 .GLCCL SNATCH
13 .PCALL .REGDEF,..V2..
14 ..V2..
15 .REGDEF
16 00000 106427 SNATCH: MTPS #340 ;INHIBIT INTERRUPTS
17 00000 000340
18 00004 016501 MCV 2(R5),R1 ;STORE ADDRESS OF INRAD INTO R1
19 00000 000002
20 00010 016721 MCV CHR,(R1)+ ;GET LOWER ORDER WORD OF HR SUM
21 00000 000020
22 00014 016711 MCV CHR1,(R1) ;GET HIGHER ORDER WORD OF HR SUM
23 00000 000030
24 00020 016501 MCV 4(R5),R1 ;STORE ADDRESS OF INLAD INTO R1
25 00000 000004
26 00024 016721 MCV CHL,(R1)+ ;GET LOWER ORDER WORD OF WL SUM
27 00000 000076
28 00030 016711 MCV CHL1,(R1) ;GET HIGHER ORDER WORD OF WL SUM
29 00000 000100
30 00034 016775 MCV CHRCNT,26(R5) ;GET NO. OF SAMPLES FOR HR
31 00000 000020
32 00006 000006
33 00042 016775 MCV CHLCNT,210(R5) ;GET NO. OF SAMPLES FOR WL
34 00000 000102
35 00010 000010
36 00050 005067 CLR CHR ;ZERO THE WORDS IN PREPARATION
37 00000 000026
38 00054 005067 CLR CHR1 ; FOR MORE DATA
39 00000 000030
40 00060 005067 CLR CHL
41 00000 000076
42 00064 005067 CLR CHL1
43 00000 000100
44 00070 005067 CLR CHRCNT
45 00000 000020
46 00074 005067 CLR CHLCNT
47 00000 000102
48 00100 106427 MTPS #0 ;ALLOW ADISR TO COLLECT MORE DATA
49 00000 000000
50 00104 000207 RTS PC ;RETURN TO WRKADJ ROUTINE
51 00000 000000 .CSECT ADCCM1 ;COMMON AREA WHERE ADISR STORES A/D DATA
52
53 .ELKW 0.
54 00020 000000 CHRCNT: .ACRD
55 .ELKW 2
56 00026 000000 CHR: .ACRD
57 00030 000000 CHR1: .ACRD
58 .ELKW 18.
59 00076 000000 CHL: .ACRD

```

SNATCH RT-11 MACRO VM02-12 29-FEB-80 15:07:10 PAGE 1+

40 00100 000000 CALL: .ACRO
41 00102 000000 CALCNT: .ACRO
42 000001 .END

SNATCH RT-11 MACRO VM02-12 29-FEB-80 15:07:10 PAGE 1+
SYMBOL TABLE

CHR	000026R	002	CHRCNT	00002CR	002	CHR1	000030R	002
CHL	000076R	002	CHLCNT	000102R	002	CHL1	000100R	002
PC	=2000007		R0	=2000000		R1	=2000001	
R2	=2000002		R3	=2000003		R4	=2000004	
R5	=2000005		SNATCH	00000CRG		SP	=2000006	

...V2 = 000001

.ABS. 000000 000

000106 001

ADCOM1 000104 002

ERRORS DETECTED: 0

FREE CCRE: 17884. ACROS

.LP:=SNATCH

```

1          .GLCBL TRSBUF
2          .NCALL .REGDEF,...V2..
3 000000   ..V2..
4 000000   .REGDEF
5 000000   106427 TRSBUF: MTPS      #200
           000200
6 000004   010146   MCV      R1,-(SP)
7 000006   010246   MCV      R2,-(SP)
8 000010   010346   MCV      R3,-(SP)
9 000012   012701   MCV      #AHRCNT,R1
           000022
10 00016   012702   MCV      #AHRC,R2
           000000
11 00022   012703   MCV      #20.,R3
           000024
12 00026   012122 LOOP:  MCV      (R1)+,(R2)+
13 00030   005303   DEC      R3
14 00032   003375   BGT      LOOP
15 00034   012701   MCV      #V02,R1
           000000
16 00040   012702   MCV      #V02C,R2
           000000
17 00044   012703   MCV      #7.,R3
           000007
18 00050   012122 LOOP0: MCV      (R1)+,(R2)+
19 00052   005303   DEC      R3
20 00054   003375   BGT      LOOP0
21 00056   012701   MCV      #AHRCNT,R1
           000022
22 00062   012703   MCV      #20.,R3
           000024
23 00066   005021 LOOP1: CLR      (R1)+
24 00070   005303   DEC      R3
25 00072   003375   BGT      LOOP1
26 00074   016767   MCV      AHRTC,AHRT
           000002
           000024
27 00102   016767   MCV      CHRC,CHR
           000004
           000026
28 00110   016767   MCV      CHR1C,CHR1
           000006
           000030
29 00116   016767   MCV      BPTC,BPTIME
           000022
           000044
30 00124   016767   MCV      BPCT1C,BPCNT1
           000030
           000052
31 00132   012767   MCV      #10.,BPCNT
           000012
           000064
32 00140   016767   MCV      SECSC,SECS
           000044
           000066
33 00146   012703   MCV      #7.,R3
           000007

```

```

34 00152 012701      MCV      #V02,R1
      000000'
35 00156 005021 LOOP2: CLR      (R1)+
36 00160 005303      DEC      R3
37 00162 003375      BGT      LOOP2
38 00164 012603      MCV      (SP)+,R3
39 00166 012602      MCV      (SP)+,R2
40 00170 012601      MCV      (SP)+,R1
41 00172 106427      MTPS     #0
      000000
42 00176 000207      RTS      PC
43
44      000000'      .CSECT  ADCCM1
45      .BLKW      9.
46 00022 000000 AHRCNT: .WORD   0
47 00024 000000 AHRT:  .WORD   0
48 00026 000000 CHR:   .WORD   0
49 00030 000000 CHR1:  .WORD   0
50      .BLKW      5.
51 00044 001046 BPTIME: .WORD  550.
52      .BLKW      2.
53 00052 000000 BPCNT1: .WORD   0.
54      .BLKW      4.
55 00064 000012 BPCNT:  .WORD  10.
56 00066 000000 SECS:   .WORD   0
57      .BLKW      6
58      000000'      .CSECT  GAS30
59 00000      V02:      .BLKW   7.
60      000000'      .CSECT  GAS30C
61 00000      V02C:     .BLKW   7.
62      000000'      .CSECT  ADCCM2
63 00000 000000 AHRC:   .WORD   0.
64 00002 000000 AHRTC:  .WORD   0
65 00004 000000 CHRC:   .WORD   0
66 00006 000000 CHR1C: .WORD   0
67      .BLKW      5.
68 00022 000000 BPTC:   .WORD   0
69      .BLKW      2.
70 00030 000000 BPCT1C: .WORD   0
71      .BLKW      5.
72 00044 000000 SECSC:  .WORD   0
73      .BLKW      1
74      000001'      .ENC

```

TRSBUF

MAIN: RT-11 MACRO VM02-12 29-FEB-80 15:10:48 PAGE 1+
SYMBOL TABLE

AMRC	000000R	005	AMRCNT	000022R	002	AHRT	000024R	002
AHRTC	000002R	005	BPCAT	000064R	002	BPCNT1	000052R	002
BPCTIC	000030R	005	BPTC	000022R	005	BPTIME	000044R	002
CHR	000026R	002	CHRC	0000C4R	005	CHR1	000030R	002
CHR1C	000006R	005	LCCF	000026R		LCOP0	000050R	
LOOP1	000066R		LCOP2	000156R		PC	=X000007	
R0	=X000000		R1	=X000001		R2	=X000002	
R3	=X000003		R4	=X000004		R5	=X000005	
SECS	000066R	002	SECSC	000044R	005	SP	=X000006	
TRSBUF	000000RG		VC2	000000R	003	VO2C	000000R	004

...V2 = 000001
. ABS. 000000 000
000200 001
ADCOM1 000104 002
GAS30 000016 003
GAS30C 000016 004
ADCOM2 000050 005

ERRORS DETECTED: 0

FREE CORE: 17813. WORDS

,LP:=TRSBUF

0001

SUBROUTINE BREATH

```

C-----
C
C   VERSION: 2.0
C   AUTHOR:  CHARLES MARN
C   DATE:    SEPTEMBER 18, 1978
C   MODIFIED BY:      WILLIAM G. CROSIER
C   DATE:          10 SEP 79
C   PURPOSE: TO DETERMINE THE OXYGEN CONSUMPTION AND THE CARBON
C             DIOXIDE PRODUCTION IN EACH BREATH.
C-----

```

```

0002   INTEGER GASBUF(2000),CALIB(1536),TIME(6)
0003   REAL N2,N2SLP,N2INT,N2OUT,N2AMB
0004   COMMON /GASCON/ GASBUF
0005   COMMON /PDATA/ CALIB
0006   COMMON /GAS30/ VO2,VCO2,VCL,IBRTH
0007   COMMON /TINCNT/ TIME
0008   COMMON /ADCOM1/ IECPRT,ICLNMY(33) !END OF BREATH POINTER IN GASBUF
0009   EQUIVALENCE (CALIB(73),STPD)
0010   EQUIVALENCE (CALIB(75),BTFS)
0011   EQUIVALENCE (CALIB(77),C2AMB)
0012   EQUIVALENCE (CALIB(79),N2AMB)
0013   EQUIVALENCE (CALIB(81),CC2AMB)
0014   EQUIVALENCE (CALIB(85),VCLSLP)
0015   EQUIVALENCE (CALIB(87),VCLINT)
0016   EQUIVALENCE (CALIB(89),G2SLP)
0017   EQUIVALENCE (CALIB(91),C2INT)
0018   EQUIVALENCE (CALIB(93),N2SLP)
0019   EQUIVALENCE (CALIB(95),N2INT)
0020   EQUIVALENCE (CALIB(97),CC2SLP)
0021   EQUIVALENCE (CALIB(99),CC2INT)
0022   EQUIVALENCE (TIME(6),IECRFG) !END OF BREATH FLAG: COUNTS BREATH
C   STORE SUBSCRIPT OF THE END OF THE BREATH SAMPLE
0023   S   IEORP=IEORPT
0024   IPOINT=IEORP-3
0025   IF(IPOINT.LT.1) IPOINT=1
CD   TYPE 1,(GASBUF(I),I=1,IECPRT)
0027   1   FORMAT(/,1X,8(1X,I5))
C   DECREMENT A/D BREATH COUNTER TO SHOW THAT THIS ONE IS BEING PROCESSED
0028   6   IEORFG=IEORFG-1
0029   IEORF=IEORFG
C   INCREMENT THE BREATH COUNTER FOR THE 30 SEC REPORT
0030   IBRTH=IBRTH+1
C   LOCATE MAXIMUM SPIROMETER VOLUME
0031   JPOINT=IPOINT-124
0032   IF(JPOINT.LT.1) JPOINT=1
0034   MAXVOL=GASBUF(JPOINT)
0035   IEND=JPOINT
0036   DO 10 I=JPOINT,IPOINT,4
0037   IF (GASBUF(I).LT.MAXVOL) GO TO 10
0039   MAXVOL=GASBUF(I)
0040   IEND=I !SUBSCRIPT OF MAX VOLUME

```



```

0041 10 CONTINUE
      C
      C AVERAGE GAS CONCENTRATIONS AT LAST 5 PCINTS CF BREATH TO REMCUE NOISE
0042      ISTART=IEND-16
0043      IF(ISTART.LT.1)ISTART=1
0044      IO2=0
0045      IN2=0
0046      ICO2=0
0047      DO 20 I=ISTART,IEND,4
0048          ICO2=ICO2+GASEUF(I+3)
0049          IO2=IO2+GASEUF(I+1)
0050          IN2=IN2+GASEUF(I+2)
0051 20 CONTINUE
0052      IDIF=((IEND-ISTART)+4)/4
0053      IF(IDIF.LT.1)ICIF=1
0054      IO2=IO2/IDIF
0055      IN2=IN2/IDIF
0056      ICO2=ICO2/IDIF
      C
      C AVERAGE LAST 3 SPIRC. VC. SAMPLES TO GET TIDAL VOLUME
0057      IF (IEND.GE.13) MAXVOL=(GASEUF(IEND-4)+GASEUF(IEND-8)
0058          +GASEUF(IEND-12))/3
      C
      C USE CALIBRATION FACTORS TO CONVERT A/D COUNTS INTO PERCENTAGES
0059      O2=O2SLP*FLOAT(IO2)+O2INT
0060      N2=N2SLP*FLCAT(IN2)+N2INT
0061      CO2=CO2SLP*FLCAT(ICO2)+CC2INT
      C
      C USE CALIBRATION FACTORS TO CONVERT A/D COUNTS INTO LITERS
0062      VOLOUT=VOLSLP*FLOAT(MAXVOL)+VOLINT
      C
      C INCREMENT MINUTE VOLUME
0063      VOL=VOL+VOLOUT*ETPS
      C
      C TYPE 100,IBRTH,C2,N2,CO2,VOL
0064 100 FORMAT(/,1X,'IBRTH  O2  N2  CO2  VOL ',I4,4(2X,FR.2))
      C
      C NORMALIZE MASS SPEC. DATA TO CORRECT FOR DRIFT
0065      ALL=O2+N2+CO2
0066      O2OUT=O2/ALL
0067      N2OUT=N2/ALL
0068      CO2OUT=CO2/ALL
      C
      C INCREMENT THE O2 CONSLMPTION AND CO2 PRODUCTION
0069      STPVOL = STPD * VOLCUT
0070      VO2 = VC2 + STPVOL * ((O2AMB/N2AMB)*N2CUT - O2OUT)
0071      VCO2 = VCO2 + STPVOL * (CC2CUT - (CC2AMB/N2AMB)*N2OUT)
      C
      C CHECK TO SEE IF THERE IS ANOTHER BREATH TO DEAL WITH
0072      IF (IEQBFG.EG.0) GO TO 900  ! NO MORE BREATHS, RETURN
0073      IF(IECBF.EG.0) GO TO 200  ! NO PRIOR BREATHS
      C
      C FIND END OF PREVIOUS BREATH
0074      IREF=(0.3-VOLINT)/VCLSLP
0075      IPOINT=IEND
0076      IPOINT=IPOINT-4
0077 150 IF (IPOINT.LT.41) GO TO 160 !BREATH .LT. 0.1 SEC. LONG
0078      IF(GASEUF(IPCINT).GT.IREF) GO TO 150
0079      GO TO 6

```

```
0086 160 IE0BFG=IE0BFG-IE0BF      !COULDN'T FIND PREVIOUS EOB(S)
0087 200 IF (IE0BPT.GT. IE0BP) GC TC 5      ! A NEW BREATH IS FINISHED
      C RESET POINTER TO THE BEGINNING OF THE GAS BUFFER
0089 900 IE0BPT=0
0090      RETURN
0091      END
```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
N2	000102	REAL*4 VARIABLE
N2OUT	000106	REAL*4 VARIABLE
IEOBP	000112	INTEGER*2 VARIABLE
IPOINT	000114	INTEGER*2 VARIABLE
IEORF	000116	INTEGER*2 VARIABLE
JPOINT	000120	INTEGER*2 VARIABLE
MAXVOL	000122	INTEGER*2 VARIABLE
IEND	000124	INTEGER*2 VARIABLE
I	000126	INTEGER*2 VARIABLE
ISTART	000130	INTEGER*2 VARIABLE
I02	000132	INTEGER*2 VARIABLE
IN2	000134	INTEGER*2 VARIABLE
IC02	000136	INTEGER*2 VARIABLE
IDIF	000140	INTEGER*2 VARIABLE
O2	000142	REAL*4 VARIABLE
FLOAT	000000	REAL*4 PROCEDURE
C02	000146	REAL*4 VARIABLE
VOLCUT	000152	REAL*4 VARIABLE
ALL	000156	REAL*4 VARIABLE
O2OUT	000162	REAL*4 VARIABLE
C02OUT	000166	REAL*4 VARIABLE
STPVOL	000172	REAL*4 VARIABLE
IREF	000176	INTEGER*2 VARIABLE

COMMON BLOCK /GASCCN/ LENGTH 007640

GASRUF 000000 INTEGER*2 ARRAY (2000)

COMMON BLOCK /PDATA/ LENGTH 006000

CALIB	000000	INTEGER*2 ARRAY (1536)
STPD	000220	REAL*4 VARIABLE
BTPS	000224	REAL*4 VARIABLE
O2AMB	000230	REAL*4 VARIABLE
N2AMB	000234	REAL*4 VARIABLE
CO2AMB	000240	REAL*4 VARIABLE
VOLSLP	000250	REAL*4 VARIABLE
VOLINT	000254	REAL*4 VARIABLE
O2SLP	000260	REAL*4 VARIABLE
O2INT	000264	REAL*4 VARIABLE
N2SLP	000270	REAL*4 VARIABLE
N2INT	000274	REAL*4 VARIABLE
CO2SLP	000300	REAL*4 VARIABLE
CO2INT	000304	REAL*4 VARIABLE

COMMON BLOCK /GAS30/ LENGTH 000016

VO2	000000	REAL*4 VARIABLE
VC02	000004	REAL*4 VARIABLE
VOL	000010	REAL*4 VARIABLE
IBWTH	000014	INTEGER*2 VARIABLE

FORTRAN IV STORAGE MAP

NAME OFFSET ATTRIBUTES

COMMON BLOCK /TIMCNT/ LENGTH 000014

TIME 000000 INTEGER*2 ARRAY (6)

LEOBFG 000012 INTEGER*2 VARIABLE

COMMON BLOCK /ADCOM1/ LENGTH 000104

LEOBPT 000000 INTEGER*2 VARIABLE

IDUMMY 000002 INTEGER*2 ARRAY (33)

C--
CC-- SUBROUTINE TO PRINT THE 30 SEC. REPCRT
C--
C-- WRITTEN BY: ROY A. REED
C-- DATE: 13-OCT-78
C--

```
0001 SUBROUTINE RPT30
0002 INTEGER NUM, INLM, IBLK
0003 INTEGER DATBUF(20)
0004 INTEGER IBRTH, CWLCNT, AHRCNT, BPCNT, BPCNT1
0005 INTEGER*4 CSPC, AHR, CWL, CSBP, CDBP
0006 REAL DDATA(640)
0007 REAL VC2, VCO2, VCL, RDATA(24)
0008 REAL HR, WL, AS, SBP, CBP
0009 COMMON /GAS30C/ VC2, VCO2, VCL, IBRTH
0010 COMMON /PDATA/ ICATA(1536)
0011 COMMON /ADCOM2/ CATBLF
0012 EQUIVALENCE (CATBUF(5), AHR), (DATBUF(7), CWL),
X(DATBUF(9), CWLCNT), (DATBUF(11), CSPC),
X(DATBUF(13), BPCNT1), (DATBUF(14), CSBP),
X(DATBUF(16), CDBP), (CATBUF(18), BPCNT),
X(DATBUF(1), AHRCNT)
0013 EQUIVALENCE (ICATA(257), CCATA(1))
0014 EQUIVALENCE (ICATA(72), NLM)
0015 EQUIVALENCE (ICATA(73), PCATA(1)), (ICATA(52), WT)
0016 DATA IBLK/2H /
0017 DATA INUM/2H /
0018 IF (IAJFLT(AHR, HR).EQ.-2) STOP
0020 IF (IAJFLT(CWL, WL).EQ.-2) STOP
0022 IF (IAJFLT(CSPC, AS).EQ.-2) STOP
0024 IF (IAJFLT(CSBP, SBP).EQ.-2) STOP
0026 IF (IAJFLT(CDBP, CBP).EQ.-2) STOP
0028 IF (HR.NE.0.) HR=(HR/FLCAT(AHRCNT))*RDATA(15)+RDATA(16)
0030 IF (WL.NE.0.) WL=(WL/FLCAT(CWLCNT))*RDATA(17)+RDATA(18)
0032 IF (AS.NE.0.) AS=(AS/FLCAT(CWLCNT))*RDATA(19)+RDATA(20)
0034 IF (SBP.NE.0.) SBP=(SBP/FLCAT(BPCNT))*RDATA(21)+RDATA(22)
0036 IF (DBP.NE.0.) CBP=(CBP/FLCAT(BPCNT))*RDATA(23)+RDATA(24)
0038 IF (NUM.EQ.1) TKL=0.
0040 IF (HR.LT.0.) HR=0.
0042 IF (WL.LT.0.) WL=0.
0044 IF (AS.LT.0.) AS=0.
0046 IF (SBP.LT.35.) SBP=0.
0048 IF (DBP.LT.35.) CBP=0.
0050 V02LM=0.
0051 VCO2LM=0.
0052 V02MKM=0.
0053 TVOL=0.
0054 RER=0.
0055 IRR=0
0056 RMV02=0.
0057 RMVCO2=0.
0058 V02HRR=0.
0059 IREM=NOC(NUM, 2)
0060 NNUM=NUM/2
```

```

0061      IF(IREM.EQ.0)ENCODE(2,51,INUM) NALM
0063  51    FORMAT(I2)
0064      VO2LM=VC2*2.
0065      VO2MKM=(VO2LM*1000.)/WT
0066      VCO2LM=VCO2*2.
0067      TVOL=VOL*2.
0068      IF(VO2LM.GT.0.)RER=VCO2LM/VC2LM
0070      IRR=IBRM*2
0071      IF(VO2LM.GT.0.)RMVC2=TVOL/VC2LM
0073      IF(VCO2LM.GT.0.)RMVCC2=TVCL/VCO2LM
0075      IF(HR.NE.0.)VC2+RB=(VO2LM*1000.)/+R
0077      IHR=IFIX(HR+.5)
0078      ISBP=IFIX(SBP+.5)
0079      IDBP=IFIX(DBP+.5)
0080      IF(IDATA(12).EQ.2)GC TO 200
0082      TWL=TWL+(WL/2.)
0083      ITWL=IFIX(TWL+.5)
0084      IWL=IFIX(WL+.5)
0085      GO TO 300
0086  200    CONTINUE
0087      IWL=IFIX(WL*10.+.5)
0088      TWL=AS
0089      ITWL=IFIX(AS+.5)
0090  300    TYPE 400,INUM,I+R,ISBP,ICBP,WL,TWL,VO2LM,VO2MKM,VCO2LM,
0091  400    XTVOL,RER,IRR,RMVC2,RMVCC2,VC2+RB
0092      IREST=ICATA(64)
0093      IEX=IDATA(68)
0094      NPT=IREST+IEX+(NUM-1)/2
0095      K=257+NPT*20
0096      KR=4+NPT*10
0097      KK=2-IREM
0098  CD     TYPE 450,IDATA(64),IDATA(68),NUM,NPT,K,KR,KK,IREM
0098  450    FORMAT(/,1X,' ICATA 64 ICATA 68 NLM NPT ',4(1X,I5),
0099          X/,1X,' K KR KK IREM ',4(1X,I5))
0101      IF(IREM.EQ.0)GC TO 550
0102      KKK=K+19
0103  500    DO 500 K2=K,KKK
0104  550    IDATA(K2)=0
0106      IF(NUM.NE.0)ICATA(K)=(NLM/2)+IREM
0108      IF(IHR.NE.0)ICATA(K+1)=(ICATA(K+1)+I+R)/KK
0110      IF(ISBP.NE.0)ICATA(K+2)=(ICATA(K+2)+ISBP)
0112      IF(IDBP.NE.0)ICATA(K+3)=(ICATA(K+3)+IDBP)
0114      IF(IWL.NE.0)ICATA(K+4)=(ICATA(K+4)+IWL)/KK
0116      IF(IRR.NE.0)ICATA(K+5)=(ICATA(K+5)+IRR)/KK
0118      IF(VO2LM.NE.0)CCATA(KR)=(CCATA(KR)+VC2LM)/FLOAT(KK)
0120      IF(VCO2LM.NE.0)COCATA(KR+1)=(COCATA(KR+1)+VCO2LM)/FLOAT(KK)
0122      IF(TVOL.NE.0)CCATA(KR+2)=(CCATA(KR+2)+TVOL)/FLOAT(KK)
0124      IF(AS.NE.0)DCATA(KR+3)=(CCATA(KR+3)+AS)/FLOAT(KK)
0126      IF(ITWL.NE.0.AND.ICATA(12).EQ.1)ICATA(6)=ITWL
0127      INUM=IBLK
0128      RETURN
0129      END

```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
INUM	000016	INTEGER*2 VARIABLE
IBLK	000014	INTEGER*2 VARIABLE
HR	000206	REAL*4 VARIABLE
WL	000212	REAL*4 VARIABLE
AS	000216	REAL*4 VARIABLE
SBP	000222	REAL*4 VARIABLE
DBP	000226	REAL*4 VARIABLE
IAJFLT	000000	INTEGER*2 PROCEDURE
FLOAT	000000	REAL*4 PROCEDURE
TWL	000232	REAL*4 VARIABLE
VOZLM	000236	REAL*4 VARIABLE
VCOZLM	000242	REAL*4 VARIABLE
VOZMKM	000246	REAL*4 VARIABLE
TVOL	000252	REAL*4 VARIABLE
RER	000256	REAL*4 VARIABLE
IRR	000262	INTEGER*2 VARIABLE
RMV02	000264	REAL*4 VARIABLE
RMV02	000270	REAL*4 VARIABLE
VOZMRP	000274	REAL*4 VARIABLE
IREM	000300	INTEGER*2 VARIABLE
MOD	000000	INTEGER*2 PROCEDURE
NNUM	000302	INTEGER*2 VARIABLE
INR	000304	INTEGER*2 VARIABLE
IFIX	000000	INTEGER*2 PROCEDURE
ISBP	000306	INTEGER*2 VARIABLE
IDBP	000310	INTEGER*2 VARIABLE
ITWL	000312	INTEGER*2 VARIABLE
IWL	000314	INTEGER*2 VARIABLE
IREST	000316	INTEGER*2 VARIABLE
IEX	000320	INTEGER*2 VARIABLE
NPT	000322	INTEGER*2 VARIABLE
K	000324	INTEGER*2 VARIABLE
KR	000326	INTEGER*2 VARIABLE
KK	000330	INTEGER*2 VARIABLE
KKK	000332	INTEGER*2 VARIABLE
KZ	000334	INTEGER*2 VARIABLE

COMMON BLOCK /GAS30C/ LENGTH 000016

VOZ	000000	REAL*4 VARIABLE
VCOZ	000004	REAL*4 VARIABLE
VOL	000010	REAL*4 VARIABLE
IBRTH	000014	INTEGER*2 VARIABLE

COMMON BLOCK /PCDATA/ LENGTH 006000

IDATA	000000	INTEGER*2 ARRAY (1536)
DDATA	001000	REAL*4 ARRAY (640)
NUM	000216	INTEGER*2 VARIABLE
RDATA	000220	REAL*4 ARRAY (24)
WT	000146	REAL*4 VARIABLE

FORTAN IV STORAGE MAP

NAME OFFSET ATTRIBUTES

COMMON BLOCK /ADCCM2/ LENGTH 000050

DATBUF	000000	INTEGER*2	ARRAY (20)
AHR	000010	INTEGER*4	VARIABLE
CWL	000014	INTEGER*4	VARIABLE
CWLCNT	000020	INTEGER*2	VARIABLE
CSPD	000024	INTEGER*4	VARIABLE
BPCNT1	000030	INTEGER*2	VARIABLE
CSBP	000032	INTEGER*4	VARIABLE
CDBP	000036	INTEGER*4	VARIABLE
BPCNT	000042	INTEGER*2	VARIABLE
AHRCNT	000000	INTEGER*2	VARIABLE

C-- SUBROUTINE RPTHCR, VERSION 2.0 29 AUG. 1979
C-- PRINTS HEADING FOR 30-SECOND SUMMARY REPORT FOR PROCES

```

0001      SUBROUTINE RPTHCR
0002      INTEGER EXSTRT,RCSTRT
0003      INTEGER AGE,UNIG,TESTRT,WLSLM
0004      REAL WT,HT,ANTEMP,AMPRES
0005      COMMON/PDATA/ICLM1,ISEX,NMSSEC(3),WLSUM,
      @UNIQ(2),NMRTST,ICPIC,IMODE,ITSTYP,
      @IDUM3(2),IHR4,IHR8,IHR12,IHR16,IDUM4,IHRMAX,
      @IDUM5,IBTHDY(3),AGE,NAME(14),ITDATE(3),
      @IDUM6(7),HT,WT,ANTEMP,
      @IDUM7,AMPRES,ICLM8(2),TESTRT(3),
      @NREST,EXSTRT(3),NEXREC,RCSTRT(3),
      @NREC
0006      COMMON/X/MD(24)
      C      DATA MD/'JA','N ','FE','B ','MA','R ','AP','R ',
      C      @'MA','Y ','JUL','N ','JUL','L ','AUG','G ','SE','P ',
      C      @'OC','T ','NOV','V ','DE','C '
      C      * * * PRINT PAGE HEADING * * *
      C      * * * USE TMX OR BX HEADER * * *
      C      ITSTYP(1)=BIKE
      C      ISTYP(2)=TMX
0007 170    IF(ITSTYP.EQ.1) GO TO 210
0009      TYPE 180
0010 180    FORMAT(48X,'EXERCISE RESPONSE TEST:  TREADMILL')
0011      GOTO 190
0012 210    TYPE 200
0013 200    FORMAT(48X,'EXERCISE RESPONSE TEST:  ERGOMETER')
0014 190    TYPE 220
0015 220    FORMAT(50X,'JSC CARDIOPULMONARY LABCRATORY')
0016      TYPE 230
0017 230    FORMAT(56X,'30 SECOND SUMMARY')
0018      IMM=2*(ITDATE(1))
0019      IMM1=IMM-1
      C      * * * PRINT SUBJ INFC * * *
0020      TYPE 240,ICPIC,(NMSSEC(J),J=1,3),
      @ITDATE(2),MD(IMM1),MD(IMM),ITDATE(3),(UNIG(I),I=1,2)
0021 240    FORMAT(/,' CPIC NO: ',I5,I2X,
      @'SS NO: ',I3,'-',I2,'-',I4,I8X,
      @'DATE OF TEST: ',I2,I1X,2A2,I2,
      @I10X,'UNIQUE AC: ',A2,I4)
0022      TYPE 250,TESTRT(1),TESTRT(2),TESTRT(3),WT,AGE
0023 250    FORMAT(1X,'START TIME OF TEST: ',I2,':',I2,':',I2,
      %14X,'WEIGHT(KG) : ',F5.1,
      %14X,'AGE (YRS) : ',I2)
      C      * * * PRINT COLUMN HEADINGS * * *
      C      IF TMX TEST, PRINT TMX HEADINGS
0024      IF(ITSTYP.EQ.1) GOTO 260
0026      TYPE 270
0027 270    FORMAT(/,'XMIN',5X,'HR',4X,'SEP',5X,'DBP',4X,'ELEV',
      @3X,'SPEED',4X,'V02 ---- V02 ----\',3X,'VCO2',5X,'M VCL',
      @3X,'RER',4X,'RESP',5X,'MV/VCO2',2X,'MV/VCO2',3X,
      @'V02/HR')

```

```

0028      TYPE 280
0029      280      FORMAT(' AVG',4X,'BPM',3X,'MM HG',3X,'MM HG',3X,
      @'X GRD',3X,'MPH',5X,'L/NA ML/KG/NA',3X,'L/MIN',
      @4X,'L/MIN',10X,'RATE',22X,'ML/BEAT')
0030      GOTO 290
0031      260      TYPE 300
0032      300      FORMAT('/',5X,'HR',4X,'SBP',5X,'DBP',4X,'POWER',
      @2X,'WL SUM',3X,'--- V02 ---',3X,'VCO2',5X,'M VOL',
      @3X,'RER',4X,'RESP',5X,'MV/VC2',2X,'MV/VC02',3X,
      @'V02/HR')
0033      TYPE 310
0034      310      FORMAT(' AVG',4X,'BPM',3X,'MM HG',3X,'MM HG',3X,
      @'WATTS',3X,'W-MIN',3X,'L/NA ML/KG/NA',3X,'L/MIN',
      @4X,'L/MIN',10X,'RATE',22X,'ML/BEAT')
0035      290      CONTINUE
0036      RETURN
0037      END
    
```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
IMM	001400	INTEGER*2 VARIABLE
IMM1	001402	INTEGER*2 VARIABLE
J	001404	INTEGER*2 VARIABLE
I	001406	INTEGER*2 VARIABLE

COMMON BLOCK /PDATA/ LENGTH 000220

IDUM1	000000	INTEGER*2 VARIABLE
ISEX	000002	INTEGER*2 VARIABLE
NMSSEC	000004	INTEGER*2 ARRAY (3)
WLSUM	000012	INTEGER*2 VARIABLE
UNIQ	000014	INTEGER*2 ARRAY (2)
NMRTST	000020	INTEGER*2 VARIABLE
ICPID	000022	INTEGER*2 VARIABLE
IMODE	000024	INTEGER*2 VARIABLE
ITSTYP	000026	INTEGER*2 VARIABLE
IDUM3	000030	INTEGER*2 ARRAY (2)
IHR4	000034	INTEGER*2 VARIABLE
IHR8	000036	INTEGER*2 VARIABLE
IHR12	000040	INTEGER*2 VARIABLE
IHR16	000042	INTEGER*2 VARIABLE
IDUM4	000044	INTEGER*2 VARIABLE
IHRMAX	000046	INTEGER*2 VARIABLE
IDUM5	000050	INTEGER*2 VARIABLE
IBTHDY	000052	INTEGER*2 ARRAY (3)
AGE	000060	INTEGER*2 VARIABLE
NAME	000062	INTEGER*2 ARRAY (14)
ITDATE	000116	INTEGER*2 ARRAY (3)
IDUM6	000124	INTEGER*2 ARRAY (7)
HT	000142	REAL*4 VARIABLE
WT	000146	REAL*4 VARIABLE
AMTEMP	000152	REAL*4 VARIABLE
IDUM7	000156	INTEGER*2 VARIABLE
AMPRES	000160	REAL*4 VARIABLE
IDUM8	000164	INTEGER*2 ARRAY (2)
TESTRT	000170	INTEGER*2 ARRAY (3)
NREST	000176	INTEGER*2 VARIABLE
EXSTRT	000200	INTEGER*2 ARRAY (3)
NEXREC	000206	INTEGER*2 VARIABLE
RCSTRT	000210	INTEGER*2 ARRAY (3)
NREC	000216	INTEGER*2 VARIABLE

COMMON BLOCK /X/ LENGTH 000060

MD	000000	INTEGER*2 ARRAY (24)
----	--------	----------------------

```
0001      SUBROUTINE TIMSET
          C      PURPOSE: SET TIMERS TO KEEP TRACK OF TIME
          C      OF DAY WITHOUT CHANGING SYSTEM CLOCK
          C      WRITTEN BY WILLIAM CROSIER
          C      DATE 21-NOV-79
          C
0002      INTEGER HR,MN,SC
0003      REAL T0,T1
0004      COMMON/TIM/T1
0005      DATA I1/0./
0006      10      TYPE 110
0007      110     FORMAT(' TIME FROM THE TIME CODE GENERATOR/READER
          X(HR:MN:SC) ?',S)
0008      ACCEPT 120,HR,MN,SC
0009      120     FORMAT(I2,1X,I2,1X,I2)
          C      GET NO. OF SECONDS SINCE SYSTEM BOOT OR CLOCK SET
0010      T0=SECNDS(0.)
0011      IF(HR.LT.0.OR.HR.GT.23.OR.MN.LT.0.OR.MN.GT.59.OR.
          XSC.LT.0.OR.SC.GT.59) GO TO 10
0013      T1=T0-(FLOAT(SC)+60.*FLOAT(MN)+3600.*FLOAT(HR))
0014      RETURN
0015      END
```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
HR	000130	INTEGER*2 VARIABLE
MN	000132	INTEGER*2 VARIABLE
SC	000134	INTEGER*2 VARIABLE
TO	000136	REAL*4 VARIABLE
SECNDS	000000	REAL*4 PROCEDURE
FLOAT	000000	REAL*4 PROCEDURE

COMMON BLOCK /TIM/ LENGTH 000004

T1	000000	REAL*4	VARIABLE
----	--------	--------	----------

```
      C
0001      SUBROUTINE TIMRC(HR,MN,SC)
      C      PURPOSE: READ TIME OF DAY, USING SYSTEM CLOCK AND
      C      PREVIOUSLY CALCLLATED TIME DIFFERENCE BETWEEN
      C      SYSTEM CLOCK TIME AND TIME CODE GEN. TIME
      C
0002      INTEGER HR,MN,SC,REM
0003      COMMON/TIM/T1
0004      T=SECNDS(0.)-T1
0005      HR=IFIX(T/3600.)
0006      REM=IFIX(T-FLCAT(HR)*3600.)
0007      MN=REM/60
0008      SC=REM-(MN*60)
0009      RETURN
0010      END
```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
------	--------	------------

HR	000014	INTEGER*2 PARAMETER VARIABLE
MN	000016	INTEGER*2 PARAMETER VARIABLE
SC	000020	INTEGER*2 PARAMETER VARIABLE
REM	000026	INTEGER*2 VARIABLE
T	000030	REAL*4 VARIABLE
SECND5	000000	REAL*4 PROCEDURE
IFIX	000000	INTEGER*2 PROCEDURE
FLQAT	000000	REAL*4 PROCEDURE

COMMON BLOCK /TIM/ LENGTH 000004

T1	000000	REAL*4 VARIABLE
----	--------	-----------------

```
0001      SUBROUTINE LIGHT(L)
          C      SUBROUTINE LIGHT -- TURNS ON APPROPRIATE PANEL LIGHTS
0002      INTEGER DROUTB
0003      DROUTB = "167772
0004      L = IPEEK(DROUTB).AND."177400 .OR. L
0005      CALL IPOKE(DROUTB,L)
          D      TYPE 1,L
          D1     FORMAT(' L=',Q7)
0006      RETURN
0007      END
```


FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
L	000014	INTEGER*2 PARAMETER VARIABLE
DROUT8	000016	INTEGER*2 VARIABLE
IPEEK	000000	INTEGER*2 PROCEDURE
IPOKE	000000	INTEGER*2 PROCEDURE


```

0001 LOGICAL*1 YES,NO, IANS
0002 LOGICAL*4 OUT,IN
0003 INTEGER NMSSEC(3),NAME(14),ITDATE(3)
0004 INTEGER EXSTRT, IDATA(1536),RCSTRT
0005 INTEGER AGE,UNIG(2),TESTRT,WLSUM
0006 REAL ADATE(3)
0007 REAL WT,RDATA(640)
0008 COMMON /PODATA/ICATA
0009 EQUIVALENCE (IDATA(257),RCATA(1))
0010 EQUIVALENCE (ICATA(2),ISEX),(IDATA(3),NMSSEC(1))
0011 EQUIVALENCE (ICATA(6),WLSUM)
0012 EQUIVALENCE (IDATA(7),LNIG(1)),(IDATA(9),NMRTST),(IDATA(10),ICPID)
0013 EQUIVALENCE (ICATA(11),IMCDE),(ICATA(12),ITSTYP),(IDATA(25),AGE)
0014 EQUIVALENCE (ICATA(26),NAME(1))
0015 EQUIVALENCE (ICATA(40),ITCATE(1)),(IDATA(50),WT),(IDATA(52),WT)
0016 DATA YES/1HY/,NO/1HN/
0017 DATA OUT/'OUT ',IN/'IN '

```

C--

C-- GET PATIENT FILE

C--

0018 3 CALL DATA(IN)

C--

CC-- VERIFY CORRECT SUBJECT READ IN

C--

```

0019 TYPE 4,(NMSSEC(I),I=1,3),LNIG(1),LNIG(2),NMRTST,ICPID
X,(ITDATE(J),J=1,3),ITSTYP
0020 4 FORMAT(/,1X,'VERIFY CORRECT SUBJECT DATA READ FROM DISK',
X//,1X,'S.S. ALN. ',I3,'-',I2,'-',I4,5X,'UNIQUE NO. ',A2,I4,
X5X,'RETEST NO. ',I5,5X,'CPIC NO. ',I4,5X,'TEST DATE ',
XI2,'/',I2,'/',I2,5X,' TEST TYPE ',I2,
X//,'S IS THIS THE CORRECT SUBJECT FILE ?? YES OR NO ?? ')

```

0021 ACCEPT 5,IANS

0022 5 FORMAT(A1)

0023 IF(IANS.NE.YES)GO TO 3

C--

C-- CHECK TO SEE IF EDIT IS REQUESTED

C--

```

0025 TYPE 6
0026 6 FORMAT(/,'S EDITOR REQUIRED ?? YES OR NO ?? ')
0027 ACCEPT 5,IANS
0028 IF(IANS.EQ.YES) CALL EDIT

```

C--

C-- GET CURRENT DATE FROM SYSTEM

C--

0030 CALL DATE(ADATE)

C--

C-- GET THE NUMBER OF REPORTS NEEDED

C--

```

0031 TYPE 12
0032 12 FORMAT(/,'S HOW MANY COPIES OF REPORT NEEDED ? ')
0033 ACCEPT 13,NUMRPT
0034 13 FORMAT(I2)
0035 IF(NUMRPT.LT.1)CALL EXIT

```

C--

```

C-- INSTRUCT THE USER TO ALIGN THE PAGE
C--
0037      TYPE 14
0038 14    FORMAT(/,1X,' PLEASE ALIGN TOP OF NEW PAGE WITH PRINTER HEAD',
           X/, ' TYPE A "RETURN" WHEN READY ')
0039      ACCEPT 5, IANS
0040      DO 530 K=1, NUMRPT
0041      IMLS=0
0042      LINE=18
C          * * * PRINT PAGE HEADING * * *
C          * * * USE TMX OR BX HEADER * * *
0043      WLSUM=ICATA(6)
0044 170    IF(ITSTYP.EQ.1) GO TO 210
0046      TYPE 180
0047 180    FORMAT(48X,'EXERCISE RESPONSE TEST:  TREADMILL')
0048      GOTO 190
0049 210    TYPE 200
0050 200    FORMAT(48X,'EXERCISE RESPONSE TEST:  ERGOMETER')
0051 190    TYPE 220
0052 220    FORMAT(50X,'JSC CARDIOPULMONARY LABCRATCHY')
0053      TYPE 230
0054 230    FORMAT(55X,'ONE MINUTE SUMMARY')
C          * * * PRINT SUBJ INFO * * *
0055      TYPE 240, ICPID, (NMSSEC(J), J=1,3),
           @ (ITDATE(J), J=1,3), (LNIG(I), I=1,2)
0056 240    FORMAT(/, ' CPID NO ', I5, 12X,
           @ 'SS NO. ', I3, '-', I2, '-', I4, 8X,
           @ 'DATE OF TEST ', I2, '/' , I2, '/' ,
           @ I2, 10X, 'UNIQUE NO ', A2, 14)
0057      TYPE 250, (IDATA(J), J=61,63), WT, AGE, WLSUM
0058 250    FORMAT(1X, 'START TIME OF TEST: ', I2, ':', I2, ':', I2,
           X14X, 'WEIGHT      (KG) : ', F5.1,
           X14X, 'AGE (YRS) : ', I2, 7X, 'MCRK LOAD SUM ', I5)
C          * * * PRINT CCLLMN HEADINGS * * *
C          IF TMX TEST, PRINT TMX HEADINGS
0059      IF(ITSTYP.EQ.1) GOTO 260
0061      TYPE 270
0062 270    FORMAT(/, 'XMIN', 5X, 'HR', 4X, 'SBP', 5X, 'DBP', 4X, 'ELEV',
           @ 3X, 'SPEED', 4X, 'V02 --- V02 ---', 3X, 'VCO2', 5X, 'M VCL',
           @ 3X, 'RER', 4X, 'RESP', 5X, 'MV/VCO2', 2X, 'MV/VCO2', 3X,
           @ 'V02/HR')
0063      TYPE 280
0064 280    FORMAT(' AVG', 4X, 'BPM', 3X, 'MM HG', 3X, 'MM HG', 3X,
           @ 'X GRD', 3X, 'MPH', 5X, 'L/MN ML/KG/MN', 3X, 'L/MIN',
           @ 4X, 'L/MIN', 10X, 'RATE', 22X, 'ML/BEAT')
0065      GOTO 290
0066 260    TYPE 300
0067 300    FORMAT(/, 'XMIN', 5X, 'HR', 4X, 'SBP', 5X, 'DBP', 4X, 'POWER',
           @ 2X, 'ML SUM', 3X, 'V02 --- V02 ---', 3X, 'VCO2', 5X, 'M VCL',
           @ 3X, 'RER', 4X, 'RESP', 5X, 'MV/VCO2', 2X, 'MV/VCO2', 3X,
           @ 'V02/HR')
0068      TYPE 310
0069 310    FORMAT(' AVG', 4X, 'BPM', 3X, 'MM HG', 3X, 'MM HG', 3X,
           @ 'WATTS', 3X, 'W-MIN', 3X, 'L/MN ML/KG/MN', 3X, 'L/MIN',

```

```

      24X,'L/MIN',10X,'RATE',22X,'WL/BEAT')
0070 290  CONTINUE
0071      IPT=64
0072      IS=257
0073      IE=IS+3
0074      ISR=4
0075      DO 520 J=1,3
0076      WLS=0.
0077      IF(J.EQ.1)TYPE 500,(IDATA(IT),IT=61,63)
0079 500  FORMAT(/,1X,'* * * REST * * *   REAL TIME: ',I2,':',I2,':',I2)
0080      IF(J.EQ.2)TYPE 501,(IDATA(IT),IT=65,67)
0082 501  FORMAT(/,1X,'* * * EXERCISE * * *   REAL TIME: ',
      X I2,':',I2,':',I2)
0083      IF(J.EQ.3) TYPE 502,(IDATA(IT),IT=69,71)
0085 502  FORMAT(/,1X,'* * * RECOVERY * * *   REAL TIME: ',
      X I2,':',I2,':',I2)
0086      IF(IDATA(IPT).EQ.0)GO TC 515
0088      DO 510 I=1,IDATA(IPT)
0089      IF(IDATA(12).EQ.2)WLS=0.
0091      WL=0.
0092      VO2MKM=0.
0093      RER=0.
0094      RVMC2=0.
0095      RVMCO2=0.
0096      VO2HRR=0.
0097      IHR=IDATA(IS+1)
0098      IF(IDATA(12).EQ.1) WL=FLCAT(IDATA(IS+4))
0100      IF(IDATA(12).EQ.2) WL=FLCAT(IDATA(IS+4))/10.
0102      IF(IDATA(12).EQ.1)WLS=WLS+FLOAT(IDATA(IS+4))
0104      IF(IDATA(12).EQ.2)WLS=RCATA(IS+3)
0106      O2=RCATA(ISR)
0107      CO2=RCATA(ISR+1)
0108      VOL=NDATA(ISR+2)
0109      IF(O2.GT.0)VO2MKM=(O2*1000.)/WT
0111      IF(O2.GT.0)RER=CC2/O2
0113      IRR=IDATA(IS+5)
0114      IF(O2.GT.0)RVMO2=VOL/O2
0116      IF(CO2.GT.0)RVMCO2=VOL/CC2
0118      IF(IHR.GT.0)VC2+RB=(O2*1000.)/FLOAT(IHR)
0120      TYPE 505,(IDATA(IP),IP=IS,IE),WL,WLS,C2,VO2MKM,CO2,VOL,RFP,IRR
      X,RVMO2,RVMCO2,VC2+RB
0121 505  FORMAT(1X,I2,3(1X,I7),F8.1,F7.1,F8.2,F8.1,2X,F8.2,F8.1,F7.2,I7,
      X 1X,3F9.1)
0122      IS=IS+20
0123      IE=IS+3
0124      ISR=ISR+10
0125 510  CONTINUE
0126 515  LINE=LINE+IDATA(IPT)
0127      IPT=IPT+4
0128 520  CONTINUE
0129      TYPE 521,(ADATE(10),ID=1,3)
0130 521  FORMAT(/,94X,'PEPCFT DATE ',2A4,A1)
0131      IF(LINE.GT.66)LINE=LINE-66
0133      LINE=66-LINE

```

```
0134      DO 526,IL=1,LINE
0135      TYPE 525
0136  525  FORMAT(1X)
0137  526  CONTINUE
0138  530  CONTINUE
0139      CALL EXIT
0140      END
```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
------	--------	------------

ADATE	000006	REAL*4 ARRAY (3)
YES	000022	LOGICAL*1 VARIABLE
NO	000023	LOGICAL*1 VARIABLE
IAN5	002612	LOGICAL*1 VARIABLE
OUT	000024	LOGICAL*4 VARIABLE
IN	000030	LOGICAL*4 VARIABLE
EXSTRT	002614	INTEGER*2 VARIABLE
RCSTRT	002616	INTEGER*2 VARIABLE
TESTRT	002620	INTEGER*2 VARIABLE
DATA	000000	REAL*4 PROCEDURE
I	002622	INTEGER*2 VARIABLE
J	002624	INTEGER*2 VARIABLE
EDIT	000000	REAL*4 PROCEDURE
DATE	000000	REAL*4 PROCEDURE
NUMRPT	002626	INTEGER*2 VARIABLE
EXIT	000000	REAL*4 PROCEDURE
K	002630	INTEGER*2 VARIABLE
IWLS	002632	INTEGER*2 VARIABLE
LINE	002634	INTEGER*2 VARIABLE
IPT	002636	INTEGER*2 VARIABLE
IS	002640	INTEGER*2 VARIABLE
IE	002642	INTEGER*2 VARIABLE
ISR	002644	INTEGER*2 VARIABLE
WLS	002646	REAL*4 VARIABLE
IT	002652	INTEGER*2 VARIABLE
WL	002654	REAL*4 VARIABLE
VO2MKM	002660	REAL*4 VARIABLE
RER	002664	REAL*4 VARIABLE
RVM02	002670	REAL*4 VARIABLE
RVMCO2	002674	REAL*4 VARIABLE
VO2HRB	002700	REAL*4 VARIABLE
IHR	002704	INTEGER*2 VARIABLE
FLOAT	000000	REAL*4 PROCEDURE
O2	002706	REAL*4 VARIABLE
CO2	002712	REAL*4 VARIABLE
VOL	002716	REAL*4 VARIABLE
IRR	002722	INTEGER*2 VARIABLE
IP	002724	INTEGER*2 VARIABLE
ID	002726	INTEGER*2 VARIABLE
IL	002730	INTEGER*2 VARIABLE

COMMON BLOCK // LENGTH 000004

COMMON 000000 REAL*4 VARIABLE

COMMON BLOCK /PDATA/ LENGTH 006000

IDATA	000000	INTEGER*2	ARRAY (1536)
RDATA	001000	REAL*4	ARRAY (640)
ISEX	000002	INTEGER*2	VARIABLE
NMSSEC	000004	INTEGER*2	ARRAY (3)
WLSUM	000012	INTEGER*2	VARIABLE

FORTRAN IV

STORAGE MAP

NAME	OFFSET	ATTRIBUTES
UNIQ	000014	INTEGER*2 ARRAY (2)
NMRTST	000020	INTEGER*2 VARIABLE
ICPID	000022	INTEGER*2 VARIABLE
IMODE	000024	INTEGER*2 VARIABLE
ITSTYP	000026	INTEGER*2 VARIABLE
AGE	000060	INTEGER*2 VARIABLE
NAME	000062	INTEGER*2 AFRAY (14)
ITDATE	000116	INTEGER*2 AFRAY (3)
HT	000142	REAL*4 VARIABLE
WT	000146	REAL*4 VARIABLE


```

C--
C-- EDITOR SUBROUTINE
C--
C-- ALLOWS THE USER TO CHANGE ANY OF THE SUBJECTS DATA
C-- STORED IN THE SUBJECT'S DATA FILE.
C--
0001 SUBROUTINE EDIT
C--
C-- SET UP VARIABLES
C--
0002 LOGICAL*1 YES,NC,IRIS
0003 LOGICAL*4 IN,CUT
0004 INTEGER IDATA(1536),IMIN(40)
0005 REAL RDATA(640)
0006 REAL*8 SECT(3)
C--
C-- SET UP COMMON BLOCK
C--
0007 COMMON /PDATA/ICATA
C--
C-- SET UP EQUIVALENCE STATEMENTS
C--
0008 EQUIVALENCE (ICATA(257),RDATA(1))
C--
C-- SET UP DATA VALUES
C--
0009 DATA IN/'IN '//,OUT/'OUT '//
0010 DATA YES/1HY//,NC/1HN/
0011 DATA SECT/'REST ','EXERCISE','RECOVERY'/
C--
C-- GET THE SECTION TO BE EDITED
C--
0012 10 TYPE 100
0013 100 FORMAT(/,1X,'SELECT SECTION TO BE EDITED ',
    &X/,5X,'1 - REST',
    &X/,5X,'2 - EXERCISE',
    &X/,5X,'3 - RECOVERY',
    &X/,'$ TYPE IN THE NUMBER OF THE SECTION TO BE EDITED ')
0014 ACCEPT 200,ISEC
0015 200 FORMAT(I3)
0016 IF(ISEC.EQ.0)GO TO 1000
0018 IF(ISEC.LT.1.OR.ISEC.GT.3)GO TO 10
0020 MREC=IDATA(64+(ISEC-1)*4)
0021 ICNT=1
0022 TYPE 300
0023 300 FORMAT(/,' TYPE IN THE NUMBER(S) OF THE MINUTE(S) TO BE EDITED ')
0024 305 TYPE 310
0025 310 FORMAT('S MINUTE ')
0026 ACCEPT 200,IMIN(ICNT)
0027 IF(IMIN(ICNT).EQ.0)GO TO 390
0029 IF(IMIN(ICNT).LE.MREC)GO TO 350
0031 TYPE 325
0032 325 FORMAT(/,1X,'MINUTE REQUESTED OUT OF RANGE; REQUEST IGNORED')
0033 GO TO 305

```

```

0034 350 ICNT=ICNT+1
0035 IF(ICNT.LT.40)GO TO 305
0037 390 ICNT=ICNT-1
0038 DO 500 IED=1,ICNT
0039 IF(ISEC.EQ.1)IPT=IMIN(IED)
0041 IF(ISEC.EQ.2)IPT=IMIN(IEC)+IDATA(64)
0043 IF(ISEC.EQ.3)IPT=IMIN(IEC)+IDATA(64)+IDATA(68)
0045 IS=257+(IPT-1)*20
0046 ISR=4+(IPT-1)*10
0047 TYPE 395,SECT(ISEC),IMIN(IEC)
0048 395 FORMAT(/,1X,'* * * ',A8,' MINUTE ',I2,' * * *',/)
0049 TYPE 400
0050 400 FORMAT('S HEART RATE ')
0051 ACCEPT 405,IRES
0052 405 FORMAT(A1)
0053 IF(IRES.NE.YES)GO TO 410
0055 TYPE 406,IDATA(IS+1)
0056 406 FORMAT('S OLC = ',I3,' NEW = ')
0057 ACCEPT 407,IDATA(IS+1)
0058 407 FORMAT(I5)
0059 410 TYPE 415
0060 415 FORMAT('S SYSTOLIC BLOOD PRESSURE ')
0061 ACCEPT 405,IRES
0062 IF(IRES.NE.YES)GO TO 420
0064 TYPE 406,IDATA(IS+2)
0065 ACCEPT 407,IDATA(IS+2)
0066 420 TYPE 425
0067 425 FORMAT('S DIASTOLIC BLOOD PRESSURE ')
0068 ACCEPT 405,IRES
0069 IF(IRES.NE.YES)GO TO 430
0071 TYPE 406,IDATA(IS+3)
0072 ACCEPT 407,IDATA(IS+3)
0073 430 IF(IDATA(12).EG.1)TYPE 435
0075 435 FORMAT('S WORK LCAC ')
0076 IF(IDATA(12).EG.2)TYPE 416
0078 416 FORMAT('S ELEVATION ')
0079 ACCEPT 405,IRES
0080 IF(IRES.NE.YES)GO TO 440
0082 ELEV=FLCAT(IDATA(IS+4))/10.
0083 TYPE 417,ELEV
0084 417 FORMAT('S OLC = ',F8.2,' NEW = ')
0085 ACCEPT 418,ELEV
0086 418 FORMAT(F8.0)
0087 IDATA(IS+4)=IFIX(ELEV*10.)
0088 440 TYPE 445
0089 445 FORMAT('S RESPIRATORY RATE ')
0090 ACCEPT 405,IRES
0091 IF(IRES.NE.YES)GO TO 450
0093 TYPE 406,IDATA(IS+5)
0094 ACCEPT 407,IDATA(IS+5)
0095 450 TYPE 455
0096 455 FORMAT('S OXYGEN CONSUMPTION ')
0097 ACCEPT 405,IRES
0098 IF(IRES.NE.YES)GO TO 460

```

```
0100      TYPE 417,RDATA(ISR)
0101      ACCEPT 418,RDATA(ISR)
0102      460      TYPE 465
0103      465      FORMAT('S  CARBON DIOXIDE PRODUCTION ')
0104      ACCEPT 405,IRES
0105      IF(IRES.NE.YES)GO TO 470
0107      TYPE 417,RDATA(ISR+1)
0108      ACCEPT 418,RDATA(ISR+1)
0109      470      TYPE 475
0110      475      FORMAT('S  MINUTE VOLUME ')
0111      ACCEPT 405,IRES
0112      IF(IRES.NE.YES)GO TO 480
0114      TYPE 417,RDATA(ISR+2)
0115      ACCEPT 418,RDATA(ISR+2)
0116      480      TYPE 485
0117      485      FORMAT('S  SPEED  MPH/RPM ')
0118      ACCEPT 405,IRES
0119      IF(IRES.NE.YES)GO TO 500
0121      TYPE 417,RDATA(ISR+3)
0122      ACCEPT 418,RDATA(ISR+3)
0123      500      CONTINUE
0124      GO TO 10
0125      1000     CALL DATA(OUT)
0126      RETURN
0127      END
```

FORTRAN IV

STORAGE MAP

NAME	OFFSET	ATTRIBUTES
IMIN	000014	INTEGER*2 ARRAY (40)
SECT	000134	REAL*8 ARRAY (3)
YES	000174	LOGICAL*1 VARIABLE
NO	000175	LOGICAL*1 VARIABLE
IRES	001310	LOGICAL*1 VARIABLE
IN	000164	LOGICAL*4 VARIABLE
OUT	000170	LOGICAL*4 VARIABLE
ISEC	001312	INTEGER*2 VARIABLE
MREC	001314	INTEGER*2 VARIABLE
ICNT	001316	INTEGER*2 VARIABLE
IED	001320	INTEGER*2 VARIABLE
IPT	001322	INTEGER*2 VARIABLE
IS	001324	INTEGER*2 VARIABLE
ISR	001326	INTEGER*2 VARIABLE
ELEV	001330	REAL*4 VARIABLE
FLOAT	000000	REAL*4 PROCEDURE
IFIX	000000	INTEGER*2 PROCEDURE
DATA	000000	REAL*4 PROCEDURE

COMMON BLOCK /PDATA/ LENGTH 006000

IDATA	000000	INTEGER*2 ARRAY (1536)
RDATA	001000	REAL*4 ARRAY (640)

[illegible]

```

C-----
C      PLOT
C-----
C
C      AUTHOR:  CHARLES MANN
C      DATE:    OCTOBER 19, 1978
C      PURPOSE: TO DRIVE THE PLOTTING ROUTINES FOR COAS
C-----
C
0001      INTEGER IDATA(1536),CPID,ISS(3),IUNG(2),ITDATE(3)
0002      INTEGER IHR(49),SBP(49),CBP(49),ICOPY
0003      INTEGER WL(49)
0004      REAL ICATE(3)
0005      REAL VO2(49),VCC2(49),RDATA(640)
0006      INTEGER IYHR(49),IYVOL(49),IYWL(49),IYSBP(49)
0007      INTEGER IYDBP(49),IYVO2(49),IYVCC2(49),IYRESP(49)
0008      INTEGER IY1(49),IY2(49),IYC2WT(49)
0009      INTEGER*4 IN
0010      COMMON /PDATA/ICATA
0011      EQUIVALENCE(ICATA(257),RCATA(1))
0012      EQUIVALENCE(ICATA(10),CPID)
0013      EQUIVALENCE(ICATA(3),ISS(1))
0014      EQUIVALENCE(ICATA(7),IUNG(1))
0015      EQUIVALENCE(ICATA(12),IMCCE)
0016      EQUIVALENCE(ICATA(40),ITCATE(1))
0017      EQUIVALENCE(ICATA(52),WT)
0018      EQUIVALENCE(ICATA(64),IRST)
0019      EQUIVALENCE(ICATA(68),IEXR)
0020      EQUIVALENCE(ICATA(72),IREC)
0021      DATA IN/'IN '/
C      ZERO OUT PLOTTING BUFFER
0022      DO 10 I=1,49
0023          IHR(I)=0
0024          SBP(I)=0
0025          DBP(I)=0
0026          WL(I)=0
0027          VO2(I)=0.
0028          VCO2(I)=0.
0029      10  CONTINUE
C      GET CURRENT DATE FROM THE SYSTEM
0030      CALL DATE(IDATE)
C      READ SUBJECT DATA IN FROM THE FILE
0031      CALL DATA(IN)
0032      ITOT=IRST+IEXR+IREC
0033      IF (ITOT.GT.49) ITOT=49
0034      DO 20 I=1,ITOT
0035          IHR(I)=IDATA(258+(I-1)*20)
0036          SBP(I)=IDATA(259+(I-1)*20)
0037          DBP(I)=IDATA(260+(I-1)*20)
0038          WL(I)=IDATA(261+(I-1)*20)
0039          VO2(I)=RDATA(4+(I-1)*10)
0040          VCO2(I)=RDATA(5+(I-1)*10)
0041      20  CONTINUE

```

```

C      TYPE 800,WL,ICLN,VC2,RDUM,VC02,RDUM
0800  FORMAT(5(10(1X,I7),/),/,/,5(10(1X,F7.2),/),
C      *      //,5(10(1X,F7.2),/))
C  PLACE DATA INTO ARRAYS SCALED FROM 0 TO 50
0043  DO 100 I=1,49
0044      IYHR(I)=(IDATA(258+(I-1)*20)+2)/4
0045      IYVOL(I)=IFIX((RCATA(6+(I-1)*10)/4.+.5)
0046      IYWL(I)=IFIX((FLCAT(IDATA(261+(I-1)*20))/6.+.5)
0047      IYSBP(I)=(IDATA(259+(I-1)*20)+3)/6
0048      IYDBP(I)=(IDATA(260+(I-1)*20)+3)/6
0049      IYVC2(I)=IFIX((RCATA(4+(I-1)*10))*10.)
0050      IY02WT(I)=IFIX(((RCATA(4+(I-1)*10)*1000.)/WT)+.5)
0051      IYVC02(I)=IFIX((RCATA(5+(I-1)*10))*10.)
0052      IYRESP(I)=(ICATA(262+(I-1)*20)+1)/2
0053  100  CONTINUE
C  PRINT THE CHOICES
0054  110  TYPE 105
0055  105  _FORMAT('0 0--NO PLOTS (STCP)',T30,' 7--SBP VS WL',/,
      *      ' 1--REGULAR 4 PLOTS',T30,' 8--SBP VS HR',/,
      *      ' 2--MIN. VOL. VS TIME',T30,' 9--SBP VS V02',/,
      *      ' 3--RESP. RATE VS TIME',T30,'10--VC02 VS VC2',/,
      *      ' 4--V02 VS TIME(ML/KG-MIN)',T30,'11--MIN. VOL. VS V02',/,
      *      ' 5--V02 VS WL',T30,'12--MIN. VOL. VS VC02',/,
      *      ' 6--HR VS WL',/,/, ' MAKE A SELECTION--',S)
C  READ THE 1ST CHOICE
0056      ACCEPT 115,IPICK1
0057  115  FORMAT(I2)
0058      IF (IPICK1.LE.0) STCP
0060      IF (IPICK1.GT.12) GO TO 110
0062      IF (IPICK1.EQ.1) GO TO 300
C  ASK FOR SECOND CHOICE
0064      TYPE 125
0065  125  FORMAT('0MAKE A SECCND SELECTION--',S)
0066      ACCEPT 115, IPICK2
C  ASK FOR THE NUMBER OF COPIES WANTED
0067      TYPE 305
0068      ACCEPT 315,ICOPY
0069      TYPE 25
0070      ACCEPT 26,I
C--THIS SECTION MOVES THE CORRECT DATA INTO THE FIRST PLOT ARRAY, IY1
0071      GO TO (110,200,210,220,230,240,250,260,270,280,290,298),IPICK1
0072  200  CALL SORT(0,IYVCL,IYVOL,IY1)
0073      GO TO 400
0074  210  CALL SORT(0,IYRESP,IYRESP,IY1)
0075      GO TO 400
0076  220  CALL SORT(0,IYC2WT,IY02WT,IY1)
0077      GO TO 400
0078  230  CALL SORT(1,IYWL,IYV02,IY1)
0079      GO TO 400
0080  240  CALL SORT(1,IYWL,IYHR,IY1)
0081      GO TO 400
0082  250  CALL SORT(1,IYWL,IYSBP,IY1)
0083      GO TO 400
0084  260  CALL SORT(1,IYHR,IYSBP,IY1)

```

```

0085      GO TO 400
0086 270   CALL SORT(1,IYVC2,IYSBP,IY1)
0087      GO TO 400
0088 280   CALL SORT(1,IYVC2,IYVC02,IY1)
0089      GO TO 400
0090 290   CALL SORT(1,IYVC2,IYVOL,IY1)
0091      GO TO 400
0092 298   CALL SORT(1,IYVC02,IYVOL,IY1)
0093 400   GO TO (110,406,410,420,430,440,450,460,470,480,490,49A),IPICK2
0094 406   CALL SORT(0,IYVCL,IYVOL,IY1)
0095      GO TO 500
0096 410   CALL SORT(0,IYRESP,IYRESP,IY2)
0097      GO TO 500
0098 420   CALL SORT(0,IYQZWT,IYQZWT,IY2)
0099      GO TO 500
0100 430   CALL SORT(1,IYWL,IYV02,IY2)
0101      GO TO 500
0102 440   CALL SORT(1,IYWL,IYFR,IY2)
0103      GO TO 500
0104 450   CALL SORT(1,IYWL,IYSBP,IY2)
0105      GO TO 500
0106 460   CALL SORT(1,IYFR,IYSBP,IY2)
0107      GO TO 500
0108 470   CALL SORT(1,IYVC2,IYSBP,IY2)
0109      GO TO 500
0110 480   CALL SORT(1,IYVC2,IYVC02,IY2)
0111      GO TO 500
0112 490   CALL SORT(1,IYVC2,IYVOL,IY2)
0113      GO TO 500
0114 498   CALL SORT(1,IYVC02,IYVOL,IY2)
0115 500   CONTINUE
0116 510   CALL GRAPH(IPICK1,IPICK2,IY1,IY2,ICPIC,ISS,
*          ILNG,ITDATE,IDATE)
0117      TYPE 35
0118      ICOPY=ICOPY-1
0119      IF (ICOPY.GT.0) GO TO 510
0121      GO TO 110
C  FIND OUT HOW MANY COPIES ARE WANTED
0122 300   TYPE 305
0123 305   FORMAT('HOW MANY COPIES DO YOU WANT?--',S)
0124      ACCEPT 315,ICOPY
0125 315   FORMAT(I2)
0126      TYPE 25
0127 25    FORMAT(' ALIGN NEW PAGE WITH PRINT HEAD.',
*             ' HIT RETURN KEY TO PROCEED.')
0128      ACCEPT 26,I
0129 26    FORMAT(I1)
0130 30    CONTINUE
0131      CALL ONE(CPIC,ISS,ILNG,ITDATE,IHR,SEP,
*             DBP,IDATE,ICUM)
0132      TYPE 35
0133 35    FORMAT('0')
0134      CALL TWO(CPIC,ISS,ILNG,ITDATE,WL,VC2,
*             VC02,IDATE,ICUM,IMODE,WT)

```


FORTRAN IV

V01C-03A

FRI 29-FEB-80 14:55:30

PAGE 004

0135 ICOPY=ICOPY-1
0136 TYPE 35
0137 IF (ICOPY.GT.0) GO TO 30
0139 GO TO 110
0140 END

FORTRAN IV

STORAGE MAP

NAME	OFFSET	ATTRIBUTES
IHR	000006	INTEGER*2 ARRAY (49)
SBP	000150	INTEGER*2 ARRAY (49)
DBP	000312	INTEGER*2 ARRAY (49)
WL	000454	INTEGER*2 ARRAY (49)
IDATE	000616	REAL*4 ARRAY (3)
VO2	000632	REAL*4 ARRAY (49)
VCO2	001136	REAL*4 ARRAY (49)
IYHR	001442	INTEGER*2 ARRAY (49)
IYVOL	001604	INTEGER*2 ARRAY (49)
IYWL	001746	INTEGER*2 ARRAY (49)
IYSBP	002110	INTEGER*2 ARRAY (49)
IYDBP	002252	INTEGER*2 ARRAY (49)
IYVO2	002414	INTEGER*2 ARRAY (49)
IYVCO2	002556	INTEGER*2 ARRAY (49)
IYRESP	002720	INTEGER*2 ARRAY (49)
IY1	003062	INTEGER*2 ARRAY (49)
IY2	003224	INTEGER*2 ARRAY (49)
IY02WT	003366	INTEGER*2 ARRAY (49)
ICOPY	004442	INTEGER*2 VARIABLE
IN	003530	INTEGER*4 VARIABLE
I	004444	INTEGER*2 VARIABLE
DATE	000000	REAL*4 PROCEDURE
DATA	000000	REAL*4 PROCEDURE
ITOT	004446	INTEGER*2 VARIABLE
IFIX	000000	INTEGER*2 PROCEDURE
FLOAT	000000	REAL*4 PROCEDURE
IPICK1	004450	INTEGER*2 VARIABLE
IPICK2	004452	INTEGER*2 VARIABLE
SORT	000000	REAL*4 PROCEDURE
GRAPH	000000	REAL*4 PROCEDURE
ICPID	004454	INTEGER*2 VARIABLE
ONE	000000	REAL*4 PROCEDURE
IDUM	004456	INTEGER*2 VARIABLE
TWO	000000	REAL*4 PROCEDURE

COMMON BLOCK /PDATA/ LENGTH 006000

IDATA	000000	INTEGER*2 ARRAY (1536)
RDATA	001000	REAL*4 ARRAY (640)
CPID	000022	INTEGER*2 VARIABLE
ISS	000004	INTEGER*2 ARRAY (3)
IUNG	000014	INTEGER*2 ARRAY (2)
IMODE	000026	INTEGER*2 VARIABLE
ITDATE	000116	INTEGER*2 ARRAY (3)
WT	000146	REAL*4 VARIABLE
IRST	000176	INTEGER*2 VARIABLE
IEXR	000206	INTEGER*2 VARIABLE
IREC	000216	INTEGER*2 VARIABLE

```
0001      SUBROUTINE DATA(ICSTAT)
0002      LOGICAL*1 IFILE(12)
0003      INTEGER FILE(4)
      C
0004      INTEGER*4 IOSTAT,CUT
      C
      C
0005      COMMON /PDATA/ IPATNT(1536) ! THIS IS PATIENT DATA.
      C
      C
0006      DATA OUT/'OUT '
0007      DATA IFILE/'D','X','1',' ' ' ' ' ' ' ' ' ' ' ' ' ' ' 'D','A','T'/
      C
0008      IF(IQSET(2).NE.0) STOP 'CLEUE ALLOCATION FALIURE'
0010      IF(IOSTAT.EQ.CUT)GO TO 10
0012      TYPE 1025
0013 1025  FORMAT('S PLEASE TYPE PATIENT DATA FILE NAME (UP TO 6 SYMBCLS):')
0014      ACCEPT 1030,(IFILE(I),I=4,9)
0015 1030  FORMAT(6A1)
      C
      C
      C
0016 10  CONTINUE
      C
      C--CONVERT ASCII FILE NAME TO RADIX 50
      C
0017      IF(IOSTAT.NE.CUT) CALL IRAC50(12,IFILE,FILE)
      C
      C--OPEN FILE:
      C
0019      ICHAN =IGETC()
0020      IF (ICHAN .LT. 0) STOP 'NO AVAILABLE CHANNEL'
0022      IF (LOOKUP(ICHAN,FILE) .LT. 0) STOP 'FILE ALLOCATION FAILURE'
      C
      C--NOW READS OR WRITES DATA:
      C
0024      IF(IOSTAT.EQ.CUT) GO TO 1050
0026      IF(IREADW(1536,IPATNT,0,ICHAN).LT.0) STOP 'DISK READ ERROR'
0028      GO TO 1060
0029 1050  CONTINUE
0030      IF(IWRITE(1536,IPATNT,0,ICHAN).LT.0) STOP 'DISK WRITE ERROR'
0032 1060  CONTINUE
0033      CALL CLCSEC(ICHAN)
0034      RETURN
0035      END
```

FORTRAN IV

STORAGE MAP

NAME	OFFSET	ATTRIBUTES
IFILE	000016	LOGICAL*1 ARRAY (12)
FILE	000032	INTEGER*2 ARRAY (4)
IOSTAT	000014	INTEGER*4 PARAMETER VARIABLE
OUT	000042	INTEGER*4 VARIABLE
IQSET	000000	INTEGER*2 PROCEDURE
I	000330	INTEGER*2 VARIABLE
IRAD50	000000	INTEGER*2 PROCEDURE
ICHAN	000332	INTEGER*2 VARIABLE
IGETC	000000	INTEGER*2 PROCEDURE
LOOKUP	000000	INTEGER*2 PROCEDURE
IREADW	000000	INTEGER*2 PROCEDURE
IWRITW	000000	INTEGER*2 PROCEDURE
CLOSEC	000000	REAL*4 PROCEDURE

COMMON BLOCK /PDATA/ LENGTH 006000

IPATNT 000000 INTEGER*2 ARRAY (1536)

```

C*****
C
0001      SUBROUTINE ONE(CPID,SS,UNQ,ITDATE,IHR,ISBP,
          @      ICBP,IDATE,ICOPY)
C
C      THIS SUBPROGRAM TAKES ARGUMENTS PASSED
C      FROM MAIN PROGRAM "PLOT" AND GENERATES
C      GRAPHS OF HEART RATE VS TIME AND BLOOD
C      PRESSURES VS TIME ON ONE PAGE OF 15
C      INCH PRINTOUT PAPER.
C
C      TIME IS A CONSTANT 50 MINUTES.
C      ARGUMENTS PASSED:
C      CPID      CARDIOPULMONARY ID NO.
C      SS      SOCIAL SECURITY NUMBER
C      UNQ      UNIQUE TEST NUMBER
C      ITDATE      DATE OF TEST
C      IHR      HEART RATE (BPM)
C      SBP      SYSTOLIC BLOOD PRES.
C      DBP      DIASTOLIC BLOOD PRES.
C      ICOPY      COPIES OF GRAPH REQ'D
C      IDATE      REPORT (TODAY'S) DATE
C
C      RUTT      JUN 78
C
C*****
0002      DIMENSION IDATE(3),IYLABL(11),IBAR(11),
          @IHR(49),I1DASH(49),I1BLNK(49),I2DASH(49),
          @I2BLNK(49),ITIME(50),ISBP(49),ICBP(49),
          @ILIN(49)
0003      REAL IDATE
0004      INTEGER IYLBL2(11),SS(3),ITDATE(3),UNQ(2),CPID
0005      DATA ITIME/18*1H,1HT,1H,1FI,1H,
          @1HM,1H,1HE,25*1H /
0006      DATA ICHAR1/1H-/
0007      DATA ICHAR2/1H /
0008      DATA IBAR/11*1H\ /
0009      DATA I1BLNK/49*' ' /
0010      DATA I2BLNK/49*' ' /
0011      DATA I1DASH/49*'-' /
0012      DATA I2DASH/49*'-' /
0013      DATA IO/1HO/
0014      DATA IX/1HX/
0015      DATA ISTAR/1H* /
0016      DATA IYLABL/0,30,60,90,120,150,180,
          @210,240,270,300/
0017      DATA IYLBL2/0,20,40,60,80,100,
          @120,140,160,180,200/
0018      DATA ILIN/' ',' ',' ',' ',' ',' ',' ',' ',' '5',
          @' ',' ',' ',' ',' ',' ',' ',' ',' '10',
          @' ',' ',' ',' ',' ',' ',' ',' ',' '15',
          @' ',' ',' ',' ',' ',' ',' ',' ',' '20',
          @' ',' ',' ',' ',' ',' ',' ',' ',' '25',

```

```

@:  :  :  :  :  :  : 30:
@:  :  :  :  :  :  : 35:
@:  :  :  :  :  :  : 40:
@:  :  :  :  :  :  : 45:
@:  :  :  :  :  :  :

```

```

0019 10 ICOPY=ICOPY-1
      C   * * * PRINT PAGE HEADING * * *
0020      TYPE 951
0021 951  FORMAT(15X,'JSC CARDIOPULMONARY LABORATORY',
      @45X,'JSC CARDIOPULMONARY LABORATORY')
      C   * * * PRINT GRAPH TITLES * * *
0022      TYPE 900
0023 900  FORMAT(/,18X,'SBP(X) & DBP(C) VS TIME',51X,
      @'HEART RATE VS TIME',/)
0024      TYPE 910
0025 910  FORMAT(18X,'BLOOD PRESSURE (MM HG)',52X,
      @'HEART RATE (BPM)')
      C   * * * PRINT Y LABELS * * *
0026      TYPE 920,IYLAB1,IYLAB2
0027 920  FORMAT(6X,I1,4X,I2,3X,I2,3X,I2,
      @2X,I3,2X,I3,2X,I3,2X,I3,2X,I3,
      @2X,I3,19X,I2,4X,I2,3X,I2,
      @3X,I2,2X,I3,2X,I3,2X,I3,2X,I3,
      @2X,I3,2X,I3,2X,I3)
      C   * * * PRINT Y LABEL MARKERS * * *
0028      TYPE 927,IBAR,IPAR
0029 927  FORMAT(2X,11(4X,1A1),17X,11(4X,1A1))
      C   * * * PRINT TOP LINE OF GRAPHS * * *
0030      TYPE 928,IBAR(1),I1DASH,IBAR(1),
      @IBAR(1),I1DASH,IBAR(1)
0031 928  FORMAT(5X,'0',A1,49A1,A1,
      @20X,'0',A1,49A1,A1)
      C   * * * SET UP THE PLOT * * *
0032      DO 100 I=1,49
0033      IMULT=0
      C   * * * CHECK IF LINE MULTIPLE OF 5 * * *
      C   DIVIDE BY FIVE AND TEST.
      C   EVEN MULTIPLE. THEN USE DASHES IN
      C   LINE ARRAY INSTEAD OF BLANKS.
      C   IF REMAINDER EQUAL 0, LINE IS
0034      XLIN=FLOAT(I)/5.
0035      IXLIN=XLIN
0036      XLDIF=XLIN-IXLIN
0037      IF(XLDIF.GT.0.01) GOTO 140
      C   * * * SET LINE ARRAYS=DASHES * * *
0039      DO 200 JJ=1,49
0040      I1BLNK(JJ)=IC-AR1
0041 200  I2BLNK(JJ)=IC-AR1
      C   * * * SCALE HR TO FIT 50 SPACES * * *
      C   AND ROUND OFF.
      C   IF IHR=0 DO NOT PRINT DATA POINT
      C   PRINT 49 BLANKS OR 49 DASHES
      C   OTHERWISE INSERT AN X TO REPRESENT
      C   HR AVG FOR THAT MINUTE.

```

```
C      * * * TEST IF IHR=0 * * *
0042 140  IF(IHR(I).LT.0) GOTO 150
0044      XDIS=(IHR(I)+2)/4
0045      IXDIS=XDIS
0046      XDIF=XDIS-IXDIS
0047      IF(XDIF.GT.0.5) IXCIS=IXCIS+1
0049      J=IXDIS
0050      I1BLNK(J)=IX

C      * * * SCALE SBP TO FIT 49 SPACES * * *
C      ROUNDOFF. IF SBP=0 DO NOT PRINT A DATA POINT.
C      PRINT 49 BLANKS OR 49 DASHES.
C      INSERT AN "X" IN THE LINE TO REPRESENT
C      SBP FOR THAT MINUTE.
C      * * * TEST FOR SBP=0 * * *?
0051 150  IF(ISBP(I).LT.0) GOTO 160
0053      XDIS1=(ISBP(I)+3)/6
0054      IXD1=XDIS1
0055      XDF1=XDIS1-IXD1
0056      IF(XDF1.GT.0.5) IXC1=IXC1+1
0058      J1=IXD1
0059      I2BLNK(J1)=IX

C      * * * SCALE DBP TO FIT 49 SPACES * * *
C      ROUNDOFF. IF DBP=0 DO NOT PRINT A DATA POINT.
C      PRINT 49 DASHES OR 49 BARS.
C      INSERT AN X IN THE LINE TO REPRESENT
C      DBP FOR THAT MINUTE.
C      * * * TEST FOR DBP=0 * * *
0060 160  IF(ICBP(I).LT.0) GOTO 180
0062      XDIS2=(IDBP(I)+3)/6
0063      IXD2=XDIS2
0064      XDF2=XDIS2-IXD2
0065      IF(XDF2.GT.0.5) IXC2=IXC2+1
0067      J2=IXD2
0068      I2BLNK(J2)=IO

C      * * * TEST IF SBP=DBP * * *
0069      IF(J1.NE.J2) GOTO 180
0071      I2BLNK(J2)=ISTAR

C      * * * LETS PRINT A DATA LINE * * *
0072 180  TYPE 190,ITIME(I),ILIN(I),IBAR(1),I2BLNK,
        @IBAR(1),ITIME(I),ILIN(I),IBAR(1),I1BLNK,IBAR(1)
0073 190  FORMAT(1X,A1,2X,A2,A1,49A1,A1,16X,A1,2X,
        @A2,A1,49A1,A1)

C      * * * RESTORE LINE ARRAYS * * *
0074      DO 300 JJJ=1,49
0075      I1BLNK(JJJ)=ICBAR2
0076 300  I2BLNK(JJJ)=ICBAR2
0077 100  CONTINUE

C      * * * NOW ADD THE BOTTOM LINE * * *
0078      TYPE 195,IBAR(1),I1CASH,IBAR(1),
        @IBAR(1),I1DASH,IBAR(1)
0079 195  FORMAT(4X,'50',A1,49A1,A1,19X,'50',
        @A1,49A1,A1)

C      * * * CK . . . TACK ON THE SURJ INFO * * *
0080      TYPE 800,CPIC,CPIC,SS,SS,UNG,UNG,ITCATE,IDATE
```

```
      & ,ITDATE,ICATE
0081 800  FORMAT(/,5X,'CPID      ',I7,53X,
      @'CPID      ',I7,/,5X,'SS NUMBER : ',
      @I3,'-',I2,'-',I4,48X,'SS NUMBER : ',I3,'-',I2,'-',I4,
      @/,5X,'UNIQUE NC.: ',A2,I4,53X,'UNIQUE
      @ NO.: ',A2,I4,/,5X,'TEST DATE : ',I2,'/',I2,'/',I2,
      @8X,'REPORT DATE : ',2A4,A1,
      @19X,'TEST DATE : ',I2,'/',I2,'/',I2,
      @8X,'REPORT DATE : ',2A4,A1)
0082      TYPE 810
0083 810  FORMAT(20X,'NOTE: AN ASTERISK (*)
      @PLOTTED ON ANY GRAPH INDICATES BOTH
      @ VARIABLES EQUAL AT THAT MINUTE')
0084      RETURN
0085      END
```


FORTRAN IV	STORAGE MAP
NAME	OFFSET ATTRIBUTES
IDATE	000032 REAL*4 PARAMETER ARRAY (3)
IYLABL	000036 INTEGER*2 ARRAY (11)
IBAR	000064 INTEGER*2 ARRAY (11)
IHR	000024 INTEGER*2 PARAMETER ARRAY (49)
I1DASH	000112 INTEGER*2 ARRAY (49)
I1BLNK	000254 INTEGER*2 ARRAY (49)
I2DASH	000416 INTEGER*2 ARRAY (49)
I2BLNK	000560 INTEGER*2 ARRAY (49)
ITIME	000722 INTEGER*2 ARRAY (50)
ISBP	000026 INTEGER*2 PARAMETER ARRAY (49)
IDBP	000030 INTEGER*2 PARAMETER ARRAY (49)
ILIN	001066 INTEGER*2 ARRAY (49)
IYLBL2	001230 INTEGER*2 ARRAY (11)
SS	000016 INTEGER*2 PARAMETER ARRAY (3)
ITDATE	000022 INTEGER*2 PARAMETER ARRAY (3)
UNQ	000020 INTEGER*2 PARAMETER ARRAY (2)
CPID	000014 INTEGER*2 PARAMETER VARIABLE
ICOPY	000034 INTEGER*2 PARAMETER VARIABLE
ICHAR1	001256 INTEGER*2 VARIABLE
ICHAR2	001260 INTEGER*2 VARIABLE
IO	001262 INTEGER*2 VARIABLE
IX	001264 INTEGER*2 VARIABLE
ISTAR	001266 INTEGER*2 VARIABLE
I	002556 INTEGER*2 VARIABLE
IMULT	002560 INTEGER*2 VARIABLE
XLIN	002562 REAL*4 VARIABLE
FLOAT	000000 REAL*4 PROCEDURE
IXLIN	002566 INTEGER*2 VARIABLE
XLDIF	002570 REAL*4 VARIABLE
JJ	002574 INTEGER*2 VARIABLE
XDIS	002576 REAL*4 VARIABLE
IXDIS	002602 INTEGER*2 VARIABLE
XDIF	002604 REAL*4 VARIABLE
J	002610 INTEGER*2 VARIABLE
XDIS1	002612 REAL*4 VARIABLE
IXD1	002616 INTEGER*2 VARIABLE
XDF1	002620 REAL*4 VARIABLE
J1	002624 INTEGER*2 VARIABLE
XDIS2	002626 REAL*4 VARIABLE
IXD2	002632 INTEGER*2 VARIABLE
XDF2	002634 REAL*4 VARIABLE
J2	002640 INTEGER*2 VARIABLE
JJJ	002642 INTEGER*2 VARIABLE

```

C*****
C
0001      SUBROUTINE TWC(CPID,SS,UNG,ITDATE,WL,VO2,
          @          VCO2,IDATE,ICOPY,IMODE,WT)
C
C      THIS SUBPROGRAM TAKES ARGUMENTS PASSED
C      FROM MAIN PROGRAM "PLOT" AND GENERATES
C      GRAPHS OF WORK LOAD VS TIME AND GASES
C      VO2 & VCO2 VS TIME ON ONE PAGE OF 15
C      INCH PRINTOUT PAPER.
C
C      TIME IS A CONSTANT 50 MINUTES.
C      ARGUMENTS PASSED:
C      CPID          CARDIOPULMONARY ID NO.
C      SS            SOCIAL SECURITY NUMBER
C      UNG           UNIQUE TEST NUMBER
C      ITDATE        DATE OF TEST
C      WL            WORK LOAD (WATTS)
C      VO2           OXYGEN CONSUMPTION (L/MIN)
C      VCO2          CARBON DIOXIDE PRODUCTION (L/MIN)
C      ICOPY         COPIES OF GRAPH REQ'D
C      IDATE         REPORT(TODAY'S) DATE
C      IMODE         BICYCLE OR TREADMILL
C
C      RUTT   JUN 78
C
C*****
0002      DIMENSION IDATE(3),IYLABL(11),IBAR(11),
          @WL(49),I1DASH(49),I1BLNK(49),I2DASH(49),
          @I2BLNK(49),ITIME(50),VO2(49),VCO2(49),
          @ILIN(49),IFIV(22)
0003      INTEGER IYLABX(11)
0004      REAL IDATE
0005      INTEGER WL,SS(3),UNG(2),ITDATE(3),CPID
0006      DATA ITIME/18*1H ,1HT,1H ,1HI,1H ,
          @1HM,1H ,1HE,25*1H /
0007      DATA ICHAR1/1H-/
0008      DATA ICHAR2/1H /
0009      DATA IBAR/11*1H\ /
0010      DATA I1BLNK/49*' ' /
0011      DATA I2BLNK/49*' ' /
0012      DATA I1CASH/49*'-' /
0013      DATA I2CASH/49*'-' /
0014      DATA IO/1HO/
0015      DATA IX/1HX/
0016      DATA ISTAR/1H* /
0017      DATA IYLABL/0,30,60,90,120,150,180,
          @210,240,270,300/
0018      DATA IYLABX/0,8,16,24,32,40,48,
          @56,64,72,80/
0019      DATA IFIV/'0.', '0.', '0.', '5.',
          @'1.', '0.', '1.', '5.', '2.',
          @'0.', '2.', '5.', '3.', '0.'

```



```

      C      LINE ARRAY INSTEAD OF BLANKS.
      C      IF REMAINDER EQUAL 0, LINE IS
0049      XLIN=FLOAT(I)/5.
0050      IXLIN=XLIN
0051      XLDIF=XLIN-IXLIN
0052      IF(XLDIF.GT.0.01) GOTO 140
      C      * * * SET LINE ARRAYS=DASHES * * *
0054      DO 200 JJ=1,49
0055      I1BLNK(JJ)=ICAR1
0056      I2BLNK(JJ)=ICAR1
      C      * * * SCALE WL TO FIT 50 SPACES * * *
      C      AND ROUND OFF
      C      IF WL=0 DO NOT PRINT DATA POINT
      C      PRINT 49 BLANKS OR 49 DASHES
      C      OTHERWISE INSERT AN X TO REPRESENT
      C      WL AVG FOR THAT MINUTE.
      C      * * * TEST IF WL=0 * * *
0057      140  IF(WL(I).LT.0) GOTO 150
0059      IF(IMODE.EQ.1)XCIS=(WL(I)+3)/6
0061      IF(IMODE.EQ.2)XCIS=((VC2(I)*1000.)/WT)+.8)/1.6
0063      XDIS=XDIS
0064      XDIF=XDIS-IXDIS
0065      IF(XDIF.GT.0.5) IXDIS=IXDIS+1
0067      J=IXDIS
0068      I1BLNK(J)=IX
      C      * * * SCALE VO2 TO FIT 49 SPACES * * *
      C      ROUNDOFF. IF VC2=0 DO NOT PRINT A DATA POINT.
      C      PRINT 49 BLANKS OR 49 DASHES.
      C      INSERT AN "X" IN THE LINE TO REPRESENT
      C      VO2 FOR THAT MINUTE.
      C      * * * TEST FOR VC2=0 * * *
0069      150  IF(VO2(I).LT.0.01) GOTO 160
0071      XDIS1=VO2(I)*10.
0072      IXD1=XDIS1
0073      J1=IXD1
0074      I2BLNK(J1)=IX
      C      * * * SCALE VCO2 TO FIT 49 SPACES * * *
      C      ROUNDOFF. IF VCC2=0 DO NOT PRINT A DATA POINT.
      C      PRINT 49 DASHES OR 49 BARS.
      C      INSERT AN X IN THE LINE TO REPRESENT
      C      VCO2 FOR THAT MINUTE.
      C      * * * TEST FOR VCO2=0 * * *
0075      160  IF(VCO2(I).LT.0.01) GOTO 180
0077      XDIS2=VCO2(I)*10.
0078      IXD2=XDIS2
0079      J2=IXD2
0080      I2BLNK(J2)=IO
      C      * * * TEST IF VC2=VCC2 * * *
0081      IF(J1.NE.J2) GOTO 180
0083      I2BLNK(J2)=ISTAR
      C      * * * LETS PRINT A DATA LINE * * *
0084      180  TYPE 190,ITIME(I),ILIN(I),IBAR(1),I2BLNK,
      @IBAR(1),ITIME(I),ILIN(I),IBAR(1),I1BLNK,IBAR(1)
0085      190  FORMAT(1X,A1,2X,A2,A1,49A1,A1,16X,A1,2X,

```

```

      @A2,A1,49A1,A1)
      C      * * * RESTORE LINE ARRAYS * * *
0086      DO 300 JJJ=1,49
0087      I1BLNK(JJJ)=ICAR2
0088      300  I2BLNK(JJJ)=ICAR2
0089      100  CONTINUE
      C      * * * NOW ADD THE BOTTOM LINE * * *
0090      TYPE 195,IBAR(1),I1CASH,IBAR(1),
      @IBAR(1),I1DASH,IBAR(1)
0091      195  FORMAT(4X,'50',A1,49A1,A1,19X,'50',
      @A1,49A1,A1)
      C      * * * OK . . . TACK ON THE SUBJ INFO * * *
0092      TYPE 800,CPID,CFID,SS,SS,LNC,UNG,ITCATE,IDATE
      & ,ITCATE,IDATE
0093      800  FORMAT(/,5X,'CPID      :',I7,53X,
      @'CPID      :',I7,/,5X,'SS NUMBER : ',
      @I3,'-',I2,'-',I4,48X,'SS NUMBER : ',I3,'-',I2,'-',I4,
      @/,5X,'UNIQUE NO.: ',A2,I4,53X,'UNIQUE
      @ NO.: ',A2,I4,/,5X,'TEST DATE : ',I2,/,I2,/,I2,
      @8X,'REPORT DATE : ',2A4,A1,
      @19X,'TEST DATE : ',I2,/,I2,/,I2,/,I2,
      @8X,'REPORT DATE : ',2A4,A1)
0094      TYPE 810
0095      810  FORMAT(20X,'NOTE: AN ASTERISK [*]
      @PLOTTED ON ANY GRAPH INDICATES BOTH
      @ VARIABLES EQUAL AT THAT MINUTE')
0096      RETURN
0097      END

```

FORTRAN IV

STORAGE MAP

NAME	OFFSET	ATTRIBUTES
IDATE	000032	REAL*4 PARAMETER ARRAY (3)
IYLABL	000042	INTEGER*2 ARPAY (11)
IBAR	000070	INTEGER*2 ARWAY (11)
WL	000024	INTEGER*2 PARAMETER ARRAY (49)
I1DASH	000116	INTEGER*2 ARRAY (49)
I1BLNK	000260	INTEGER*2 ARRAY (49)
I2DASH	000422	INTEGER*2 ARRAY (49)
I2BLNK	000564	INTEGER*2 ARRAY (49)
ITIME	000726	INTEGER*2 ARRAY (50)
VO2	000026	REAL*4 PARAMETER ARRAY (49)
VC02	000030	REAL*4 PARAMETER ARRAY (49)
ILIN	001072	INTEGER*2 ARRAY (49)
IFIV	001234	INTEGER*2 ARRAY (22)
IYLABX	001310	INTEGER*2 ARRAY (11)
SS	000016	INTEGER*2 PARAMETER ARRAY (3)
UNQ	000020	INTEGER*2 PARAMETER ARRAY (2)
ITDATE	000022	INTEGER*2 PARAMETER ARRAY (3)
CPID	000014	INTEGER*2 PARAMETER VARIABLE
ICOPY	000034	INTEGER*2 PARAMETER VARIABLE
IMODE	000036	INTEGER*2 PARAMETER VARIABLE
WT	000040	REAL*4 PARAMETER VARIABLE
ICHAR1	001336	INTEGER*2 VARIABLE
ICHAR2	001340	INTEGER*2 VARIABLE
IC	001342	INTEGER*2 VARIABLE
IX	001344	INTEGER*2 VARIABLE
ISTAR	001346	INTEGER*2 VARIABLE
I	003002	INTEGER*2 VARIABLE
IMULT	003004	INTEGER*2 VARIABLE
XLIN	003006	REAL*4 VARIABLE
FLOAT	000000	REAL*4 PROCEDURE
IXLIN	003012	INTEGER*2 VARIABLE
XLDIF	003014	REAL*4 VARIABLE
JJ	003020	INTEGER*2 VARIABLE
XDIS	003022	REAL*4 VARIABLE
IXDIS	003026	INTEGER*2 VARIABLE
XDIF	003030	REAL*4 VARIABLE
J	003034	INTEGER*2 VARIABLE
XDIS1	003036	REAL*4 VARIABLE
IXD1	003042	INTEGER*2 VARIABLE
J1	003044	INTEGER*2 VARIABLE
XDIS2	003046	REAL*4 VARIABLE
IXD2	003052	INTEGER*2 VARIABLE
J2	003054	INTEGER*2 VARIABLE
JJJ	003056	INTEGER*2 VARIABLE

CC

C
C SUBROUTINE TO SCRT THE X ARRAY DATA, AVERAGE THE CORRESPONDING
C / ARRAY DATA AND THEN STCR THE AVERAGE VALUE IN THE CORRECT
C WORD OF THE OUT ARRAY.

C AUTHOR : ROY A. REED DATE : 27-OCT-78

CC

0001 SUBROUTINE SORT(SFLG,IXA,IYA,IOA)
0002 INTEGER SFLG,IXA(49),IYA(49),IOA(49)

C
C
D TYPE 800,IXA,SFLG,IYA,SFLG,IOA,SFLG
D200 FORMAT(' ***** SORT ***** ---IXA---',/,
D * S(10(1X,I6)/),/,', ' ---IYA---',/,
D * S(10(1X,I6)/),/,5X,' ---IOA---',/,
C * S(10(1X,I6)/))

0003 DO 100 IZ=1,49 ! LOOP TO ZERO
0004 100 IOA(IZ)=0 ! THE OUTPUT ARRAY

0005 IF(SFLG.EQ.1) GC TC 300 ! TEST TO SEE IF A SCRT IS NEEDED.

0007 DO 200 IT=1,49 ! NO SCRT NEEDED
0008 200 IOA(IT)=IYA(IT) ! DO A DIRECT TRANSFER OF DATA
D TYPE 800,IXA,SFLG,IYA,SFLG,IOA,SFLG
0009 RETURN ! FINSHED RETURN TO CALLING ROUTINE.

0010 300 DO 600 IW=1,49 ! SCRT NEEDED; SET UP LOOP TO SORT-
C ! VALUES FROM SMALLEST TO LARGEST.
0011 AVG=0. ! ZERO THE ACCUMULATOR FOR EQUAL
C ! VALUES.
0012 ICNT=0 ! ZERO THE COUNTER FOR EQUAL VALUES
C ! IN THE X ARRAY

0013 DO 500 IS=1,49 ! LOOP TO GET ALL EQUAL VALUES IN
C ! INCREASING ORDER.
0014 IF(IXA(IS).NE.IW) GC TO 500 ! IF THE X VALUE IS GREATER THAN
C ! THE TEST VALLE BRANCH.
0016 AVG=AVG+FLOAT(IYA(IS)) ! SUMATE ALL EQUAL VALUFS
0017 ICNT=ICNT+1 ! INCREMENT THE COUNTER FOR EQUAL
C ! VALUES.

0018 500 CONTINUE

0019 IF (ICNT.EQ.0) GC TO 600
0021 IOA(IW)=IFIX((AVG/FLOAT(ICNT))+.5) ! CALCULATE THE AVERAGE VALUE
0022 600 CONTINUE

C
C
C
C! AND STORE THE AVERAGE IN
! THE CORRECT WORD IN OUT
! ARRAY.0023
0024RETURN
END

! ALL DONE RETURN TO CALLING ROUTINE

FORTRAN IV

STORAGE MAP

NAME	OFFSET	ATTRIBUTES
IXA	000016	INTEGER*2 PARAMETER ARRAY (49)
IYA	000020	INTEGER*2 PARAMETER ARRAY (49)
IOA	000022	INTEGER*2 PARAMETER ARRAY (49)
SFLG	000014	INTEGER*2 PARAMETER VARIABLE
IZ	000024	INTEGER*2 VARIABLE
IT	000026	INTEGER*2 VARIABLE
IW	000030	INTEGER*2 VARIABLE
AVG	000032	REAL*4 VARIABLE
ICNT	000036	INTEGER*2 VARIABLE
IS	000040	INTEGER*2 VARIABLE
FLOAT	000000	REAL*4 PROCEDURE
IFIX	000000	INTEGER*2 PROCEDURE

```

C
C-----
0001      SUBROUTINE GRAPH(IPICK1,IPICK2,IY1,IY2,ICPID,
*          ISS,IUNG,ITCATE,IDATE)
C-----
C
C      AUTHOR:  CHARLES MANN
C      DATE:    OCTOBER 31, 1978
C      PURPOSE: THIS IS A GENERAL PURPOSE ROUTINE FOR PLOTTING WHATEVER
C      DATA IS PASSED TO IT FOR THE CDAS PLOTS.
C-----
C
0002      INTEGER IY1(49),IY2(49),ISS(3),IUNG(2),ITDATE(3)
0003      INTEGER IHEAD1(12),IHEAD2(12),IXLAB1(22),IXLAB2(22)
0004      INTEGER IYLAB1(22),IYLAB2(22),ICASH1(49),ICASH2(49),IX
0005      INTEGER IBLNK1(49),IBLNK2(49),ICASH,IBLANK
0006      REAL IDATE(3)
0007      DATA IDASH1/49*'-',//,IDASH2/49*'-',//,IBLNK1/49*' '//,IBLNK2/49*' '//
0008      DATA IX/'X',//,ICASH/'-',//,IBLANK/' '//
0009      DATA IBAR/1H\//
D      TYPE 800,IY1,IPICK1,IY2,IPICK2
D800    FORMAT(' ***** GRAPH ***** ',/,5(10(1X,I6),/),/,
C      *      '      '      '---- IY2 ---',/,
D      *      5(10(1X,I6),/))
C      PRINT CLINIC TITLE
0010      TYPE 25
0011      25    FORMAT('0',T19,'JSC CARDIOPULMONARY LABORATORY',
*              T83,'JSC CARDIOPULMONARY LABORATORY',/)
C      GET THE LABELS FOR THE PLOTS
0012      CALL TITLES(IPICK1,IHEAD1,IXLAB1,IYLAB1)
0013      CALL TITLES(IPICK2,IHEAD2,IXLAB2,IYLAB2)
C      PRINT THE PLOT TITLES AND Y AXES
0014      TYPE 35,IHEAD1,IHEAD2,IYLAB1,IYLAB2,IXLAB1(1),IXLAB1(2),
*          IDASH1,IXLAB2(1),IXLAB2(2),ICASH2
0015      35    FORMAT('0',T22,I2A2,T88,I2A2,/,,'0',2X,11(1X,2A2),
*          T66,11(1X,2A2),/,,' ',11(4X,'\'),T67,11(4X,'\'),/,,' ',
*          I2A2,'\ ',49A1,'\ ',T67,2A2,'\ ',49A1,'\ ')
C-----PRINT THE PLOTS-----
0016      DO 40 LINE=1,49
0017      IF (MOD(LINE,5).NE.0) GO TO 20
C      USE A DASHED LINE
0019      IF (IY1(LINE).LE.0) GO TO 5
0021      IF (IY1(LINE).GT.49) IY1(LINE)=49
0023      IDASH1(IY1(LINE))=IX
0024      5    IF (IY2(LINE).LE.0) GO TO 10
0026      IF (IY2(LINE).GT.49) IY2(LINE)=49
0028      IDASH2(IY2(LINE))=IX
0029      10    J=(LINE/5)*2+1
0030      TYPE 15,IXLAB1(J),IXLAB1(J+1),ICASH1,
*          IXLAB2(J),IXLAB2(J+1),ICASH2
0031      15    FORMAT(2X,2A2,'\ ',49A1,'\ ',T67,2A2,'\ ',49A1,'\ ')
0032      IDASH1(IY1(LINE))=ICASH
0033      IDASH2(IY2(LINE))=ICASH
0034      GO TO 40

```

```
C USE A BLANK LINE
0035 20 IF (IY1(LINE).LE.0) GO TO 22
0037 IF (IY1(LINE).GT.49) IY1(LINE)=49
0039 IBLNK1(IY1(LINE))=IX
0040 22 IF (IY2(LINE).LE.0) GO TO 23
0042 IF (IY2(LINE).GT.49) IY2(LINE)=49
0044 IBLNK2(IY2(LINE))=IX
0045 23 TYPE 28, IBLNK1,IBLNK2
0046 28 FORMAT(6X,'\\',49A1,'\\',T71,'\\',49A1,'\\')
0047 IBLNK1(IY1(LINE))=IBLANK
0048 IBLNK2(IY2(LINE))=IBLANK
0049 40 CONTINUE
C PRINT THE LAST LINE AND THE SUBJECT/TEST IDENTIFICATION
C * * * NOW ADD THE BOTTOM LINE * * *
0050 TYPE 195,IXLAB1(21),IXLAB1(22),IBAR,IDASH1,IBAR,
@IXLAB2(21),IXLAB2(22),IBAR,IDASH2,IBAR
0051 195 FORMAT(2X,2A2,51A1,T67,
@2A2,51A1)
C * * * OK . . . TACK ON THE SUBJ INFO * * *
0052 TYPE 50,ICPID,ICPID,ISS,ISS,IUNG,ILNG,ITDATE,IDATE
& ,ITDATE,IC/TE
0053 50 FORMAT(/,5X,'CPID :',I7,53X,
@'CPID :',I7,/,5X,'SS NUMBER : ',
@I3,'-',I2,'-',I4,48X,'SS NUMBER : ',I3,'-',I2,'-',I4,
@/,5X,'UNIQUE NO.: ',A2,I4,53X,'UNIQUE
@ NO.: ',A2,I4,/,5X,'TEST DATE : ',I2,'/',I2,'/',I2,
@8X,'REPORT DATE : ',2A4,A1,
@19X,'TEST DATE : ',I2,'/',I2,'/',I2,
@8X,'REPORT DATE : ',2A4,A1)
0054 RETURN
0055 END
```

FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
IY1	000020	INTEGER*2 PARAMETER ARRAY (49)
IY2	000022	INTEGER*2 PARAMETER ARRAY (49)
ISS	000026	INTEGER*2 PARAMETER ARRAY (3)
IUNG	000030	INTEGER*2 PARAMETER ARRAY (2)
ITDATE	000032	INTEGER*2 PARAMETER ARRAY (3)
IHEAD1	000036	INTEGER*2 AFRAY (12)
IHEAD2	000066	INTEGER*2 ARRAY (12)
IXLAB1	000116	INTEGER*2 ARRAY (22)
IXLAB2	000172	INTEGER*2 ARRAY (22)
IYLAB1	000246	INTEGER*2 ARRAY (22)
IYLAB2	000322	INTEGER*2 ARRAY (22)
IDASH1	000376	INTEGER*2 ARRAY (49)
IDASH2	000540	INTEGER*2 ARRAY (49)
IBLNK1	000702	INTEGER*2 ARRAY (49)
IBLNK2	001044	INTEGER*2 ARRAY (49)
IDATE	000034	REAL*4 PARAMETER AFRAY (3)
IPICK1	000014	INTEGER*2 PARAMETER VARIABLE
IPICK2	000016	INTEGER*2 PARAMETER VARIABLE
ICPID	000024	INTEGER*2 PARAMETER VARIABLE
IX	001206	INTEGER*2 VARIABLE
IDASH	001210	INTEGER*2 VARIABLE
IBLANK	001212	INTEGER*2 VARIABLE
IBAR	001214	INTEGER*2 VARIABLE
TITLES	000000	REAL*4 PROCECURE
LINE	002134	INTEGER*2 VARIABLE
MOD	000000	INTEGER*2 PROCECURE
J	002136	INTEGER*2 VARIABLE

```

0001 SUBROUTINE TITLES(ICODE,IHEAD,IXLABEL,IYLABEL)
C-----
C
C AUTHOR: CHARLES MANN
C DATE OCTOBER 31, 1978
C PURPOSE: TO STORE THE TITLE AND AXIS LABELS IN
C THE APPROPRIATE ARRAYS
C-----
C
0002 INTEGER IHEAD(12,11),IXLABL(22),IYLABEL(22),IYXVAL(2,11)
0003 INTEGER ILABEL(22,5),IHEAC(12)
0004 DATA IHEAD1/
1      'MI','N','VO','L','VS','T','IM','E'
2      'R','ES','P','RA','TE','V','S','TI','ME'
3      'VO','2','VS','T','IM','E','(M','L','KG','-M','IN')
4      'V','C2','V','S','WL'
5      'H','V','S','WL'
6      'S','BP','V','S','WL'
7      'S','BP','V','S','HR'
8      'S','BP','V','S','VO','2'
9      'V','CO','2','VS','V','O2'
A      'M','IN','V','OL','V','S','VO','2'
B      'MI','N','VO','L','VS','V','CO','2'
0005 DATA ILABEL/
1      '0','5','10','15','20','25'
&      '30','35','40','45','50'
2      '0','20','40','60','80','1','00'
&      '1','20','1','40','1','60','1','80','2','00'
3      '0','30','60','90','1','20','1','50'
&      '1','80','2','10','2','40','2','70','3','00'
4      '0','0','0','5','1','0','1','5','2','0','2','.5'
&      '3','0','3','5','4','0','4','5','5','0'
5      '0','10','20','30','40','50'
&      '60','70','80','90','1','00'
0006 DATA IYXVAL/2,1,5,1,1,1,4,3,2,3,3,3,
*          3,2,3,4,4,4,2,4,2,4/
D TYPE 800,ICODE
C800 FORMAT(' ***** TITLES ***** ICODE=',I3)
0007 JCODE=ICODE-1
0008 DO 100 I=1,12
0009 100 IHEAD(I)=IHEAD1(I,JCODE)
0010 DO 200 I=1,22
0011 IXLABL(I)=ILABEL(I,IYXVAL(2,JCODE))
0012 IYLABEL(I)=ILABEL(I,IYXVAL(1,JCODE))
0013 200 CONTINUE
0014 RETURN
0015 END

```

FORTRAN IV

STORAGE MAP

NAME	OFFSET	ATTRIBUTES
IHEAD1	000024	INTEGER*2 ARRAY (12,11) VECTORED
IXLABL	000020	INTEGER*2 PARAMETER ARRAY (22)
IYLABL	000022	INTEGER*2 PARAMETER ARRAY (22)
IYXVAL	000434	INTEGER*2 ARRAY (2,11)
ILABEL	000510	INTEGER*2 ARRAY (22,5) VECTORED
IHEAD	000016	INTEGER*2 PARAMETER ARRAY (12)
ICODE	000014	INTEGER*2 PARAMETER VARIABLE
JCODE	001104	INTEGER*2 VARIABLE
I	001106	INTEGER*2 VARIABLE

```

*****
*****
*
*
*          TECHNOLOGY INCORPORATED
*          LIFE SCIENCE DIVISION
*
*****
*
* PROGRAM NAME:..... EDICON
* AUTHOR:..... RCY A. REED
* DATE:..... 2/NOV/78
*
-----
*
* COMPUTER SYSTEM:..... DEC PDP 1103 "LSI-11"
* OPERATING SYSTEM:..... RT-11  VO1C-03A
*
-----
*
* COMPILING SEQUENCE:
*
*   .R FORTRA <CR>
*   *EDICON=EDICON <CR>
*   **<CR>
*
-----
*
* RUN MODUAL LINKING SEQUENCE:
*
*   .R LINK <CR>
*   *EDICON=EDICON,DX0:SYSLIB/F <CR>
*   **<CR>
*
-----
*
* CALLING SEQUENCE:
*
*   .P EDICON <CR>
*
-----
*
* PURPOSE:
*
*   ALLOWS THE USER TO MODIFY THE CONSTANTS USED BY THE
*   CALIBRATION ROUTINE IN THE MAIN PROGRAM PROCES.
*   THE CONSTANTS ARE STORED IN THE DISK FILE CALCON.DAT.
*
*****
*****
CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC

```

```

0001 LOGICAL*1 IANS,YES,NO
0002 INTEGER FILE(4)
0003 REAL CALBUF(256)
0004 DATA YES/1HY/,NO/1HN/
0005 CALL IRAD50(12,'DXOCALCONCAT',FILE)
0006 ICHAN=IGETC()
0007 IF(LOOKUP(ICCHAN,FILE).LT.0) STOP 'LOOKUP ERROR'
0009 IF(IREADW(256,CALBUF,1,ICCHAN).LT.0) STOP 'DISK READ ERROR'
0011 TYPE 2
0012 2 FORMAT(/,'$ LIST OF CURRENT VALUES NEEDED ?? YES OR NO ??')
0013 ACCEPT 3,IANS
0014 3 FORMAT(A1)
0015 IF(IANS.NE.YES)GO TO 8
0017 DO 5 I=1,22
0018 TYPE 4,I,CALBUF(I)
0019 4 FORMAT(/,1X,'CALIBRATION FOR WORD ',I2,' IS ',F10.2)
0020 5 CONTINUE
0021 8 TYPE 6
0022 6 FORMAT(/,'$ TYPE IN NUMBER OF WORD TO EDIT ')
0023 ACCEPT 7,I
0024 7 FORMAT(I3)
0025 IF(I.EQ.0)GO TO 20
0027 TYPE 10,I,CALBUF(I)
0028 10 FORMAT(/,'$ WORD ',I4,' ) OLD = ',F10.2,' NEW = ',)
0029 ACCEPT 15,CALBUF(I)
0030 15 FORMAT(F10.0)
0031 GO TO 8
0032 20 CONTINUE
0033 IF(IWRITW(256,CALBUF,1,ICCHAN).LT.0)STOP 'DISK WRITE ERROR'
0035 CALL CLOSEC(ICCHAN)
0036 STOP
0037 END

```


FORTRAN IV STORAGE MAP

NAME	OFFSET	ATTRIBUTES
FILE	000006	INTEGER*2 ARRAY (4)
CALBUF	000016	REAL*4 ARRAY (256)
IAN5	002412	LOGICAL*1 VARIABLE
YES	002016	LOGICAL*1 VARIABLE
NO	002017	LOGICAL*1 VARIABLE
IRAD50	000000	INTEGER*2 PROCEDURE
ICHAN	002414	INTEGER*2 VARIABLE
IGETC	000000	INTEGER*2 PROCEDURE
LOOKUP	000000	INTEGER*2 PROCEDURE
IREADW	000000	INTEGER*2 PROCEDURE
I	002416	INTEGER*2 VARIABLE
IWRITW	000000	INTEGER*2 PROCEDURE
CLOSEC	000000	REAL*4 PROCEDURE

APPENDIX I
FORMAT OF THE FLOPPY DISK FILE

PRECEDING PAGE BLANK NOT FILMED

CDAS PATIENT DATA FILE - FIRST BLOCK - PRETEST DATA

DUM 1	SUBJECT SEX		SOCIAL		SECURITY		NUMBER		WLSUM		UNIQUE		NUMBER		RETEST		CPID		TEST MODE 1=LIVE 2=TAPE		TYPE OF TEST B=1, T=2	
	1	A	2	I	3	I	4	I	5	I	6	A	7	I	8	I	9	I	10	I	11	I

MAJOR REV. #	MINOR REV. #		TARGET HR @4		TARGET HR @8		TARGET HR @12		TARGET HR @16		DUM 4		MAX TARGET HR		DUM 5		DATE OF MONTH		BIRTH DAY		YEAR	
	13	I	14	I	15	I	16	I	17	I	18	I	19	I	20	I	21	I	22	I	23	I

AGE	SUBJECT NAME (14 WORDS, 28 CHARACTERS)	
	25	A
26	A	
27	A	
28	A	
29	A	
30	A	
31	A	
32	A	
33	A	
34	A	
35	A	
36	A	

SUBJECT	NAME CONT'D		TEST DATE		CALIBRATION		DUM 6 (16 WORDS)																
	37	A	38	A	39	I	40	I	41	I	42	I	43	I	44	I	45	I	46	I	47	I	48

DUM 6	HEIGHT (CM)		WEIGHT (KG)		AMBIENT TEMP. (°C)		DUM 7		AMBIENT (mm Hg)		PRES.		DUM 8 (2 WORDS)										
	49	F	50	F	51	F	52	F	53	F	54	F	55	F	56	F	57	F	58	F	59	F	60

TIME BEGIN REST	NO. OF REST RECORDS		TIME BEGIN EXERCISE		NO. OF EXERCISE RECORDS		TIME BEGIN RECOVERY		NO. OF RECOVERY RECORDS														
	61	I	62	I	63	I	64	I	65	I	66	I	67	I	68	I	69	I	70	I	71	I	72

FIRST BLOCK (CONTINUED) - CALIBRATION & PFT DATA

STPD FACTOR		BTPS	FACTOR	AMBIENT	O ₂ %	AMBIENT	N ₂ %	AMBIENT	CO ₂ %	DUM 9 (2 WORDS)													
F	73	F	74	F	75	F	76	F	77	F	78	F	79	F	80	F	81	F	82	F	83	F	84

[illegible]

CO2 SLOPE	CO2 Y-INT	HR-SLOPE	HR-Y-INT	WL/ELEV SLOPE	WL/ELEV	Y-INT
F 97	F 98	F 100	F 102	F 104	F 105	F 106

[illegible]

	FVC	FEV ₁
127	128	130
126	129	131
125	128	130
124	127	129
123	126	128
122	125	127
121	124	126

PEFR	MWFR	MEFR	% FEV1/FVC	
F 133	F 134	F 136	F 140	F 143
F 135	F 137	F 138	F 142	F 144

[illegible][illegible][illegible][illegible][illegible][illegible]

Abstract

[illegible][illegible][illegible]

End of last record →

TRAILER (all zeros-words 241-256)	
241	242
243	244
245	246
247	248
249	250
251	252
253	254
255	256

TRAILER CONT.					
249		247	248	249	250
					251
					252
					253
					254
					255
					256